



Surgical Technique

about us

Auxin Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements



Guidelines

This publication sets forth detailed recommended procedures for using Auxein Medical devices and instruments.

It offers guidance that needs to be heeded. However, with any such technical guide, each surgeon must consider the unique needs of each patient and make appropriate adjustments when and as required.

A workshop training under DAIS Academy by Auxein will provide assistance prior to first surgery. It is vital to know that all non-sterile devices must be cleaned and sterilized before use.

Moreover, multi-component instruments must be disassembled for cleaning. The surgeon must discuss all relevant risks, including the finite lifetime of the device, with the patient, when necessary.

Please NOTE that all the bone screws referenced in this document here are not approved for screw attachment or fixation in the areas not mentioned in this publication.

Warning:

This description is not sufficient for immediate application of the instrumentation. Instruction by a surgeon experienced in handling this instrumentation is highly recommended.



Contents

Product Overview	01
AO Principles	02
Plate Feature	03
Variable Angle Locking	04
Indications	05
Surgical Techniques Steps	07-11
Surgical Implants Composition	12-15
Screws	16-17
Implant Box Composition	18-28
Instruments	29-38

INTRODUCTION

AUXEIN MEDICAL'S Variable Angle System consists of variety of Bone plates & fixed angle Bone Screws. This system is a single use implantable device for long term duration (intended for continuous use for more than 30 days) contacting radius and Humerus bone and its surrounding tissues. Variable Angle System plates are having fixed angle and variable angle holes. Variable locking Holes allow up to 15° off-axis screw angulation in all directions in order to address the individual fracture patterns. Properties of fixed angle plates enable their successful using even in less quality and osteoporotic bones. It is mainly useful during intra-articular fractures treatment.




AO PRINCIPLES

Anatomic reduction

Fracture reduction and fixation to restore anatomical relationships.

Early, active mobilization

Early and safe mobilization and rehabilitation of the injured part and the patient as a whole.



Stable fixation


Fracture fixation providing absolute or relative stability, as required by the patient, the injury, and the personality of the

Preservation of blood supply

Preservation of the blood supply to soft tissue and bone by gently reduction technique and careful handling.


Plate Features

- The AV-Wiselock holes combined with variable angle locking screws allow up to 15° off-axis angulation in all directions.
- The variable angle locking screws offer a fixed-angle construct to support the articular surface
- Fixation can be achieved in osteoporotic bone
- Manufactured in stainless steel & titanium




Variable Angle Locking


- Screws can be angled anywhere within 30° cone around the central axis of the plate hole. (Figure 1)
- Four columns of threads in the Variable Angle locking hole provide four points of threaded locking between the Variable Angle Plate and the Variable Angle screw, forming a fixed-angle construct at the desired screw angle. (Figure 3)




(1)



(2)



(3)




INDICATIONS

The 2.4mm Variable Angle Distal Radius Plates are indicated for fixation of complex intra- and extra-articular fractures and osteotomies of the distal radius and other small bones in adults, skeletally mature adolescents, and the following adolescent distal radius fractures:



- Intra-articular fractures exiting the epiphysis
- Intra-articular fractures exiting the metaphysis
- Physeal crush injuries
- Any injuries which cause growth arrest to the distal radius

2.4mm Variable Angle Two-Column Volar Distal Radius Plate



1. Reduce Fracture and Position Plate.

Reduce the fracture. The reduction method will be fracture specific. Beginning with the elongated hole in the shaft of the plate, drill with the 1.8 mm drill bit using the Universal Drill Guide 2.4.

Insert a 2.4 mm Cortical screw in the elongated hole in the plate shaft. Adjust the plate position if necessary and tighten the screw.

The order of screw insertion in the shaft and metaphysis may vary depending on the fracture pattern and reduction technique.

Instrumentation:-

Codes	Names
1312-04-1.8	Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm
1312-22	Quick Coupling Handle
1312-18	Star Screwdriver Shaft, T8
1312-17	Depth Gauge with Protector measuring upto 50mm
1312-15	Self-Centering Double Drill Guide Ø1.8/2.4mm

Alternative

Codes	Names
1312-16	Self-Centering Double Drill Guide Ø2.0/2.7mm
1312-05-2.0	Drill Bit with Quick Coupling End, Ø2.0mm x Length 140mm
1312-01	Depth Gauge measuring upto 60mm

Alternatively, insert 2.7 mm Cortical screws into the shaft. Use the Universal Drill Guide 2.7 in the unthreaded part of the hole. Drill with the drill bit 2.0 mm.

Precaution: Reverse bending or use of the incorrect instrumentation for bending may weaken the plate and lead to premature plate failure (e.g. breakage). Do not bend the plate beyond what is required to match the anatomy



2. Insert Proximal Screw

Determine where 2.4 mm VA locking screws or 2.4 mm cortical screws will be used in the shaft of the plate. Insert these screws beginning with the most proximal screw.


For VA locking screws, carefully insert the threaded drill sleeve 1.8/ 2.4 in line with the hole's axis until it is seated in the desired locking hole. Use the 1.8 mm drill bit.

Read the screw length directly from the laser mark on the drill bit, or use the corresponding depth gauge to determine the screw length. To insert the screw use the T8 screwdriver shaft in combination with the torque limiter 0.8 Nm

Precaution: Use of the Torque limiting attachment is mandatory when inserting locking screws into variable angle locking holes, to ensure the adequate torque is applied.

Instrumentation:-

Codes	Names
1312-04-1.8	Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm
1312-13-1.8	Threaded Drill Sleeve, Ø1.8/2.4mm (Double Start Thread)
1312-01	Depth Gauge measuring upto 60mm
1312-18	Star Screwdriver Shaft, T8
1312-21	Star Screwdriver Shaft, T8
7-019-03	Handle for Torque Limiting Attachment
1312-22	Quick Coupling Handle



Note: For dense bone, visually inspect if the screw is countersunk after tightening with the torque limiter. If required, carefully tighten without the torque limiter until the screw head is flush with the plate surface.

For 2.4 mm Cortical Screws use the Universal Drill Guide 2.4 and drill with the 1.8 mm drill bit. Measure the screw length with the corresponding depth gauge. Insert the screw with the T8 screwdriver shaft.

3. Insert AV-Wiselock Screw



The AV-Wiselock drill guide tip is inserted and keyed into the cloverleaf design of the AV-Wiselock hole.

Use the funnel-shaped end of the VA-LCP drill guide to drill variable angle holes at the desired angle.

Precaution: Do not use the threaded LCP drill sleeve in variable angle locking holes.

Instrumentation:-

Codes	Names
1312-18	Star Screwdriver Shaft, T8
1312-04-1.8	Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm
1312-30	Universal Variable Angle Drill Guide, Ø1.8mm
1312-22	Quick Coupling Handle



4. Drilling for AV-Wiselock Screw

When the Universal Variable Angle Drill Guide, Ø1.8mm is engaged in the AV-Wiselock locking hole, use the 1.8 mm drill bit to drill to the desired depth at the desired angle.

The funnel of the drill guide allows the drill bit a total angle variation of 30°. Verify the drill bit angle under C-arm to ensure the desired angle has been achieved. If necessary, drill at a different angle and verify again under C-arm. Use the depth gauge for 2.0 mm and 2.4 mm screws to measure the correct screw length.

The fixed-angle end of the drill guide only allows the drill bit to follow the nominal trajectory of the locking hole.



5. Preliminary Screw Placement

Insert the AV-Wiselock Screw screws manually with the self retaining T8 Stardriver screwdriver shaft and Quick Coupling handle until just before the screw head is seated in the locking hole.


Do not over-tighten the screw so they can be easily removed if necessary.

Instrumentation:-

Codes	Names
1312-22	Quick Coupling Handle
1312-18	Star Screwdriver Shaft, T8

Optional Instrumentation:-

Codes	Names
1312-31	Screwdriver Star Head



6. Confirm Proper Joint Reconstruction

Confirm proper joint reconstruction, screw placement and screw length using multiple C-arm views.

Ensure that the distal screws are not in the joint using additional views such as oblique view.



7. Locking of Variable Angle Screws in AV-Wiselock Screw Holes

Use the torque limiting attachment to perform the final locking step for the VA locking screws.

Precaution: Use of the torque limiter is mandatory when inserting locking screws into variable angle locking holes, to ensure the adequate torque is applied.

Note: For dense bone, visually inspect if the screw is countersunk after tightening with the torque limiter. If required, carefully tighten without the torque limiter until the screw head is flush with the plate surface



Instrumentation:-

Codes	Names
1312-21	Torque Limiting Attachment, 0.8Nm
7-019-03	Handle for Torque Limiting Attachment
1312-18	Star Screwdriver Shaft, T8
1312-22	Quick Coupling Handle

Use the appropriate method for surgical closure of the incision. Closing of Incision

8. Implant Removal

To remove locking screws, first unlock all screws from the plate; then remove the screws completely from the bone.

This prevents rotation of the plate when removing the last locking screw.



2.4mm

Distal Radius System

- AV-Wiselock Dorsal System**
- AV-Wiselock Volar Extra-articular System**
- Variable Angle Volar Column System**
- Variable Angle Volar Rim System**


2.4mm AV-Wiselock Dorsal Distal Radius Plate, Radial Column

Holes	Stainless Steel	Titanium
5	11-045-05SS	11-045-05TI
6	11-045-06SS	11-045-06TI




2.4mm AV-Wiselock Dorsal Distal Radius Plate, Intermediate Column, (2Head Holes)

Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	11-046-03LSS	11-046-03LTI	11-046-03RSS	11-046-03RTI
4	11-046-04LSS	11-046-04LTI	11-046-04RSS	11-046-04RTI




2.4mm AV-Wiselock "L" Dorsal Distal Radius Plate, (2Head Holes)

Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	11-047-03LSS	11-047-03LTI	11-047-03RSS	11-047-03RTI
5	11-047-05LSS	11-047-05LTI	11-047-05RSS	11-047-05RTI




2.4mm AV-Wiselock "L" Dorsal Distal Radius Plate, (3 Head Holes)

Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	11-048-03LSS	11-048-03LTI	11-048-03RSS	11-048-03RTI
5	11-048-05LSS	11-048-05LTI	11-048-05RSS	11-048-05RTI



2.4mm AV-Wiselock "L" Dorsal Distal Radius Plate, Oblique (3 Head Holes)

Holes	Left Direction		Right Direction	
	Stainless Steel	Titanium	Stainless Steel	Titanium
3	11-049-03LSS	11-049-03LTI	11-049-03RSS	11-049-03RTI
5	11-049-05LSS	11-049-05LTI	11-049-05RSS	11-049-05RTI



2.4mm AV-Wiselock "T" Dorsal Distal Radius Plate, (3 Head Holes)

Holes	Stainless Steel	Titanium
3	11-050-03SS	11-050-03TI
5	11-050-05SS	11-050-05TI



2.4mm AV-Wise Distal Radius Volar Plate, Extra-Articular, (5 Head Holes)

Stainless Steel	Titanium	Head Hole	Shaft Hole	Direction
11-043-03LSS	11-043-03LTI	5	3	Left
11-043-03RSS	11-043-03RTI	5	3	Right
11-043-05LSS	11-043-05LTI	5	5	Left
11-043-05RSS	11-043-05RTI	5	5	Right




2.4mm AV-Wise Distal Radius Volar Plate, Extra-Articular, (4 Head Holes)

Stainless Steel	Titanium	Head Hole	Shaft Hole	Direction
11-044-03LSS	11-044-03LTI	4	3	Left
11-044-03RSS	11-044-03RTI	4	3	Right
11-044-05LSS	11-044-05LTI	4	5	Left
11-044-05RSS	11-044-05RTI	4	5	Right




2.4mm Variable Angle Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Dia	Length
549.006	TI-549.006	Ø2.4	6mm
549.008	TI-549.008	Ø2.4	8mm
549.010	TI-549.010	Ø2.4	10mm
549.012	TI-549.012	Ø2.4	12mm
549.014	TI-549.014	Ø2.4	14mm
549.016	TI-549.016	Ø2.4	16mm
549.018	TI-549.018	Ø2.4	18mm
549.020	TI-549.020	Ø2.4	20mm
549.022	TI-549.022	Ø2.4	22mm
549.024	TI-549.024	Ø2.4	24mm
549.026	TI-549.026	Ø2.4	26mm
549.028	TI-549.028	Ø2.4	28mm
549.030	TI-549.030	Ø2.4	30mm




2.4mm Cortical Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Dia	Length
1144.06	TI-1144.06	Ø2.4	6mm
1144.08	TI-1144.08	Ø2.4	8mm
1144.10	TI-1144.10	Ø2.4	10mm
1144.12	TI-1144.12	Ø2.4	12mm
1144.14	TI-1144.14	Ø2.4	14mm
1144.16	TI-1144.16	Ø2.4	16mm
1144.18	TI-1144.18	Ø2.4	18mm
1144.20	TI-1144.20	Ø2.4	20mm
1144.22	TI-1144.22	Ø2.4	22mm
1144.24	TI-1144.24	Ø2.4	24mm
1144.26	TI-1144.26	Ø2.4	26mm
1144.28	TI-1144.28	Ø2.4	28mm
1144.30	TI-1144.30	Ø2.4	30mm




2.7mm Cortical Screws, Self-Tapping (Star Head)

Stainless Steel	Titanium	Dia	Length
1242.06	TI-1242.06	Ø2.7	6mm
1242.08	TI-1242.08	Ø2.7	8mm
1242.10	TI-1242.10	Ø2.7	10mm
1242.12	TI-1242.12	Ø2.7	12mm
1242.14	TI-1242.14	Ø2.7	14mm
1242.16	TI-1242.16	Ø2.7	16mm
1242.18	TI-1242.18	Ø2.7	18mm
1242.20	TI-1242.20	Ø2.7	20mm
1242.22	TI-1242.22	Ø2.7	22mm
1242.24	TI-1242.24	Ø2.7	24mm
1242.26	TI-1242.26	Ø2.7	26mm
1242.28	TI-1242.28	Ø2.7	28mm
1242.30	TI-1242.30	Ø2.7	30mm



2.4mm Wise-Lock Screw, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Dia	Length
1145.06	TI-1145.06	Ø2.4	6mm
1145.08	TI-1145.08	Ø2.4	8mm
1145.10	TI-1145.10	Ø2.4	10mm
1145.12	TI-1145.12	Ø2.4	12mm
1145.14	TI-1145.14	Ø2.4	14mm
1145.16	TI-1145.16	Ø2.4	16mm
1145.18	TI-1145.18	Ø2.4	18mm
1145.20	TI-1145.20	Ø2.4	20mm
1145.22	TI-1145.22	Ø2.4	22mm
1145.24	TI-1145.24	Ø2.4	24mm
1145.26	TI-1145.26	Ø2.4	26mm
1145.28	TI-1145.28	Ø2.4	28mm
1145.30	TI-1145.30	Ø2.4	30mm



2.4mm

Implant Box Composition

**11-051 : Implant Box with cover for 2.4mm AV-Wiselock
Distal Radius Dorsal and Extra-articular System**



2.4mm AV-Wiselock Dorsal Distal Radius Plate, Radial Column

Stainless Steel	Titanium	Shaft Hole	Units
11-045-05SS	11-045-05TI	5	1
11-045-06SS	11-045-06TI	6	1

2.4mm AV-Wiselock Dorsal Distal Radius Plate, Intermediate Column, (2 Head Holes)

Left		Right				
Stainless Steel	Titanium	Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
11-046-03LSS	11-046-03LTI	11-046-03RSS	11-046-03RTI	3	2	1
11-046-04LSS	11-046-04LTI	11-046-04RSS	11-046-04RTI	4	2	1

2.4mm AV-Wiselock "L" Dorsal Distal Radius Plate, (2 Head Holes)

Left		Right				
Stainless Steel	Titanium	Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
11-047-03LSS	11-047-03LTI	11-047-03RSS	11-047-03RTI	3	3	1
11-047-05LSS	11-047-05LTI	11-047-05RSS	11-047-05RTI	5	3	1

2.4mm AV-Wiselock "L" Dorsal Distal Radius Plate, (3 Head Holes)

Left		Right				
Stainless Steel	Titanium	Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
11-048-03LSS	11-048-03LTI	11-048-03RSS	11-048-03RTI	3	3	1
11-048-05LSS	11-048-05LTI	11-048-05RSS	11-048-05RTI	5	3	1

2.4mm AV-Wiselock "L" Dorsal Distal Radius Plate, Oblique (3 Head Holes)

Left		Right				
Stainless Steel	Titanium	Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
11-049-03LSS	11-049-03LTI	11-049-03RSS	11-049-03RTI	3	3	1
11-049-05LSS	11-049-05LTI	11-049-05RSS	11-049-05RTI	5	3	1

2.4mm AV-Wiselock "T" Dorsal Distal Radius Plate, (3 Head Holes)

Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
11-050-03SS	11-050-03TI	3	3	1
11-050-05SS	11-050-05TI	5	3	1

2.4mm AV-Wiselock Distal Radius Volar Plate, Extra-articular (5 Head Holes)

Left		Right				
Stainless Steel	Titanium	Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
11-043-03LSS	11-043-03LTI	11-043-03RSS	11-043-03RTI	3	5	1
11-043-05LSS	11-043-05LTI	11-043-05RSS	11-043-05RTI	5	5	1

2.4mm AV-Wiselock Distal Radius Volar Plate, Extra-articular (4 Head Holes)

Left		Right				
Stainless Steel	Titanium	Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
11-044-03LSS	11-044-03LTI	11-044-03RSS	11-044-03RTI	3	4	1
11-044-05LSS	11-044-05LTI	11-044-05RSS	11-044-05RTI	5	4	1

2.4mm Variable Angle Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
549.006	TI-549.006	6mm	6
549.008	TI-549.008	8mm	6
549.010	TI-549.010	10mm	6
549.012	TI-549.012	12mm	6
549.014	TI-549.014	14mm	6
549.016	TI-549.016	16mm	6
549.018	TI-549.018	18mm	6
549.020	TI-549.020	20mm	6
549.022	TI-549.022	22mm	6
549.024	TI-549.024	24mm	6
549.026	TI-549.026	26mm	6
549.028	TI-549.028	28mm	6
549.030	TI-549.030	30mm	6

2.4mm Wise-Lock Screw, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1145.06	TI-1145.06	6mm	6
1145.08	TI-1145.08	8mm	6
1145.10	TI-1145.10	10mm	6
1145.12	TI-1145.12	12mm	6
1145.14	TI-1145.14	14mm	6
1145.16	TI-1145.16	16mm	6
1145.18	TI-1145.18	18mm	6
1145.20	TI-1145.20	20mm	6
1145.22	TI-1145.22	22mm	6
1145.24	TI-1145.24	24mm	6
1145.26	TI-1145.26	26mm	6
1145.28	TI-1145.28	28mm	6
1145.30	TI-1145.30	30mm	6

2.4mm Cortical Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1144.06	TI-1144.06	6mm	3
1144.08	TI-1144.08	8mm	3
1144.10	TI-1144.10	10mm	3
1144.12	TI-1144.12	12mm	3
1144.14	TI-1144.14	14mm	3
1144.16	TI-1144.16	16mm	3
1144.18	TI-1144.18	18mm	3
1144.20	TI-1144.20	20mm	3
1144.22	TI-1144.22	22mm	3
1144.24	TI-1144.24	24mm	3
1144.26	TI-1144.26	26mm	3
1144.28	TI-1144.28	28mm	3
1144.30	TI-1144.30	30mm	3

2.7mm Cortical Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1242.06	TI-1242.06	6mm	3
1242.08	TI-1242.08	8mm	3
1242.10	TI-1242.10	10mm	3
1242.12	TI-1242.12	12mm	3
1242.14	TI-1242.14	14mm	3
1242.16	TI-1242.16	16mm	3
1242.18	TI-1242.18	18mm	3
1242.20	TI-1242.20	20mm	3
1242.22	TI-1242.22	22mm	3
1242.24	TI-1242.24	24mm	3
1242.26	TI-1242.26	26mm	3
1242.28	TI-1242.28	28mm	3
1242.30	TI-1242.30	30mm	3

Implant Box with cover for 2.4mm AV-Wiselock Distal Radius Dorsal and Extra-articular System

Code	Units
11-051	1

**1312-33 Implant Box with cover for
2.4mm Variable Angle Distal Radius System**



2.4mm Variable Angle Two-Column Volar Distal Radius Plate, Narrow

Left		Right		Shaft Hole	Head Hole	Width (mm)	Units
Stainless Steel	Titanium	Stainless Steel	Titanium				
546.002L	TI-546.002L	546.202R	TI-546.002R	2	6	19.5	1
546.003L	TI-546.003L	546.203R	TI-546.003R	3	6	19.5	1
546.004L	TI-546.004L	546.204R	TI-546.004R	4	6	19.5	1
546.005L	TI-546.005L	546.205R	TI-546.005R	5	6	19.5	1
546.006L	TI-546.006L	546.206R	TI-546.006R	6	6	19.5	1
546.007L	TI-546.007L	546.207R	TI-546.007R	7	6	19.5	1
546.008L	TI-546.008L	546.208R	TI-546.008R	8	6	19.5	1

Left		Right		Shaft Hole	Head Hole	Width (mm)	Units
Stainless Steel	Titanium	Stainless Steel	Titanium				
546.202L	TI-546.202L	546.202R	TI-546.202R	2	6	22	1
546.203L	TI-546.203L	546.203R	TI-546.203R	3	6	22	1
546.204L	TI-546.204L	546.204R	TI-546.204R	4	6	22	1
546.205L	TI-546.205L	546.205R	TI-546.205R	5	6	22	1
546.206L	TI-546.206L	546.206R	TI-546.206R	6	6	22	1
546.207L	TI-546.207L	546.207R	TI-546.207R	7	6	22	1
546.208L	TI-546.208L	546.208R	TI-546.208R	8	6	22	1

2.4mm Variable Angle Two-Column Volar Distal Radius Plate, Wide

Left		Right		Shaft Hole	Head Hole	Width (mm)	Units
Stainless Steel	Titanium	Stainless Steel	Titanium				
548.002L	TI-548.002L	548.002R	TI-548.002R	2	7	25.5	1
548.003L	TI-548.003L	548.003R	TI-548.003R	3	7	25.5	1
548.004L	TI-548.004L	548.004R	TI-548.004R	4	7	25.5	1
548.005L	TI-548.005L	548.005R	TI-548.005R	5	7	25.5	1
548.006L	TI-548.006L	548.006R	TI-548.006R	6	7	25.5	1
548.007L	TI-548.007L	548.007R	TI-548.007R	7	7	25.5	1
548.008L	TI-548.008L	548.008R	TI-548.008R	8	7	25.5	1

Aiming Block For Variable Angle Two-Column Volar Distal Radius Plate, Narrow 6 Head Holes, 19.5mm

Left				
Stainless Steel	Titanium	Head Hole	Width	Units
650-19.5L	TI-650-19.5L	6	19.5	2

Right				
Stainless Steel	Titanium	Head Hole	Width	Units
650-19.5R	TI-650-19.5R	6	19.5	2

Left				
Stainless Steel	Titanium	Head Hole	Width	Units
650-22L	TI-650-22L	6	22	2

Right				
Stainless Steel	Titanium	Head Hole	Width	Units
650-22R	TI-650-22R	6	22	2

Aiming Block For Variable Angle Two-Column Volar Distal Radius Plate, Wide 7 Head Holes, 25.5mm

Left				
Stainless Steel	Titanium	Head Hole	Width	Units
650-25.5L	TI-650-25.5L	7	25.5	2

Right				
Stainless Steel	Titanium	Head Hole	Width	Units
650-25.5R	TI-650-25.5R	7	25.5	2

Aiming Block Drill Sleeve

Stainless Steel	Titanium	Units
650-01	TI-650-01	2

Aiming Block Screw

Stainless Steel	Titanium	Units
650-02	TI-650-02	5

2.4mm Variable Angle Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
549.006	TI-549.006	6mm	6
549.008	TI-549.008	8mm	6
549.010	TI-549.010	10mm	6
549.012	TI-549.012	12mm	6
549.014	TI-549.014	14mm	6
549.016	TI-549.016	16mm	6
549.018	TI-549.018	18mm	6
549.020	TI-549.020	20mm	6
549.022	TI-549.022	22mm	6
549.024	TI-549.024	24mm	6
549.026	TI-549.026	26mm	6
549.028	TI-549.028	28mm	6
549.030	TI-549.030	30mm	6

2.4mm Wise-Lock Screw, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1145.06	TI-1145.06	6mm	6
1145.08	TI-1145.08	8mm	6
1145.10	TI-1145.10	10mm	6
1145.12	TI-1145.12	12mm	6
1145.14	TI-1145.14	14mm	6
1145.16	TI-1145.16	16mm	6
1145.18	TI-1145.18	18mm	6
1145.20	TI-1145.20	20mm	6
1145.22	TI-1145.22	22mm	6
1145.24	TI-1145.24	24mm	6
1145.26	TI-1145.26	26mm	6
1145.28	TI-1145.28	28mm	6
1145.30	TI-1145.30	30mm	6

2.4mm Cortical Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1144.06	TI-1144.06	6mm	3
1144.08	TI-1144.08	8mm	3
1144.10	TI-1144.10	10mm	3
1144.12	TI-1144.12	12mm	3
1144.14	TI-1144.14	14mm	3
1144.16	TI-1144.16	16mm	3
1144.18	TI-1144.18	18mm	3
1144.20	TI-1144.20	20mm	3
1144.22	TI-1144.22	22mm	3
1144.24	TI-1144.24	24mm	3
1144.26	TI-1144.26	26mm	3
1144.28	TI-1144.28	28mm	3
1144.30	TI-1144.30	30mm	3

2.7mm Cortical Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1242.06	TI-1242.06	6mm	3
1242.08	TI-1242.08	8mm	3
1242.10	TI-1242.10	10mm	3
1242.12	TI-1242.12	12mm	3
1242.14	TI-1242.14	14mm	3
1242.16	TI-1242.16	16mm	3
1242.18	TI-1242.18	18mm	3
1242.20	TI-1242.20	20mm	3
1242.22	TI-1242.22	22mm	3
1242.24	TI-1242.24	24mm	3
1242.26	TI-1242.26	26mm	3
1242.28	TI-1242.28	28mm	3
1242.30	TI-1242.30	30mm	3

10-050 : Implant Box with cover for 2.4mm Variable Angle Volar Rim Distal Radius System



2.4mm Variable Angle Volar Rim Distal Radius Plate

Left		Right				
Stainless Steel	Titanium	Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
560-06-03L	TI-560-06-03L	560-06-03R	TI-560-06-03R	3	6	1
560-06-04L	TI-560-06-04L	560-06-04R	TI-560-06-04R	4	6	1
560-06-05L	TI-560-06-05L	560-06-05R	TI-560-06-05R	5	6	1
560-06-06L	TI-560-06-06L	560-06-06R	TI-560-06-06R	6	6	1
560-06-07L	TI-560-06-07L	560-06-07R	TI-560-06-07R	7	6	1
560-06-08L	TI-560-06-08L	560-06-08R	TI-560-06-08R	8	6	1

Left		Right				
Stainless Steel	Titanium	Stainless Steel	Titanium	Shaft Hole	Head Hole	Units
560-07-03L	TI-560-07-03L	560-07-03R	TI-560-07-03R	3	7	1
560-07-04L	TI-560-07-04L	560-07-04R	TI-560-07-04R	4	7	1
560-07-05L	TI-560-07-05L	560-07-05R	TI-560-07-05R	5	7	1
560-07-06L	TI-560-07-06L	560-07-06R	TI-560-07-06R	6	7	1
560-07-07L	TI-560-07-07L	560-07-07R	TI-560-07-07R	7	7	1
560-07-08L	TI-560-07-08L	560-07-08R	TI-560-07-08R	8	7	1

Aiming Block For Variable Angle Volar Rim Distal Radius Plate, 6 Head Hole

Left			
Stainless Steel	Titanium	Head Hole	Units
651-06L	TI-651-06L	6	2

Right			
Stainless Steel	Titanium	Head Hole	Units
651-06R	TI-651-06R	6	2

Aiming Block For Variable Angle Volar Rim Distal Radius Plate, 7 Head Hole

Left			
Stainless Steel	Titanium	Head Hole	Units
651-07L	TI-651-07L	7	2

Right			
Stainless Steel	Titanium	Head Hole	Units
651-07R	TI-651-07R	7	2

Aiming Block Drill Sleeve

Stainless Steel	Titanium	Units
650-01	TI-650-01	2

Aiming Block Screw

Stainless Steel	Titanium	Units
650-02	TI-650-02	5

2.4mm Variable Angle Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
549.006	TI-549.006	6mm	6
549.008	TI-549.008	8mm	6
549.010	TI-549.010	10mm	6
549.012	TI-549.012	12mm	6
549.014	TI-549.014	14mm	6
549.016	TI-549.016	16mm	6
549.018	TI-549.018	18mm	6
549.020	TI-549.020	20mm	6
549.022	TI-549.022	22mm	6
549.024	TI-549.024	24mm	6
549.026	TI-549.026	26mm	6
549.028	TI-549.028	28mm	6
549.030	TI-549.030	30mm	6

2.4mm Wise-Lock Screw, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1145.06	TI-1145.06	6mm	6
1145.08	TI-1145.08	8mm	6
1145.10	TI-1145.10	10mm	6
1145.12	TI-1145.12	12mm	6
1145.14	TI-1145.14	14mm	6
1145.16	TI-1145.16	16mm	6
1145.18	TI-1145.18	18mm	6
1145.20	TI-1145.20	20mm	6
1145.22	TI-1145.22	22mm	6
1145.24	TI-1145.24	24mm	6
1145.26	TI-1145.26	26mm	6
1145.28	TI-1145.28	28mm	6
1145.30	TI-1145.30	30mm	6

2.4mm Cortical Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1144.06	TI-1144.06	6mm	3
1144.08	TI-1144.08	8mm	3
1144.10	TI-1144.10	10mm	3
1144.12	TI-1144.12	12mm	3
1144.14	TI-1144.14	14mm	3
1144.16	TI-1144.16	16mm	3
1144.18	TI-1144.18	18mm	3
1144.20	TI-1144.20	20mm	3
1144.22	TI-1144.22	22mm	3
1144.24	TI-1144.24	24mm	3
1144.26	TI-1144.26	26mm	3
1144.28	TI-1144.28	28mm	3
1144.30	TI-1144.30	30mm	3

2.7mm Cortical Screws, Self-Tapping, (Star Head)

Stainless Steel	Titanium	Length	Units
1242.06	TI-1242.06	6mm	3
1242.08	TI-1242.08	8mm	3
1242.10	TI-1242.10	10mm	3
1242.12	TI-1242.12	12mm	3
1242.14	TI-1242.14	14mm	3
1242.16	TI-1242.16	16mm	3
1242.18	TI-1242.18	18mm	3
1242.20	TI-1242.20	20mm	3
1242.22	TI-1242.22	22mm	3
1242.24	TI-1242.24	24mm	3
1242.26	TI-1242.26	26mm	3
1242.28	TI-1242.28	28mm	3
1242.30	TI-1242.30	30mm	3

Please choose any one implant box as per requirement

11-051 Implant Box with cover for 2.4mm AV-Wiselock Distal Radius Dorsal and Extra-articular System



10-050 Implant Box with cover for 2.4mm Variable Angle Volar Rim Distal Radius System



11-077 Implant Box with cover for 2.4mm AV-Wiselock Hand and Wrist - Arthrodesis System



1312-33 Implant Box with cover for 2.4mm Variable Angle Distal Radius system



1312-33	Implant Box with cover for 2.4mm Variable Angle Distal Radius system	1
10-050	Implant Box with cover for 2.4mm Variable Angle Volar Rim Distal Radius System	1
11-051	Implant Box with cover for 2.4mm AV-Wiselock Distal Radius Dorsal & Extra-Articular System	1
11-077	Implant Box with cover for 2.4mm AV-Wiselock Hand and Wrist - Arthrodesis System	1

2.4mm

Instrumentation

1312-01 Depth Gauge measuring upto 60mm




1312-02 T-Handle With Quick Coupling



1312-03-15 Retractor, 15mm



1312-04-1.8 Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm



1312-05-2.0 Drill Bit with Quick Coupling End, Ø2.0mm x Length 140mm




1312-06-2.4 Drill Bit with Quick Coupling End, Ø2.4mm x Length 140mm



1312-07-2.7 Drill Bit with Quick Coupling End, Ø2.7mm x Length 140mm



1312-08-1.2 Guide Wire, Ø1.2mm x Length 150mm



1312-09-1.2 Guide Wire With Threaded Tip, Ø1.2mm x Thread Length 5mm x Length 150mm




1312-10-2.4 Tap, HA 2.4mm



1312-11-2.7 Tap, HA 2.7mm



1312-13-1.8 Threaded Drill Sleeve, Ø1.8/2.4mm (Double Start Thread)



1312-15 Self-Centering Double Drill Guide Φ 1.8/2.4mm




1312-16 Self-Centering Double Drill Guide Φ 2.0/2.7mm



1312-17 Depth Gauge with Protector measuring upto 50mm



1312-18 Star Screwdriver Shaft, T8



1312-21 Torque Limiting Attachment, 0.8Nm



7-019-03 Handle for Torque Limiting Attachment



1312-22-6.5

Retractor, 6.5mm



1312-23

Double Periosteal Elevator, Round



1312-24

Double Periosteal Elevator, Flat



1312-25

Reduction Forcep with Points, Small



1312-26

Reduction Forcep with Points, Large



1312-27

Reduction Forcep with Points, Serrated



1312-28

Bending Forcep

**1312-30**

Universal Variable Angle Drill Guide, Ø1.8mm

**1312-20-2.7**

Screw Holding Sleeve for 2.7mm Cortical Screw

**1312-20-2.4**

Screw Holding Sleeve for 2.4mm Cortical Screw

**1312-30-2.4**

Screw Holding Sleeve for 2.4mm Wise-lock Screw

**1312-22**

Quick Coupling Handle



1312-29

Sharp Hook




1312-31

Screwdriver Star Head




1312-35

Instrument trays for 2.4mm AV-Wiselock Distal Radius Instrument Set






1312-34

Container for 2.4mm AV-Wiselock Distal Radius Instrument



1312-000 2.4mm AV-Wiselock Distal Radius Instrument Set



Code	Set Consisting of	Units
1312-01	Depth Gauge measuring upto 60mm	1
1312-02	T-Handle With Quick Coupling	1
1312-03-15	Retractor, 15mm	1
1312-04-1.8	Drill Bit with Quick Coupling End, Ø1.8mm x Length 140mm	3
1312-05-2.0	Drill Bit with Quick Coupling End, Ø2.0mm x Length 140mm	2
1312-06-2.4	Drill Bit with Quick Coupling End, Ø2.4mm x Length 140mm	1
1312-07-2.7	Drill Bit with Quick Coupling End, Ø2.7mm x Length 140mm	1
1312-08-1.2	Guide Wire, Ø1.2mm x Length 150mm	5
1312-09-1.2	Guide Wire With Threaded Tip, Ø1.2mm x Thread Length 5mm x Length 150mm	2
1312-10-2.4	Tap, HA 2.4mm	1
1312-11-2.7	Tap, HA 2.7mm	1
1312-13-1.8	Threaded Drill Sleeve, Ø1.8/2.4mm (Double Start Thread)	3
1312-15	Self-Centering Double Drill Guide Ø1.8/2.4mm	1
1312-16	Self-Centering Double Drill Guide Ø2.0/2.7mm	1
1312-17	Depth Gauge with Protector measuring upto 50mm	1
1312-18	Star Screwdriver Shaft, T8	2
1312-21	Torque Limiting Attachment, 0.8Nm	1
7-019-03	Handle for Torque Limiting Attachment	1
1312-22-6.5	Retractor, 6.5mm	1
1312-23	Double Periosteal Elevator, Round	1
1312-24	Double Periosteal Elevator, Flat	1
1312-25	Reduction Forcep with Points, Small	1
1312-26	Reduction Forcep with Points, Large	1
1312-27	Reduction Forcep with Points, Serrated	1
1312-28	Bending Forcep	2
1312-30	Universal Variable Angle Drill Guide, Ø1.8mm	1
1312-20-2.7	Screw Holding Sleeve for 2.7mm Cortical Screw	1
1312-20-2.4	Screw Holding Sleeve for 2.4mm Cortical Screw	1
1312-30-2.4	Screw Holding Sleeve for 2.4mm Wise-lock Screw	1
1312-22	Quick Coupling Handle	1
1312-29	Sharp Hook	1
1312-31	Screwdriver Star Head	1
1312-35	Instrument trays for 2.4mm AV-Wiselock Distal Radius Instrument Set	2
1312-34	Container for 2.4mm AV-Wiselock Distal Radius Instrument	1

Please choose any one implant box as per requirement

1312-33	Implant Box with cover for 2.4mm Variable Angle Distal Radius system	1
10-050	Implant Box with cover for 2.4mm Variable Angle Volar Rim Distal Radius System	1
11-051	Implant Box with cover for 2.4mm AV-Wiselock Distal Radius Dorsal & Extra-Articular System	1
11-077	Implant Box with cover for 2.4mm AV-Wiselock Hand and Wrist - Arthrodesis System	1