

TARSUS SCREW SYSTEM

TARSUS SCREW SYSTEM



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.


Our Achievements



The Tarsus Screw is an implantable stabilization device designed to prevent hyperpronation (Talotarsal Displacement), an excessive inward rotation of the foot that occurs when the talus and the calcaneus slip out of alignment. It is consisting of the implants of various sizes and an instrumentation set.

The implant is manufactured in Titanium alloy as per ISO 5832-3 Standard.

Features:



Improved Screw thread design for better resistance to pull out force

New design Screw head that lock with Driver through Positioning Sleeve

Indications

Tarsus Screw is a stabilization device used in the treatment of hyperpronating instability of the hindfoot. The implant is designed to stabilize the talus to prevent excessive anterior, and/or medial, and/or plantarflexion of the talus, while allowing normal talotarsal joint motion.

Tarsus Screw


Dia	Titanium
6	1-048-06TI
7	1-048-07TI
8	1-048-08TI
9	1-048-09TI
10	1-048-10TI
11	1-048-11TI
12	1-048-12TI


STERILE

Dia	Titanium
6	1-048-06TI-S
7	1-048-07TI-S
8	1-048-08TI-S
9	1-048-09TI-S
10	1-048-10TI-S
11	1-048-11TI-S
12	1-048-12TI-S

7-149-01

Tarsus Universal Driver


7-149-02

Positioning Sleeve



Tarsus Trial Implant

Dia	Code
6	7-149-03
7	7-149-04
8	7-149-05
9	7-149-06
10	7-149-07
11	7-149-08
12	7-149-09


7-149-10

Curved Tenotomy Scissors



7-149-12

Guide Wire with Blunt Tip, Ø1.9mm X 250mm




7-149-11

Container for Tarsus Screw and Instrument



7-149-13

Lid for Tarsus Screw Implant and Instrument Set



7-149 : Tarsus Screw Instrument and Implant Set



7-149 : Tarsus Screw Instrument and Implant Set

Code	Set Consisting of	Qty.
7-149-01	Tarsus Universal Driver	1
7-149-02	Positioning Sleeve	1
7-149-03	Tarsus Trial Implant, Ø6mm	1
7-149-04	Tarsus Trial Implant, Ø7mm	1
7-149-05	Tarsus Trial Implant, Ø8mm	1
7-149-06	Tarsus Trial Implant, Ø9mm	1
7-149-07	Tarsus Trial Implant, Ø10mm	1
7-149-08	Tarsus Trial Implant, Ø11mm	1
7-149-09	Tarsus Trial Implant, Ø12mm	1
7-149-10	Curved Tenotomy Scissors	1
7-149-12	Guide Wire with Blunt Tip, Ø1.9mm X 250mm	3
7-149-11	Container for Tarsus Screw and Instrument	1
7-149-13	Lid for Tarsus Screw Implant and Instrument Set	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com