


RADIAL HEAD SYSTEM



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements




Radial Head System

The Radial Head implant is a specialized medical device designed to replace a damaged or fractured radial head.


The Radial Head System is radial head replacement implant which is designed to facilitate replacement of the patient's radial head. For the specific radial head fractures, the surgeon can opt amongst various diameters and heights which suits the patient anatomy in efficient way.

- Contoured lateral surface enhances interface with the annular ligament.
- Medial ulnar articular zone facilitates the optimum contact with radial notch.



Various Sizes with multiple Combination
of Radial Head and Stem

Implants Available in
Cobalt Chrome & Cobalt Chrome(with HA Coating)



Radial Head , Cobalt Chrome, Sterile


Dia	Height	Code
Ø16mm	Standard	8-013-16-01CC-S
Ø18mm	Standard	8-013-18-01CC-S
Ø20mm	Standard	8-013-20-01CC-S
Ø22mm	Standard	8-013-22-01CC-S
Ø24mm	Standard	8-013-24-01CC-S
Ø16mm	+2	8-013-16-02CC-S
Ø18mm	+2	8-013-18-02CC-S
Ø20mm	+2	8-013-20-02CC-S
Ø22mm	+2	8-013-22-02CC-S
Ø24mm	+2	8-013-24-02CC-S
Ø16mm	+4	8-013-16-04CC-S
Ø18mm	+4	8-013-18-04CC-S
Ø20mm	+4	8-013-20-04CC-S
Ø22mm	+4	8-013-22-04CC-S
Ø24mm	+4	8-013-24-04CC-S

Radial Stem Standard, Cobalt Chrome, Sterile



Dia	Length	Code
Ø4.5mm	16mm	8-014-4.5-16CC-S
Ø5.5mm	16mm	8-014-5.5-16CC-S
Ø6.5mm	16mm	8-014-6.5-16CC-S
Ø7.5mm	17mm	8-014-7.5-17CC-S
Ø8.5mm	17mm	8-014-8.5-17CC-S
Ø9.5mm	17mm	8-014-9.5-17CC-S

Radial Stem Long, Cobalt Chrome, Sterile



Dia	Length	Code
Ø4.5mm	24mm	8-015-4.5-24CC-S
Ø5.5mm	24mm	8-015-5.5-24CC-S
Ø6.5mm	24mm	8-015-6.5-24CC-S
Ø7.5mm	25mm	8-015-7.5-25CC-S
Ø8.5mm	25mm	8-015-8.5-25CC-S
Ø9.5mm	25mm	8-015-9.5-25CC-S

Radial Stem Standard, Cobalt Chrome (with HA Coating), Sterile



Dia	Length	Code
Ø4.5mm	16mm	8-014-4.5-16CC-HA-S
Ø5.5mm	16mm	8-014-5.5-16CC-HA-S
Ø6.5mm	16mm	8-014-6.5-16CC-HA-S
Ø7.5mm	17mm	8-014-7.5-17CC-HA-S
Ø8.5mm	17mm	8-014-8.5-17CC-HA-S
Ø9.5mm	17mm	8-014-9.5-17CC-HA-S

Radial Stem Long, Cobalt Chrome (with HA Coating), Sterile






Dia	Length	Code
Ø4.5mm	24mm	8-015-4.5-24CC-HA-S
Ø5.5mm	24mm	8-015-5.5-24CC-HA-S
Ø6.5mm	24mm	8-015-6.5-24CC-HA-S
Ø7.5mm	25mm	8-015-7.5-25CC-HA-S
Ø8.5mm	25mm	8-015-8.5-25CC-HA-S
Ø9.5mm	25mm	8-015-9.5-25CC-HA-S

7-154 : Instrument Set for Radial Head System, Reusable
7-154-01 Awl for Radial Head


Sequential Reamer

Dia	Code
Ø4.5mm	7-154-02
Ø5.5mm	7-154-03
Ø6.5mm	7-154-04
Ø7.5mm	7-154-05
Ø8.5mm	7-154-06
Ø9.5mm	7-154-07


7-154-08 Planer for Radial Head

7-154-09 Inserter cap for Planer

7-154-10 Impactor for Radial Head


7-154-11


Head Sizer


Radial Head Trial

Dia	Height	Code
Ø16mm	Standard	7-154-12-01
Ø18mm	Standard	7-154-12-02
Ø20mm	Standard	7-154-12-03
Ø22mm	Standard	7-154-12-04
Ø24mm	Standard	7-154-12-05
Ø16mm	+2	7-154-12-06
Ø18mm	+2	7-154-12-07
Ø20mm	+2	7-154-12-08
Ø22mm	+2	7-154-12-09
Ø24mm	+2	7-154-12-10
Ø16mm	+4	7-154-12-11
Ø18mm	+4	7-154-12-12
Ø20mm	+4	7-154-12-13
Ø22mm	+4	7-154-12-14
Ø24mm	+4	7-154-12-15


Radial Stem Trial Standard

Dia	Length	Code
Ø4.5mm	16mm	7-154-13-01
Ø5.5mm	16mm	7-154-13-02
Ø6.5mm	16mm	7-154-13-03
Ø7.5mm	17mm	7-154-13-04
Ø8.5mm	17mm	7-154-13-05
Ø9.5mm	17mm	7-154-13-06



7-154 : Instrument Set for Radial Head System, Reusable

Code	Description	Qty.
7-154-01	Awl for Radial Head	1
7-154-02	Sequential Reamer, Ø4.5mm	1
7-154-03	Sequential Reamer, Ø5.5mm	1
7-154-04	Sequential Reamer, Ø6.5mm	1
7-154-05	Sequential Reamer, Ø7.5mm	1
7-154-06	Sequential Reamer, Ø8.5mm	1
7-154-07	Sequential Reamer, Ø9.5mm	1
7-154-08	Planer for Radial Head	1
7-154-09	Insert cap for Planer	1
7-154-10	Impactor for Radial Head	1
7-154-11	Head Sizer	1
7-154-12-01	Radial Head Trial, Ø16mm, Standard Height	1
7-154-12-02	Radial Head Trial, Ø18mm, Standard Height	1
7-154-12-03	Radial Head Trial, Ø20mm, Standard Height	1
7-154-12-04	Radial Head Trial, Ø22mm, Standard Height	1
7-154-12-05	Radial Head Trial, Ø24mm, Standard Height	1
7-154-12-06	Radial Head Trial, Ø16mm, (+2) Height	1
7-154-12-07	Radial Head Trial, Ø18mm, (+2) Height	1
7-154-12-08	Radial Head Trial, Ø20mm, (+2) Height	1
7-154-12-09	Radial Head Trial, Ø22mm, (+2) Height	1
7-154-12-10	Radial Head Trial, Ø24mm, (+2) Height	1
7-154-12-11	Radial Head Trial, Ø16mm, (+4) Height	1
7-154-12-12	Radial Head Trial, Ø18mm, (+4) Height	1
7-154-12-13	Radial Head Trial, Ø20mm, (+4) Height	1
7-154-12-14	Radial Head Trial, Ø22mm, (+4) Height	1
7-154-12-15	Radial Head Trial, Ø24mm, (+4) Height	1
7-154-13-01	Radial Stem Trial Standard, Ø4.5mm x Length 16mm	1
7-154-13-02	Radial Stem Trial Standard, Ø5.5mm x Length 16mm	1
7-154-13-03	Radial Stem Trial Standard, Ø6.5mm x Length 16mm	1
7-154-13-04	Radial Stem Trial Standard, Ø7.5mm x Length 17mm	1
7-154-13-05	Radial Stem Trial Standard, Ø8.5mm x Length 17mm	1
7-154-13-06	Radial Stem Trial Standard, Ø9.5mm x Length 17mm	1

7-154-S : Instrument Set for Radial Head System, Sterile, Single Use

Code	Description	Qty.
7-154-01-S	Awl for Radial Head, Sterile	1
7-154-02-S	Sequential Reamer, Ø4.5mm, Sterile	1
7-154-03-S	Sequential Reamer, Ø5.5mm, Sterile	1
7-154-04-S	Sequential Reamer, Ø6.5mm, Sterile	1
7-154-05-S	Sequential Reamer, Ø7.5mm, Sterile	1
7-154-06-S	Sequential Reamer, Ø8.5mm, Sterile	1
7-154-07-S	Sequential Reamer, Ø9.5mm, Sterile	1
7-154-08-S	Planer for Radial Head, Sterile	1
7-154-09-S	Insert cap for Planer, Sterile	1
7-154-10-S	Impactor for Radial Head, Sterile	1
7-154-11-S	Head Sizer, Sterile	1
7-154-12-01-S	Radial Head Trial, Ø16mm, Standard Height, Sterile	1
7-154-12-02-S	Radial Head Trial, Ø18mm, Standard Height, Sterile	1
7-154-12-03-S	Radial Head Trial, Ø20mm, Standard Height, Sterile	1
7-154-12-04-S	Radial Head Trial, Ø22mm, Standard Height, Sterile	1
7-154-12-05-S	Radial Head Trial, Ø24mm, Standard Height, Sterile	1
7-154-12-06-S	Radial Head Trial, Ø16mm, (+2) Height, Sterile	1
7-154-12-07-S	Radial Head Trial, Ø18mm, (+2) Height, Sterile	1
7-154-12-08-S	Radial Head Trial, Ø20mm, (+2) Height, Sterile	1
7-154-12-09-S	Radial Head Trial, Ø22mm, (+2) Height, Sterile	1
7-154-12-10-S	Radial Head Trial, Ø24mm, (+2) Height, Sterile	1
7-154-12-11-S	Radial Head Trial, Ø16mm, (+4) Height, Sterile	1
7-154-12-12-S	Radial Head Trial, Ø18mm, (+4) Height, Sterile	1
7-154-12-13-S	Radial Head Trial, Ø20mm, (+4) Height, Sterile	1
7-154-12-14-S	Radial Head Trial, Ø22mm, (+4) Height, Sterile	1
7-154-12-15-S	Radial Head Trial, Ø24mm, (+4) Height, Sterile	1
7-154-13-01-S	Radial Stem Trial Standard, Ø4.5mm x Length 16mm, Sterile	1
7-154-13-02-S	Radial Stem Trial Standard, Ø5.5mm x Length 16mm, Sterile	1
7-154-13-03-S	Radial Stem Trial Standard, Ø6.5mm x Length 16mm, Sterile	1
7-154-13-04-S	Radial Stem Trial Standard, Ø7.5mm x Length 17mm, Sterile	1
7-154-13-05-S	Radial Stem Trial Standard, Ø8.5mm x Length 17mm, Sterile	1
7-154-13-06-S	Radial Stem Trial Standard, Ø9.5mm x Length 17mm, Sterile	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com