


auxilock® 



Biobraid

AUXILOCK[®] BIOBRAID SUTURES



AUXILOCK[®] BioBraid Sutures are braided sterile sutures prepared from Ultra High Molecular Weight Polyethylene (UHMWPE). They can be used in soft tissue approximations along with allograft tissues in arthroscopy procedures. BioBraid sutures are non-absorbable and do not impose any significant changes in tensile strength retention known to occur in vivo. BioBraid sutures are available in a variety of colour combinations for easy suture management in complicated repairs.

Ordering Information

Code	Product Description	Suture Colour	Needle
6-011-03	AUXILOCK [®] Two #2 BioBraid: White/Blue & White/Black, 36in Total Length		
6-011-04	AUXILOCK [®] #2 BioBraid: Blue, With Needles: MO-6 & CE, 36in Total Length		
6-011-05	AUXILOCK [®] #2 BioBraid: White, With Needle: MO-6, 36in Total Length		
6-011-06	AUXILOCK [®] #2 BioBraid: White/Black, With Needle: MO-6, 36in Total Length		
6-011-08	AUXILOCK [®] #2 BioBraid: White/Blue, With Needle: MO-6, 36in Total Length		
6-011-09	AUXILOCK [®] #2 BioBraid: White/Black, 36in Total Length		
6-011-10	AUXILOCK [®] #2 BioBraid: White/Green, With Needle: MO-6, 36in Total Length		
6-011-11	AUXILOCK [®] #2 BioBraid: White/Green, 36in Total Length		
6-011-12	AUXILOCK [®] #2 BioBraid: Blue, 36in Total Length		
6-011-13	AUXILOCK [®] #2 BioBraid: White, 36in Total Length		
6-011-17	AUXILOCK [®] #2 BioBraid: White/Blue, 36in Total Length		
6-011-32	AUXILOCK [®] #2 BioBraid: White, With Needle: CCS, 36in Total Length		
6-011-07	AUXILOCK [®] #5 BioBraid: White, With Needle: CCS, 36in Total Length		
6-011-33	AUXILOCK [®] #2-0 BioBraid: White/Blue, 36in Total Length		

Note: Each box consist of 10 Pieces


MO-6 : - 1/2 Circle Taper Point Long Taper Ratio

CE : - 1/2 Circle Reverse Cutting

CCS : - 1/2 Circle Trocar point







STERILE EO

AUXILOCK[®]
BIOBRAID INFINITY
LOOP NEEDLE



AUXILOCK[®] BioBraid Infinity Loop is a continuous loop of BioBraid with a straight needle. BioBraid Infinity Loop utilises the whipstitch technique for graft preparation and reduces the time spent. It also compresses the graft uniformly and improves strength.


Ordering Information

Code	Product Description	Suture Colour	Needle
6-012-01	AUXILOCK [®] #2 BioBraid Infinity Loop: White, With Needle: Straight, 20in Loop Length, 24in Total Length		
6-012-02	AUXILOCK [®] #2 BioBraid Infinity Loop: White/Blue, With Needle: Straight, 20in Loop Length, 24in Total Length		
6-012-03	AUXILOCK [®] #5 BioBraid Infinity Loop: White/Blue, With Needle: Straight, 20in Loop Length, 24in Total Length		

Note: Each box consist of 10 Pieces








STERILE 

AUXILOCK[®] BIOBRAID SUTURE TAPE



AUXILOCK[®] BioBraid Suture Tape made from UHMWPE. The suture tape design accommodates the tape in the centre with a suture on the either sides of the tape. It can be used in Knee, Shoulder, and Foot & Ankle surgeries including AC joint reconstruction, knotless double-row rotator cuff repair, patella fracture repair, MPFL Reconstruction etc. The broad footprint of the BioBraid Suture Tape is appropriate for repairs in degenerative cuff tissue where tissue pull-through may be a concern.

Ordering Information

Code	Product Description	Needle
6-013-01	AUXILOCK [®] 1.8mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-02	AUXILOCK [®] 1.4mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-03	AUXILOCK [®] 1.4mm BioBraid Suture Tape: White/Black, 39in Total length	
6-013-04	AUXILOCK [®] 1.6mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-05	AUXILOCK [®] 2.0mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-06	AUXILOCK [®] 1.0mm BioBraid Suture Tape: White/Blue, 39in Total length	
6-013-07	AUXILOCK [®] 1.4mm BioBraid Suture Tape: Solid Black, 39in Total length	

Note: Each box consist of 10 Pieces

STERILE EO

AUXILOCK PROTEAN[®] CANNULA

The silicon cannula enhances both the inside and outside of the arthroscopy workspace visibility and movement. It establishes a stable portal that eliminates the risk of cannula loss during the insertion and removal of equipment. They are simply inserted into the incision, and numerous surgeries involving the shoulder, knee, hip, and elbow can be completed with the silicon Cannula.

- The double dam has been intricately designed to successfully control fluid.
- Translucent hue for enhanced visibility and stiff flange to mitigate unintentional removal of cannula.
- 5 mm spacer ring for gradual length.



AUXILOCK Protean[®] Cannula

Code	Set Consisting of
6-060-01	AUXILOCK Protean [®] Cannula, Ø6.0 X 20mm
6-060-02	AUXILOCK Protean [®] Cannula, Ø6.0 X 30mm
6-060-03	AUXILOCK Protean [®] Cannula, Ø6.0 X 40mm
6-060-04	AUXILOCK Protean [®] Cannula, Ø6.0 X 50mm
6-060-05	AUXILOCK Protean [®] Cannula, Ø8.0 X 20mm
6-060-06	AUXILOCK Protean [®] Cannula, Ø8.0 X 30mm
6-060-07	AUXILOCK Protean [®] Cannula, Ø8.0 X 40mm
6-060-08	AUXILOCK Protean [®] Cannula, Ø8.0 X 50mm
6-060-09	AUXILOCK Protean [®] Cannula, Ø8.0 X 60mm
6-060-10	AUXILOCK Protean [®] Cannula, Ø10.0 X 20mm
6-060-11	AUXILOCK Protean [®] Cannula, Ø10.0 X 30mm
6-060-12	AUXILOCK Protean [®] Cannula, Ø10.0 X 40mm
6-060-13	AUXILOCK Protean [®] Cannula, Ø10.0 X 50mm





USA

Auxein Inc.

1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172

Tel: +1 305 395 6062

E Fax: +1 305 395 6262

Email: USoffice@auxein.com

www.auxein.com

MEXICO

Auxein México, S.A. de C.V.

Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760

Tel: +521 55 7261 0318

Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.

Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197

Email: info@auxein.com