


**auxilock<sup>®</sup>**   
HIP Arthroscopy System



**HIP**  
Arthroscopy System

**AUXILOCK<sup>®</sup>**  
**KNOTLESS**  
**PEEK CF PUSH-IN**  
**SUTURE ANCHOR**



**Ordering Information**

| Code     | Product Description                                                                                                 |
|----------|---------------------------------------------------------------------------------------------------------------------|
| 6-038-02 | AUXILOCK <sup>®</sup> 2.8mm Knotless PEEK CF Push-In Suture Anchor With Hip Length Driver, With Suture Passer (Hip) |

**Recommended Instrumentation**

| Code     | Product Description                     | Compatible Implants |
|----------|-----------------------------------------|---------------------|
| 7-138-01 | 8 Point Drill Guide, Long               | 6-038-02            |
| 7-138-02 | Ø2.5mm X 360mm Knotless Drill Bit, Long | 6-038-02            |
| 7-138-03 | Trocar Tip Obturator, Long              | 6-038-02            |


STERILE EO

AUXILOCK<sup>®</sup>

KNOTLESS

PEEK OPTIMA PUSH-IN

SUTURE ANCHOR



Ordering Information

| Code     | Product Description                                                                                                     |
|----------|-------------------------------------------------------------------------------------------------------------------------|
| 6-039-02 | AUXILOCK <sup>®</sup> 2.8mm Knotless PEEK OPTIMA Push-In Suture Anchor With Hip Length Driver, With Suture Passer (Hip) |

Recommended Instrumentation

| Code     | Product Description                     | Compatible Implants |
|----------|-----------------------------------------|---------------------|
| 7-138-01 | 8 Point Drill Guide, Long               | 6-039-02            |
| 7-138-02 | Ø2.5mm X 360mm Knotless Drill Bit, Long | 6-039-02            |
| 7-138-03 | Trocar Tip Obturator, Long              | 6-039-02            |

STERILE EO

7-138-01  
8 Point Drill Guide, Long



7-138-02  
Ø2.5mm X 360mm Knotless Drill Bit, Long  
(for 2.8mm Knotless Anchor)




7-138-03  
Trocar Tip Obturator, Long



## Optional Instruments

7-138-04  
Tissue Grasper Long





#### **USA**

##### **Auxein Inc.**

1500 Nw 89th Court, Suite 107-108  
Doral, Florida 33172

Tel: +1 305 395 6062

E Fax: +1 305 395 6262

Email: [USoffice@auxein.com](mailto:USoffice@auxein.com)

#### **MEXICO**

##### **Auxein México, S.A. de C.V.**

Tepic 139 int 801, Colonia Roma Sur,  
Alcaldía Cuauhtémoc, CDMX,  
México, C.P. 06760

Tel: +521 55 7261 0318

Email: [info@auxein.mx](mailto:info@auxein.mx)

#### **INDIA**

##### **Auxein Medical Pvt. Ltd.**

Plot No. 168-169-170, Phase-4,  
Kundli Industrial Area,  
HSIIDC, Sector-57, Sonapat - 131028, Haryana

Tel: +91 99106 43638 | Fax: +91 86077 70197

Email: [info@auxein.com](mailto:info@auxein.com)

[www.auxein.com](http://www.auxein.com)