


CRANIO MAXILLOFACIAL SYSTEM



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements



Orthopaedic Implants

Contents

1.0mm	1 - 8	
1.5mm	9-23	
2.0mm	24-49	
2.4mm	50-60	
Lag screw Kirschner Screw IMF Screw Anchor Screw	61-64	
Instrumentation	65-83	

Neuro | Facial | Mandible | Orthognathic



CMF System a plating platform for internal fixation of the cranio maxillo facial skeleton addressing neuro, craniofacial, mandibular, and orthognathic surgery. It is simple, yet comprehensive system that offers flexibility and ease of use.

- All screws work with all plates within each CMF System.
- Reduced plate/screw profile.
- Standardized instrumentation.

Indications:

Auxein's CMF system is intended for oral, maxillofacial surgery:

- Trauma
- Reconstruction surgery
- Orthognathic Surgery (surgical correction of dentofacial deformities).

Bone Plate & Screw



1.0mm

1.5mm

2.0mm

2.4mm



Neuro
Facial
Mandible
Orthognathic

Screw System	Matched Drill Bit	Pin Emergency Screw	Bone Plate (Thickness)	Indication
1.0mm	0.8mm	1.3mm	0.4mm Mini Plate	Repair and fixation of fracture of cranial base, orbital floor, nasofrontal region and maxillofacial zone of children
1.5mm	1.2mm	1.9mm	0.6mm Mini Plate	Fixation of fracture of zygoma and zygomatic arch, eye pit, maxillary
1.5mm	1.2mm	1.9mm	1.0mm Micro Plate	Correction and fixation of osteotomy of maxillary and mandible
2.0mm	1.7mm	2.3mm	0.6mm Low-incisura Plate	Fracture of Midface, correction and fixation of maxillary
2.0mm	1.7mm	2.3mm	0.8mm Micro Plate	Fracture of maxillary, fixation of orthogathia
2.0mm	1.7mm	2.3mm	1.0mm Micro Plate	Fracture of maxillary and mandible, fixation of orthogathia
2.0mm Locking	1.7mm	2.3mm	1.0mm Locking Plate	Fracture of maxillary and mandible, fixation of orthogathia and osteotomy
2.0mm Locking	1.7mm	2.3mm	1.2mm Locking Plate	Fracture of mandible, fixation of orthogathia and osteotomy
2.0mm	1.7mm	2.3mm	1.5mm One-Way Compression Plate	Reconstruction and fixation of defect of mandible
2.4mm	2.0mm	2.7mm	2.0mm General Reconstruction Plate	Fracture of mandible, fixation of osteotomy
2.4mm Locking	2.0mm	2.7mm	2.0mm Locking Reconstruction Plate	Reconstruction and fixation of defect of mandible
2.4mm Locking	2.0mm	2.7mm	2.0mm Two-Way Compression & Locking Reconstruction Plate	Reconstruction and fixation of defect of mandible
2.4mm	2.0mm	2.7mm	2.5mm General Reconstruction Plate	Fracture of mandible, fixation of osteotomy
2.4mm Locking	2.0mm	2.7mm	2.5mm Locking Reconstruction Plate	Reconstruction and fixation of defect of mandible
2.4mm Locking	2.0mm	2.7mm	2.5mm Two-Way Compression & Locking Reconstruction Plate	Reconstruction and fixation of defect of mandible




1.0mm Straight Plate



Holes	Codes
4	01.04
6	01.06

1.0mm L-Plate, 90°



Direction	Codes
Left	02.03L
Right	02.03R

1.0mm Double Y Plate



Holes	Codes
6	03.06

1.0mm Y-Plate with Bridge



Holes	Codes
4	04.04

1.0mm Y-Plate



Holes	Codes
5	05.05

1.0mm T-Plate



Holes	Codes
6	06.06

1.0mm Square Plate




Hole	Codes
2x2	07.01
3x3	07.02

1.0mm Rectangular Grill Plate, Type-I



Holes	Codes
2x4	08.01

1.0mm Rectangular Grill Plate, Type-II



Holes	Codes
4x6	09.01


1.0mm Rectangular Grill Plate, Type-III



Holes	Codes
2x2	10.01
2x3	10.02
2x4	10.03
4x5	10.04

STERILE
1.0mm Straight Plate


Holes	Codes
4	01.04-S
6	01.06-S

1.0mm L-Plate, 90°


Direction	Codes
Left	02.03L-S
Right	02.03R-S

1.0mm Double Y Plate


Holes	Codes
6	03.06-S

1.0mm Y-Plate with Bridge


Holes	Codes
4	04.04-S

1.0mm Y-Plate


Holes	Codes
5	05.05-S

1.0mm T-Plate



Holes	Codes
6	06.06-S

1.0mm Square Plate


Hole	Codes
2x2	07.01-S
3x3	07.02-S

1.0mm Rectangular Grill Plate, Type-I



Holes	Codes
2x4	08.01-S

1.0mm Rectangular Grill Plate, Type-II



Holes	Codes
4x6	09.01-S

1.0mm Rectangular Grill Plate, Type-III


Holes	Codes
2x2	10.01-S
2x3	10.02-S
2x4	10.03-S
4x5	10.04-S


1.0mm Orbital Floor Plate, Type-I


Codes	11.01
-------	-------

1.0mm Orbital Floor Plate, Type-II



Codes	12.01
-------	-------

1.0mm Orbital Floor Plate, Type-III


Codes	13.01
-------	-------

1.0mm Orbital Floor Plate, Type-IV


Codes	13.11
-------	-------

1.0mm Orbital Floor Plate, Type-V



Codes	14.01
-------	-------

1.0mm Orbital Floor Plate, Type-VI



Direction	Codes
Left	15.01L
Right	15.01R

1.0mm Orbital Floor Plate, Type-VII



Direction	Codes
Left	16.01L
Right	16.01R

1.0mm Orbital Floor Plate, Type-VIII


Direction	Codes
Left	17.01L
Right	17.01R

STERILE
1.0mm Orbital Floor Plate, Type-I


Codes	11.01-S
-------	---------

1.0mm Orbital Floor Plate, Type-II



Codes	12.01-S
-------	---------

1.0mm Orbital Floor Plate, Type-III


Codes	13.01-S
-------	---------

1.0mm Orbital Floor Plate, Type-IV


Codes	13.11-S
-------	---------

1.0mm Orbital Floor Plate, Type-V



Codes	14.01-S
-------	---------

1.0mm Orbital Floor Plate, Type-VI


Direction	Codes
Left	15.01L-S
Right	15.01R-S

1.0mm Orbital Floor Plate, Type-VII


Direction	Codes
Left	16.01L-S
Right	16.01R-S

1.0mm Orbital Floor Plate, Type-VIII


Direction	Codes
Left	17.01L-S
Right	17.01R-S

1.0mm Self Tapping Screws


Length	Codes
2mm	18.02
3mm	18.03
4mm	18.04
5mm	18.05
6mm	18.06
7mm	18.07

STERILE

Length	Codes
2mm	18.02-S
3mm	18.03-S
4mm	18.04-S
5mm	18.05-S
6mm	18.06-S
7mm	18.07-S


1.0mm Self Tapping & Drilling Screws


Length	Codes
2mm	19.02
3mm	19.03
4mm	19.04
5mm	19.05
6mm	19.06
7mm	19.07

STERILE

Length	Codes
2mm	19.02-S
3mm	19.03-S
4mm	19.04-S
5mm	19.05-S
6mm	19.06-S
7mm	19.07-S


1.3mm Emergency Screws, Self Tapping


Length	Codes
2mm	20.02
3mm	20.03
4mm	20.04
5mm	20.05
6mm	20.06
7mm	20.07

STERILE

Length	Codes
2mm	20.02-S
3mm	20.03-S
4mm	20.04-S
5mm	20.05-S
6mm	20.06-S
7mm	20.07-S



7-134-11

Implant and Instrument Set for 1.0mm CMF System



Code	Item	Qty.
01.04	1.0mm Straight Plate, 4 Holes	1
01.06	1.0mm Straight Plate, 6 Holes	1
02.03L	1.0mm L-Plate, 90° - Left	1
02.03R	1.0mm L-Plate, 90° - Right	1
03.06	1.0mm Double Y Plate, 6 Holes	1
04.04	1.0mm Y-Plate with Bridge, 4 Holes	1
05.05	1.0mm Y-Plate , 5 Holes	1
06.06	1.0mm T-Plate, 6 Holes	1
07.01	1.0mm Square Plate, 2 x 2 Holes	1
07.02	1.0mm Square Plate, 2 x 3 Holes	1
08.01	1.0mm Rectangular Grill Plate, Type-I, 2 x 4 Holes	1
09.01	1.0mm Rectangular Grill Plate, Type-II, 4 x 6 Holes	1
10.01	1.0mm Rectangular Grill Plate, Type-III, 2 x 2 Holes	1
10.02	1.0mm Rectangular Grill Plate, Type-III, 2 x 3 Holes	1
10.03	1.0mm Rectangular Grill Plate, Type-III, 2 x 4 Holes	1
10.04	1.0mm Rectangular Grill Plate, Type-III, 4 x 5 Holes	1
11.01	1.0mm Orbital Floor Plate, Type-I	1
12.01	1.0mm Orbital Floor Plate, Type-II	1
13.01	1.0mm Orbital Floor Plate, Type-III	1
13.11	1.0mm Orbital Floor Plate, Type-IV	1
14.01	1.0mm Orbital Floor Plate, Type-V	1
15.01L	1.0mm Orbital Floor Plate, Type-VI - Left	1
15.01R	1.0mm Orbital Floor Plate, Type-VI - Right	1
16.01L	1.0mm Orbital Floor Plate, Type-VII - Left	1
16.01R	1.0mm Orbital Floor Plate, Type-VII - Right	1
17.01L	1.0mm Orbital Floor Plate, Type-VIII - Left	1
17.01R	1.0mm Orbital Floor Plate, Type-VIII - Right	1
Code	1.0mm Bone Screw System	Qty.
18.02	1.0mm Self Tapping Screws, Length 2mm	10
18.03	1.0mm Self Tapping Screws, Length 3mm	10
18.04	1.0mm Self Tapping Screws, Length 4mm	10
18.05	1.0mm Self Tapping Screws, Length 5mm	10
18.06	1.0mm Self Tapping Screws, Length 6mm	10
18.07	1.0mm Self Tapping Screws, Length 7mm	10
19.02	1.0mm Self Tapping & Drilling Screws, Length 2mm	20
19.03	1.0mm Self Tapping & Drilling Screws, Length 3mm	20
19.04	1.0mm Self Tapping & Drilling Screws, Length 4mm	20
19.05	1.0mm Self Tapping & Drilling Screws, Length 5mm	20
19.06	1.0mm Self Tapping & Drilling Screws, Length 6mm	20
19.07	1.0mm Self Tapping & Drilling Screws, Length 7mm	20
20.02	1.3mm Emergency Screws, Self Tapping, Length 2mm	10
20.03	1.3mm Emergency Screws, Self Tapping, Length 3mm	10
20.04	1.3mm Emergency Screws, Self Tapping, Length 4mm	10
20.05	1.3mm Emergency Screws, Self Tapping, Length 5mm	10
20.06	1.3mm Emergency Screws, Self Tapping, Length 6mm	10
20.07	1.3mm Emergency Screws, Self Tapping, Length 7mm	10

Bone Plate & Screw




1.0mm

1.5mm

2.0mm

2.4mm




Neuro
Facial
Mandible
Orthognathic


1.5mm Plate, 2 Holes


Space	Codes
6mm	21.06
8mm	21.08
10mm	21.10
12mm	21.12

1.5mm Straight Plate


Holes	Codes
4	22.04
6	22.06
8	22.08
12	22.12
14	22.14

1.5mm Plate, 4 Holes



Space	Codes
4.5mm	23.01
6.0mm	23.02
7.5mm	23.03
9.0mm	23.04
12.0mm	23.05

1.5mm Plate, 6 Holes



Space	Codes
6.0mm	24.01
7.5mm	24.02
9.0mm	24.03
12.0mm	24.04

1.5mm Continuous Plate



Holes	Codes
6	25.06
8	25.08
16	25.16

1.5mm Orbital Plate, 4 Holes



Space	Codes
6mm	26.01
8mm	26.02

1.5mm Orbital Plate, 6 Holes


Space	Codes
8mm	27.01
12mm	27.02
14mm	27.03

1.5mm Orbital Plate, 7 Holes


Space	Codes
7.5mm	28.01

1.5mm Orbital Plate, 8 Holes


Space	Codes
8mm	29.01

1.5mm Continous Orbital Plate



Holes	Codes
6	30.06
8	30.08
10	30.10

STERILE
1.5mm Plate, 2 Holes



Space	Codes
6mm	21.06-S
8mm	21.08-S
10mm	21.10-S
12mm	21.12-S

1.5mm Straight Plate


Holes	Codes
4	22.04-S
6	22.06-S
8	22.08-S
12	22.12-S
14	22.14-S

1.5mm Plate, 4 Holes


Space	Codes
4.5mm	23.01-S
6.0mm	23.02-S
7.5mm	23.03-S
9.0mm	23.04-S
12.0mm	23.05-S

1.5mm Plate, 6 Holes



Space	Codes
6.0mm	24.01-S
7.5mm	24.02-S
9.0mm	24.03-S
12.0mm	24.04-S

1.5mm Continuous Plate



Holes	Codes
6	25.06-S
8	25.08-S
16	25.16-S

1.5mm Orbital Plate, 4 Holes



Space	Codes
6mm	26.01-S
8mm	26.02-S

1.5mm Orbital Plate, 6 Holes


Space	Codes
8mm	27.01-S
12mm	27.02-S
14mm	27.03-S

1.5mm Orbital Plate, 7 Holes


Space	Codes
7.5mm	28.01-S

1.5mm Orbital Plate, 8 Holes



Space	Codes
8mm	29.01-S

1.5mm Continous Orbital Plate



Holes	Codes
6	30.06-S
8	30.08-S
10	30.10-S


1.5mm L-Plate, 90°



Holes	Direction	Codes
2+2	Left	31.22L
2+2	Right	31.22R
2+3	Left	31.23L
2+3	Right	31.23R
2+4	Left	31.24L
2+4	Right	31.24R
2+5	Left	31.25L
2+5	Right	31.25R
3+5	Left	31.35L
3+5	Right	31.35R

1.5mm L-Plate with Bridge, 90°


Holes	Direction	Space	Codes
2+2	Left	5.5mm	32.01L
2+2	Right	5.5mm	32.01R
2+2	Left	7.5mm	32.02L
2+2	Right	7.5mm	32.02R
2+3	Left	9.0mm	32.03L
2+3	Right	9.0mm	32.03R
2+3	Left	15.0mm	32.04L
2+3	Right	15.0mm	32.04R

1.5mm L-Plate, 130°



Holes	Direction	Codes
2+2	Left	33.01L
2+2	Right	33.01R
2+4	Left	33.02L
2+4	Right	33.02R

1.5mm L-Plate with Bridge, 100°



Holes	Direction	Space	Codes
2+2	Left	5.5mm	34.01L
2+2	Right	5.5mm	34.01R
2+2	Left	7.0mm	34.02L
2+2	Right	7.0mm	34.02R
2+2	Left	9.0mm	34.03L
2+2	Right	9.0mm	34.03R
2+3	Left	9.0mm	34.04L
2+3	Right	9.0mm	34.04R
3+4	Left	5.0mm	34.05L
3+4	Right	5.0mm	34.05R

STERILE
1.5mm L-Plate, 90°



Holes	Direction	Codes
2+2	Left	31.22L-S
2+2	Right	31.22R-S
2+3	Left	31.23L-S
2+3	Right	31.23R-S
2+4	Left	31.24L-S
2+4	Right	31.24R-S
2+5	Left	31.25L-S
2+5	Right	31.25R-S
3+5	Left	31.35L-S
3+5	Right	31.35R-S

1.5mm L-Plate with Bridge, 90°



Holes	Direction	Space	Codes
2+2	Left	5.5mm	32.01L-S
2+2	Right	5.5mm	32.01R-S
2+2	Left	7.5mm	32.02L-S
2+2	Right	7.5mm	32.02R-S
2+3	Left	9.0mm	32.03L-S
2+3	Right	9.0mm	32.03R-S
2+3	Left	15.0mm	32.04L-S
2+3	Right	15.0mm	32.04R-S

1.5mm L-Plate, 130°


Holes	Direction	Codes
2+2	Left	33.01L-S
2+2	Right	33.01R-S
2+4	Left	33.02L-S
2+4	Right	33.02R-S

1.5mm L-Plate with Bridge, 100°



Holes	Direction	Space	Codes
2+2	Left	5.5mm	34.01L-S
2+2	Right	5.5mm	34.01R-S
2+2	Left	7.0mm	34.02L-S
2+2	Right	7.0mm	34.02R-S
2+2	Left	9.0mm	34.03L-S
2+2	Right	9.0mm	34.03R-S
2+3	Left	9.0mm	34.04L-S
2+3	Right	9.0mm	34.04R-S
3+4	Left	5.0mm	34.05L-S
3+4	Right	5.0mm	34.05R-S


1.5mm Double X-Plate, Small


Holes	Space	Codes
4+2	4.0mm	35.01
4+2	5.5mm	35.02
4+2	7.0mm	35.03
4+2	8.5mm	35.04

1.5mm Double X-Plate, Medium



Holes	Codes
4+3	36.01


1.5mm Double X-Plate, Large


Holes	Codes
8+3	37.01

1.5mm X-Plate



Holes	Codes
2+2	38.01

1.5mm Double X-Plate with Support



Holes	Space	Codes
8+2	10mm	39.01
8+2	12mm	39.02
8+2	15mm	39.03

1.5mm Cross Bend Plate


Holes	Codes
2+3	40.01

1.5mm Y-Plate with Bridge, Type-I


Holes	Space	Codes
2+2	8.0mm	41.01
2+2	12.0mm	41.02
2+2	15.0mm	41.03

1.5mm Y-Plate with Bridge, Type-II



Holes	Space	Codes
2+3	5.5mm	42.01
2+3	7.0mm	42.02
2+3	8.5mm	42.03

1.5mm Y-Plate


Holes	Codes
2+3	43.01
2+4	43.02
4+7	43.03

1.5mm H-Plate



Holes	Codes
6+1	44.01
8+1	44.02
10+1	44.03
10+2	44.04

STERILE
1.5mm Double X-Plate, Small


Holes	Space	Codes
4+2	4.0mm	35.01-S
4+2	5.5mm	35.02-S
4+2	7.0mm	35.03-S
4+2	8.5mm	35.04-S

1.5mm Double X-Plate, Medium


Holes	Codes
4+3	36.01-S


1.5mm Double X-Plate, Large

Holes	Codes
8+3	37.01-S

1.5mm X-Plate



Holes	Codes
2+2	38.01-S


1.5mm Double X-Plate with Support

Holes	Space	Codes
8+2	10mm	39.01-S
8+2	12mm	39.02-S
8+2	15mm	39.03-S

1.5mm Cross Bend Plate


Holes	Codes
2+3	40.01-S


1.5mm Y-Plate with Bridge, Type-I


Holes	Space	Codes
2+2	8.0mm	41.01-S
2+2	12.0mm	41.02-S
2+2	15.0mm	41.03-S

1.5mm Y-Plate with Bridge, Type-II

Holes	Space	Codes
2+3	5.5mm	42.01-S
2+3	7.0mm	42.02-S
2+3	8.5mm	42.03-S



1.5mm Y-Plate


Holes	Codes
2+3	43.01-S
2+4	43.02-S
4+7	43.03-S

1.5mm H-Plate

Holes	Codes
6+1	44.01-S
8+1	44.02-S
10+1	44.03-S
10+2	44.04-S





1.5mm H-Plate, Big

Holes	Codes
8	45.08
16	45.16

1.5mm H-Plate, Curve


Holes	Codes
10+1	46.01
10+2	46.02


1.5mm I-Plate


Space	Codes
8mm	47.01
10mm	47.02
12mm	47.03

1.5mm T-Plate

Holes	Codes
3+3	48.01
3+4	48.02
3+7	48.03
7+9	48.04


1.5mm T-Plate with Bridge


Holes	Space	Codes
2+2	9mm	49.01
2+3	9mm	49.02
2+3	12mm	49.03

1.5mm T-Plate, Long


Holes	Codes
2+4	50.01


1.5mm T-Plate, Long-Oblique


Holes	Direction	Codes
2+4	Left	51.01L
2+4	Right	51.01R

1.5mm T-Plate, Oblique

Holes	Direction	Codes
3+3	Left	52.01L
3+3	Right	52.01R



1.5mm Z-Plate


Space	Direction	Codes
9mm	Left	53.01L
9mm	Right	53.01R
11mm	Left	53.02L
11mm	Right	53.02R
13mm	Left	53.03L
13mm	Right	53.03R

1.5mm Square Plate

Holes	Codes
2x2	54.01
2x3	54.02




STERILE

1.5mm H-Plate, Big

Holes	Codes
8	45.08-S
16	45.16-S

1.5mm H-Plate, Curve


Holes	Codes
10+1	46.01-S
10+2	46.02-S


1.5mm I-Plate


Space	Codes
8mm	47.01-S
10mm	47.02-S
12mm	47.03-S

1.5mm T-Plate

Holes	Codes
3+3	48.01-S
3+4	48.02-S
3+7	48.03-S
7+9	48.04-S


1.5mm T-Plate with Bridge


Holes	Space	Codes
2+2	9mm	49.01-S
2+3	9mm	49.02-S
2+3	12mm	49.03-S

1.5mm T-Plate, Long


Holes	Codes
2+4	50.01-S


1.5mm T-Plate, Long-Oblique


Holes	Direction	Codes
2+4	Left	51.01L-S
2+4	Right	51.01R-S

1.5mm T-Plate, Oblique

Holes	Direction	Codes
3+3	Left	52.01L-S
3+3	Right	52.01R-S


1.5mm Z-Plate


Space	Direction	Codes
9mm	Left	53.01L-S
9mm	Right	53.01R-S
11mm	Left	53.02L-S
11mm	Right	53.02R-S
13mm	Left	53.03L-S
13mm	Right	53.03R-S

1.5mm Square Plate

Holes	Codes
2x2	54.01-S
2x3	54.02-S





1.5mm Rectangular Grill Plate, Type-I


Holes	Codes
4X2	55.01
4X6	55.02

1.5mm Rectangular Grill Plate, Type-II


Holes	Codes
2X2	56.01
2X3	56.02
2X4	56.03
4X6	56.04


1.5mm Triangular Plate


Holes	Codes
3	57.01

1.5mm Rectangular Grill Plate, Type-III

Holes	Codes
2X2	58.01
2X3	58.02
2X4	58.03
4X6	58.04


1.5mm Rectangular Grill Plate, Type-IV



Holes	Codes
4X2	55.01
4X6	55.02

1.5mm Rectangular Grill Plate, Type-V



Holes	Codes
11X17	60.1
2X21	60.2
2X11	60.3

1.5mm T-Gap Plate, Angled

Space	Direction	Codes
13.5mm	Left	62.01L
13.5mm	Right	62.01R
15.5mm	Left	62.02L
15.5mm	Right	62.02R
17.5mm	Left	62.03L
17.5mm	Right	62.03R


1.5mm T-Gap Plate, Straight

Space	Codes
13.5mm	61.01
15.5mm	61.02
17.5mm	61.03




STERILE
1.5mm Rectangular Grill Plate, Type-I


Holes	Codes
4X2	55.01-S
4X6	55.02-S

1.5mm Rectangular Grill Plate, Type-II


Holes	Codes
2X2	56.01-S
2X3	56.02-S
2X4	56.03-S
4X6	56.04-S


1.5mm Triangular Plate


Holes	Codes
3	57.01-S

1.5mm Rectangular Grill Plate, Type-III

Holes	Codes
2X2	58.01-S
2X3	58.02-S
2X4	58.03-S
4X6	58.04-S


1.5mm Rectangular Grill Plate, Type-IV



Holes	Codes
4X2	55.01-S
4X6	55.02-S

1.5mm Rectangular Grill Plate, Type-V



Holes	Codes
11X17	60.1-S
2X21	60.2-S
2X11	60.3-S

1.5mm T-Gap Plate, Angled

Space	Direction	Codes
13.5mm	Left	62.01L-S
13.5mm	Right	62.01R-S
15.5mm	Left	62.02L-S
15.5mm	Right	62.02R-S
17.5mm	Left	62.03L-S
17.5mm	Right	62.03R-S



1.5mm T-Gap Plate, Straight

Space	Codes
13.5mm	61.01-S
15.5mm	61.02-S
17.5mm	61.03-S




1.5mm Orbital Floor Plate, Type-III


Codes	65.01
-------	-------

1.5mm Orbital Floor Plate, Type-IV



Codes	66.01
-------	-------

1.5mm Orbital Floor Plate, Type-VI



Direction	Codes
Left	67.01L
Right	67.01R

1.5mm Orbital Floor Plate, Type-VII




Direction	Codes
Left	68.01L
Right	68.01R

1.5mm Orbital Floor Plate, Type-VIII


Direction	Codes
Left	69.01L
Right	69.01R

1.5mm Maxilla Plate, Space



Space	Direction	Codes
3mm	Left	64.01L
3mm	Right	64.01R
5mm	Left	64.02L
5mm	Right	64.02R
7mm	Left	64.03L
7mm	Right	64.03R
9mm	Left	64.04L
9mm	Right	64.04R
11mm	Left	64.05L
11mm	Right	64.05R


1.5mm Chin Plate


Space	Codes
2mm	63.01
4mm	63.02
6mm	63.03
8mm	63.04
10mm	63.05

STERILE
1.5mm Orbital Floor Plate, Type-III


Codes	65.01-S
-------	---------

1.5mm Orbital Floor Plate, Type-IV



Codes	66.01-S
-------	---------

1.5mm Orbital Floor Plate, Type-VI



Direction	Codes
Left	67.01L-S
Right	67.01R-S

1.5mm Orbital Floor Plate, Type-VII




Direction	Codes
Left	68.01L-S
Right	68.01R-S

1.5mm Orbital Floor Plate, Type-VIII


Direction	Codes
Left	69.01L-S
Right	69.01R-S

1.5mm Maxilla Plate, Space


Space	Direction	Codes
3mm	Left	64.01L-S
3mm	Right	64.01R-S
5mm	Left	64.02L-S
5mm	Right	64.02R-S
7mm	Left	64.03L-S
7mm	Right	64.03R-S
9mm	Left	64.04L-S
9mm	Right	64.04R-S
11mm	Left	64.05L-S
11mm	Right	64.05R-S


1.5mm Chin Plate


Space	Codes
2mm	63.01-S
4mm	63.02-S
6mm	63.03-S
8mm	63.04-S
10mm	63.05-S

1.5mm Self Tapping Screws


Length	Codes
3mm	71.01
4mm	71.02
5mm	71.03
6mm	71.04
7mm	71.05
8mm	71.06
9mm	71.07
11mm	71.08
13mm	71.09
15mm	71.10

STERILE

Length	Codes
3mm	71.01-S
4mm	71.02-S
5mm	71.03-S
6mm	71.04-S
7mm	71.05-S
8mm	71.06-S
9mm	71.07-S
11mm	71.08-S
13mm	71.09-S
15mm	71.10-S


1.5mm Self Tapping & Self Drilling Screws


Length	Codes
4mm	70.01
5mm	70.02
6mm	70.03
7mm	70.04

STERILE

Length	Codes
4mm	70.01-S
5mm	70.02-S
6mm	70.03-S
7mm	70.04-S


1.9mm Emergency Screws, Self Tapping


Length	Codes
3mm	72.01
4mm	72.02
5mm	72.03
6mm	72.04
7mm	72.05
8mm	72.06
9mm	72.07
11mm	72.08
13mm	72.09
15mm	72.10

STERILE

Length	Codes
3mm	72.01-S
4mm	72.02-S
5mm	72.03-S
6mm	72.04-S
7mm	72.05-S
8mm	72.06-S
9mm	72.07-S
11mm	72.08-S
13mm	72.09-S
15mm	72.10-S




7-134-21 Implant and Instrument Set for 1.5mm CMF System

Code	Item	Qty.
21.06	1.5mm Plate, 2 Holes, Space 6.0mm	1
21.08	1.5mm Plate, 2 Holes, Space 8.0mm	1
21.10	1.5mm Plate, 2 Holes, Space 10.0mm	1
21.12	1.5mm Plate, 2 Holes, Space 12.0mm	1
22.04	1.5mm Straight Plate, 4 Holes	1
22.06	1.5mm Straight Plate, 6 Holes	1
22.08	1.5mm Straight Plate, 8 Holes	1
22.12	1.5mm Straight Plate, 12 Holes	1
22.14	1.5mm Straight Plate, 14 Holes	1
23.01	1.5mm Plate, 4 Holes, Space 4.5mm	1
23.02	1.5mm Plate, 4 Holes, Space 6.0mm	1
23.03	1.5mm Plate, 4 Holes, Space 7.5mm	1
23.04	1.5mm Plate, 4 Holes, Space 9.0mm	1
23.05	1.5mm Plate, 4 Holes, Space 12.0mm	1
24.01	1.5mm Plate, 6 Holes, Space 6.0mm	1
24.02	1.5mm Plate, 6 Holes, Space 7.5mm	1
24.03	1.5mm Plate, 6 Holes, Space 9.0mm	1
24.04	1.5mm Plate, 6 Holes, Space 12.0mm	1
25.06	1.5mm Continuous Plate, 6 Holes	1
25.08	1.5mm Continuous Plate, 8 Holes	1
25.16	1.5mm Continuous Plate, 16 Holes	1
26.01	1.5mm Orbital Plate, 4 Holes, Space 6.0mm	1
26.02	1.5mm Orbital Plate, 4 Holes, Space 8.0mm	1
27.01	1.5mm Orbital Plate, 6 Holes, Space 8.0mm	1
27.02	1.5mm Orbital Plate, 6 Holes, Space 12.0mm	1
27.03	1.5mm Orbital Plate, 6 Holes, Space 14.0mm	1
29.01	1.5mm Orbital Plate, 8 Holes, Space 8.0mm	1

30.06	1.5mm Continuous Orbital Plate, 6 Holes	1
30.08	1.5mm Continuous Orbital Plate, 8 Holes	1
30.10	1.5mm Continuous Orbital Plate, 10 Holes	1
31.22L	1.5mm L-Plate, 90°, 2+2 Holes - Left	1
31.23L	1.5mm L-Plate, 90°, 2+3 Holes - Left	1
31.24L	1.5mm L-Plate, 90°, 2+4 Holes - Left	1
31.25L	1.5mm L-Plate, 90°, 2+5 Holes - Left	1
31.35L	1.5mm L-Plate, 90°, 3+5 Holes - Left	1
31.22R	1.5mm L-Plate, 90°, 2+2 Holes - Right	1
31.23R	1.5mm L-Plate, 90°, 2+3 Holes - Right	1
31.24R	1.5mm L-Plate, 90°, 2+4 Holes - Right	1
31.25R	1.5mm L-Plate, 90°, 2+5 Holes - Right	1
31.35R	1.5mm L-Plate, 90°, 3+5 Holes - Right	1
32.01R	1.5mm L-Plate with Bridge, 90°, 2+2 Holes, Space 5.5mm - Right	1
32.02R	1.5mm L-Plate with Bridge, 90°, 2+2 Holes, Space 7.5mm - Right	1
32.03R	1.5mm L-Plate with Bridge, 90°, 2+3 Holes, Space 9.0mm - Right	1
32.04R	1.5mm L-Plate with Bridge, 90°, 2+3 Holes, Space 15.0mm - Right	1
32.01L	1.5mm L-Plate with Bridge, 90°, 2+2 Holes, Space 5.5mm - Left	1
32.02L	1.5mm L-Plate with Bridge, 90°, 2+2 Holes, Space 7.5mm - Left	1
32.03L	1.5mm L-Plate with Bridge, 90°, 2+3 Holes, Space 9.0mm - Left	1
32.04L	1.5mm L-Plate with Bridge, 90°, 2+3 Holes, Space 15.0mm - Left	1
33.01R	1.5mm L-Plate, 130°, 2+2 Holes - Right,	1
33.02R	1.5mm L-Plate, 130°, 2+4 Holes - Right	1
33.01L	1.5mm L-Plate, 130°, 2+2 Holes - Left	1
33.02L	1.5mm L-Plate, 130°, 2+4 Holes - Left	1
35.01	1.5mm Double X-Plate, Small, 4+2 Holes, Space 4.0mm	1
35.02	1.5mm Double X-Plate, Small, 4+2 Holes, Space 5.5mm	1
35.03	1.5mm Double X-Plate, Small, 4+2 Holes, Space 7.0mm	1
35.04	1.5mm Double X-Plate, Small, 4+2 Holes, Space 8.5mm	1
36.01	1.5mm Double X-Plate, Medium, 4+3 Holes	1
37.01	1.5mm Double X-Plate, Large, 8+3 Holes	1
38.01	1.5mm X-Plate, 2+2 Holes	1
39.01	1.5mm Double X-Plate with Support, 8+2 Holes, Space 10.0mm	1
39.02	1.5mm Double X-Plate with Support, 8+2 Holes, Space 12.0mm	1
39.03	1.5mm Double X-Plate with Support, 8+2 Holes, Space 15.0mm	1
41.01	1.5mm Y-Plate with Bridge, Type-I, 2+2 Holes, Space 8.0mm	1
41.02	1.5mm Y-Plate with Bridge, Type-I, 2+2 Holes, Space 12.0mm	1
41.03	1.5mm Y-Plate with Bridge, Type-I, 2+2 Holes, Space 15.0mm	1
42.01	1.5mm Y-Plate with Bridge, Type-II, 2+3 Holes, Space 5.5mm	1
42.02	1.5mm Y-Plate with Bridge, Type-II, 2+3 Holes, Space 7.0mm	1
42.03	1.5mm Y-Plate with Bridge, Type-II, 2+3 Holes, Space 8.5mm	1
43.01	1.5mm Y-Plate, 2+3 Holes	1
43.02	1.5mm Y-Plate, 2+4 Holes	1
43.03	1.5mm Y-Plate, 4+7 Holes	1
44.01	1.5mm H-Plate, 6+1 Holes	1
44.02	1.5mm H-Plate, 8+1 Holes	1
44.03	1.5mm H-Plate, 10+1 Holes	1
44.04	1.5mm H-Plate, 10+2 Holes	1
45.08	1.5mm H-Plate, Big, 8 Holes	1
45.16	1.5mm H-Plate, Big, 16 Holes	1
47.01	1.5mm I-Plate, 4 Holes, Space 8.0mm	1
47.02	1.5mm I-Plate, 4 Holes, Space 10.0mm	1
47.03	1.5mm I-Plate, 4 Holes, Space 12.0mm	1

48.01	1.5mm T-Plate, 3+3 Holes	1
48.02	1.5mm T-Plate, 3+4 Holes	1
48.03	1.5mm T-Plate, 5+7 Holes	1
48.04	1.5mm T-Plate, 7+9 Holes	1
49.01	1.5mm T-Plate with Bridge, 2+2 Holes, Space 9.0mm	1
49.02	1.5mm T-Plate with Bridge, 2+3 Holes, Space 9.0mm	1
49.03	1.5mm T-Plate with Bridge, 2+3 Holes, Space 12.0mm	1
50.01	1.5mm T-Plate, Long, 2+4 Holes	1
54.01	1.5mm Square Plate, 2 x 2 Holes	1
54.02	1.5mm Square Plate, 2 x 3 Holes	1
58.01	1.5mm Rectangular Grill Plate, Type-III, 2 x 2 Holes	1
58.02	1.5mm Rectangular Grill Plate, Type-III, 2 x 3 Holes	1
58.03	1.5mm Rectangular Grill Plate, Type-III, 2 x 4 Holes	1
58.04	1.5mm Rectangular Grill Plate, Type-III, 4 x 6 Holes	1
60.01	1.5mm Rectangular Grill Plate, Type-V, 11 x 17 Holes	1
60.02	1.5mm Rectangular Grill Plate, Type-V, 2 x 21 Holes	1
60.03	1.5mm Rectangular Grill Plate, Type-V, 2 x 11 Holes	1
65.01	1.5mm Orbital Floor Plate, Type-III	1
66.01	1.5mm Orbital Floor Plate, Type-IV	1
67.01R	1.5mm Orbital Floor Plate, Type-VI - Right	1
67.01L	1.5mm Orbital Floor Plate, Type-VI - Left	1
68.01R	1.5mm Orbital Floor Plate, Type-VII - Right	1
68.01L	1.5mm Orbital Floor Plate, Type-VII - Left	1
69.01R	1.5mm Orbital Floor Plate, Type-VIII - Right	1
69.01L	1.5mm Orbital Floor Plate, Type-VIII - Left	1
Code	1.5mm Screw System	Qty.
70.01	1.5mm Self Tapping & Self Drilling Screws, Length 4mm	15
70.02	1.5mm Self Tapping & Self Drilling Screws, Length 5mm	15
70.03	1.5mm Self Tapping & Self Drilling Screws, Length 6mm	15
70.04	1.5mm Self Tapping & Self Drilling Screws, Length 7mm	15
71.01	1.5mm Self Tapping Screws, Length 3mm	15
71.02	1.5mm Self Tapping Screws, Length 4mm	15
71.03	1.5mm Self Tapping Screws, Length 5mm	15
71.04	1.5mm Self Tapping Screws, Length 6mm	15
71.05	1.5mm Self Tapping Screws, Length 7mm	15
71.06	1.5mm Self Tapping Screws, Length 8mm	15
71.07	1.5mm Self Tapping Screws, Length 9mm	15
71.08	1.5mm Self Tapping Screws, Length 11mm	15
71.09	1.5mm Self Tapping Screws, Length 13mm	15
71.10	1.5mm Self Tapping Screws, Length 15mm	15
72.01	1.9mm Emergency Screws, Self Tapping, Length 3mm	15
72.02	1.9mm Emergency Screws, Self Tapping, Length 4mm	15
72.03	1.9mm Emergency Screws, Self Tapping, Length 5mm	15
72.04	1.9mm Emergency Screws, Self Tapping, Length 6mm	15
72.05	1.9mm Emergency Screws, Self Tapping, Length 7mm	15
72.06	1.9mm Emergency Screws, Self Tapping, Length 8mm	15
72.07	1.9mm Emergency Screws, Self Tapping, Length 9mm	15
72.08	1.9mm Emergency Screws, Self Tapping, Length 11mm	15
72.09	1.9mm Emergency Screws, Self Tapping, Length 13mm	15
72.10	1.9mm Emergency Screws, Self Tapping, Length 15mm	15

Bone Plate & Screw




1.0mm

1.5mm

2.0mm


2.4mm




Neuro
Facial
Mandible
Orthognathic


2.0mm Straight Plate



Holes	Codes (0.6mm)
4	73.04
5	73.05
6	73.06
14	73.14
16	73.16

2.0mm Straight Plate with Bridge, Type-I


Holes	Space	Codes (0.6mm)
4	9mm	74.01
4	12mm	74.02
6	9mm	74.03
6	12mm	74.04

2.0mm Straight Plate 2 Holes, Type-I


Space	Codes (0.6mm)
8mm	75.01
10mm	75.02
12mm	75.03

2.0mm Straight Plate with Bridge, Type-II



Holes	Space	Codes (0.6mm)
4	6mm	76.01
4	8mm	76.02
4	10mm	76.03
4	12mm	76.04
8	8mm	76.05
8	10mm	76.06
8	12mm	76.07

2.0mm Straight Plate 5 Holes, Type-I



Codes (0.6mm)	77.01
---------------	-------

2.0mm Straight Plate 6 Holes, Type-I


Codes (0.6mm)	78.01
---------------	-------

2.0mm L-Plate with bridge, Type-I, 90°



Holes	Direction	Space	Codes (0.6mm)
2+2	Left	8mm	79.01L
2+2	Right	8mm	79.01R
2+2	Left	10mm	79.02L
2+2	Right	10mm	79.02R
2+2	Left	12mm	79.03L
2+2	Right	12mm	79.03R
2+2	Left	15mm	79.04L
2+2	Right	15mm	79.04R
2+3	Left	9mm	79.05L
2+3	Right	9mm	79.05R
2+3	Left	13mm	79.06L
2+3	Right	13mm	79.06R

2.0mm L-Plate, Type-I, 90°



Holes	Direction	Codes (0.6mm)
2+2	Left	80.01L
2+2	Right	80.01R
2+4	Left	80.02L
2+4	Right	80.02R

STERILE
2.0mm Straight Plate



Holes	Codes (0.6mm)
4	73.04-S
5	73.05-S
6	73.06-S
14	73.14-S
16	73.16-S

2.0mm Straight Plate with Bridge, Type-I


Holes	Space	Codes (0.6mm)
4	9mm	74.01-S
4	12mm	74.02-S
6	9mm	74.03-S
6	12mm	74.04-S

2.0mm Straight Plate 2 Holes, Type-I


Space	Codes (0.6mm)
8mm	75.01-S
10mm	75.02-S
12mm	75.03-S

2.0mm Straight Plate with Bridge, Type-II



Holes	Space	Codes (0.6mm)
4	6mm	76.01-S
4	8mm	76.02-S
4	10mm	76.03-S
4	12mm	76.04-S
8	8mm	76.05-S
8	10mm	76.06-S
8	12mm	76.07-S

2.0mm Straight Plate 5 Holes, Type-I



Codes (0.6mm)	77.01-S
---------------	---------

2.0mm Straight Plate 6 Holes, Type-I



Codes (0.6mm)	78.01-S
---------------	---------

2.0mm L-Plate with bridge, Type-I, 90°



Holes	Direction	Space	Codes (0.6mm)
2+2	Left	8mm	79.01L-S
2+2	Right	8mm	79.01R-S
2+2	Left	10mm	79.02L-S
2+2	Right	10mm	79.02R-S
2+2	Left	12mm	79.03L-S
2+2	Right	12mm	79.03R-S
2+2	Left	15mm	79.04L-S
2+2	Right	15mm	79.04R-S
2+3	Left	9mm	79.05L-S
2+3	Right	9mm	79.05R-S
2+3	Left	13mm	79.06L-S
2+3	Right	13mm	79.06R-S

2.0mm L-Plate, Type-I, 90°


Holes	Direction	Codes (0.6mm)
2+2	Left	80.01L-S
2+2	Right	80.01R-S
2+4	Left	80.02L-S
2+4	Right	80.02R-S


2.0mm L-Plate with Bridge, Type-II, 90°





Holes	Direction	Space	Codes (0.6mm)
2+2	Left	6mm	81.01L
2+2	Right	6mm	81.01R
2+2	Left	8mm	81.02L
2+2	Right	8mm	81.02R
2+2	Left	10mm	81.03L
2+2	Right	10mm	81.03R
3+3	Left	10mm	81.04L
3+3	Right	10mm	81.04R

2.0mm L-Oblique Plate with Bridge, 90°



Holes	Direction	Space	Codes (0.6mm)
2+2	Left	6mm	82.01L
2+2	Right	6mm	82.01R
2+2	Left	8mm	82.02L
2+2	Right	8mm	82.02R
2+2	Left	10mm	82.03L
2+2	Right	10mm	82.03R
3+3	Left	10mm	82.04L
3+3	Right	10mm	82.04R

2.0mm H-Plate, 11 Holes

Codes (0.6mm)	83.01
---------------	-------


2.0mm Straight Plate



Holes	Codes (0.8mm)	Codes (1.0mm)
4	84.01	84.09
5	84.02	84.10
6	84.03	84.11
8	84.04	84.12
10	84.05	84.13
12	84.06	84.14
14	84.07	84.15
16	84.08	84.16

2.0mm Straight Plate with Bridge, Type-III


Holes	Space	Codes (0.8mm)	Codes (1.0mm)
4	6mm	85.01	85.08
4	8mm	85.02	85.09
4	10mm	85.03	85.10
4	12mm	85.04	85.11
6	8mm	85.05	85.12
6	10mm	85.06	85.13
6	12mm	85.07	85.14


STERILE

2.0mm L-Plate with Bridge, Type-II, 90°



Holes	Direction	Space	Codes (0.6mm)
2+2	Left	6mm	81.01L-S
2+2	Right	6mm	81.01R-S
2+2	Left	8mm	81.02L-S
2+2	Right	8mm	81.02R-S
2+2	Left	10mm	81.03L-S
2+2	Right	10mm	81.03R-S
3+3	Left	10mm	81.04L-S
3+3	Right	10mm	81.04R-S




2.0mm L-Oblique Plate with Bridge, 90°



Holes	Direction	Space	Codes (0.6mm)
2+2	Left	6mm	82.01L-S
2+2	Right	6mm	82.01R-S
2+2	Left	8mm	82.02L-S
2+2	Right	8mm	82.02R-S
2+2	Left	10mm	82.03L-S
2+2	Right	10mm	82.03R-S
3+3	Left	10mm	82.04L-S
3+3	Right	10mm	82.04R-S

2.0mm H-Plate, 11 Holes

Codes (0.6mm)	83.01-S
---------------	---------




2.0mm Straight Plate



Holes	Codes (0.8mm)	Codes (1.0mm)
4	84.01-S	84.09-S
5	84.02-S	84.10-S
6	84.03-S	84.11-S
8	84.04-S	84.12-S
10	84.05-S	84.13-S
12	84.06-S	84.14-S
14	84.07-S	84.15-S
16	84.08-S	84.16-S

2.0mm Straight Plate with Bridge, Type-III



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
4	6mm	85.01-S	85.08-S
4	8mm	85.02-S	85.09-S
4	10mm	85.03-S	85.10-S
4	12mm	85.04-S	85.11-S
6	8mm	85.05-S	85.12-S
6	10mm	85.06-S	85.13-S
6	12mm	85.07-S	85.14-S



2.0mm Straight Plate, 5 Holes, Type-II



Codes (0.8mm)	Codes (1.0mm)
86.01	86.02

2.0mm Straight Plate, 6 Holes, Type-II




Codes (0.8mm)	Codes (1.0mm)
87.01	87.02

2.0mm Straight Plate, 2 Holes, Type-II




Space	Codes (0.8mm)	Codes (1.0mm)
8mm	88.01	88.04
10mm	88.02	88.05
12mm	88.03	88.06

2.0mm Straight Plate with Bridge, Type IV




Holes	Space	Codes (0.8mm)	Codes (1.0mm)
4	9mm	89.01	89.05
4	12mm	89.02	89.06
6	9mm	89.03	89.07
6	12mm	89.04	89.08

2.0mm Orbital Plate with Bridge, Type-I



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
4	8mm	90.01	90.06
4	12mm	90.02	90.07
6	8mm	90.03	90.08
6	12mm	90.04	90.09
6	15mm	90.05	90.10

2.0mm Orbital Plate with Bridge, Type-II




Holes	Space	Codes (0.8mm)	Codes (1.0mm)
5	7.5mm	91.01	91.03
7	7.5mm	91.02	91.04

2.0mm Orbital Plate



Holes	Codes (0.8mm)	Codes (1.0mm)
6	92.01	92.04
8	92.02	92.05
12	92.03	92.06

2.0mm L-Plate, Type-II, 90°



Holes	Direction	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	93.01	93.05
2+2	Right	93.02	93.06
2+4	Left	93.03	93.07
2+4	Right	93.04	93.08

STERILE
2.0mm Straight Plate, 5 Holes, Type-II



Codes (0.8mm)	Codes (1.0mm)
86.01-S	86.02-S

2.0mm Straight Plate, 6 Holes, Type-II



Codes (0.8mm)	Codes (1.0mm)
87.01-S	87.02-S

2.0mm Straight Plate, 2 Holes, Type-II



Space	Codes (0.8mm)	Codes (1.0mm)
8mm	88.01-S	88.04-S
10mm	88.02-S	88.05-S
12mm	88.03-S	88.06-S

2.0mm Straight Plate with Bridge, Type IV


Holes	Space	Codes (0.8mm)	Codes (1.0mm)
4	9mm	89.01-S	89.05-S
4	12mm	89.02-S	89.06-S
6	9mm	89.03-S	89.07-S
6	12mm	89.04-S	89.08-S

2.0mm Orbital Plate with Bridge, Type-I



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
4	8mm	90.01-S	90.06-S
4	12mm	90.02-S	90.07-S
6	8mm	90.03-S	90.08-S
6	12mm	90.04-S	90.09-S
6	15mm	90.05-S	90.10-S

2.0mm Orbital Plate with Bridge, Type-II



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
5	7.5mm	91.01-S	91.03-S
7	7.5mm	91.02-S	91.04-S

2.0mm Orbital Plate



Holes	Codes (0.8mm)	Codes (1.0mm)
6	92.01-S	92.04-S
8	92.02-S	92.05-S
12	92.03-S	92.06-S

2.0mm L-Plate, Type-II, 90°



Holes	Direction	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	93.01-S	93.05-S
2+2	Right	93.02-S	93.06-S
2+4	Left	93.03-S	93.07-S
2+4	Right	93.04-S	93.08-S


2.0mm L-Plate with Bridge, Type-III, 90°



Holes	Direction	Space	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	8mm	94.01	94.15
2+2	Right	8mm	94.02	94.16
2+2	Left	10mm	94.03	94.17
2+2	Right	10mm	94.04	94.18
2+2	Left	12mm	94.05	94.19
2+2	Right	12mm	94.06	94.20
2+2	Left	15mm	94.07	94.21
2+2	Right	15mm	94.08	94.22
3+3	Left	9mm	94.09	94.23
3+3	Right	9mm	94.10	94.24
3+3	Left	13mm	94.11	94.25
3+3	Right	13mm	94.12	94.26
2+3	Left	12mm	94.13	94.27
2+3	Right	12mm	94.14	94.28

2.0mm L-Plate with Bridge, 100°



Holes	Direction	Space	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	6mm	95.01	95.13
2+2	Right	6mm	95.02	95.14
2+2	Left	8mm	95.03	95.15
2+2	Right	8mm	95.04	95.16
2+2	Left	12mm	95.05	95.17
2+2	Right	12mm	95.06	95.18
3+3	Left	10mm	95.07	95.19
3+3	Right	10mm	95.08	95.20
3+4	Left	12mm	95.09	95.21
3+4	Right	12mm	95.10	95.22
4+6	Left	12mm	95.11	95.23
4+6	Right	12mm	95.12	95.24

2.0mm L-Plate, 130°



Holes	Direction	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	96.01	96.05
2+2	Right	96.02	96.06
2+4	Left	96.03	96.07
2+4	Right	96.04	96.08

2.0mm L-Plate with Bridge, Type-IV, 90°



Holes	Direction	Space	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	6mm	97.01	97.07
2+2	Right	6mm	97.02	97.08
2+2	Left	8mm	97.03	97.09
2+2	Right	8mm	97.04	97.10
2+2	Left	10mm	97.05	97.11
2+2	Right	10mm	97.06	97.12

STERILE
2.0mm L-Plate with Bridge, Type-III, 90°



Holes	Direction	Space	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	8mm	94.01-S	94.15-S
2+2	Right	8mm	94.02-S	94.16-S
2+2	Left	10mm	94.03-S	94.17-S
2+2	Right	10mm	94.04-S	94.18-S
2+2	Left	12mm	94.05-S	94.19-S
2+2	Right	12mm	94.06-S	94.20-S
2+2	Left	15mm	94.07-S	94.21-S
2+2	Right	15mm	94.08-S	94.22-S
3+3	Left	9mm	94.09-S	94.23-S
3+3	Right	9mm	94.10-S	94.24-S
3+3	Left	13mm	94.11-S	94.25-S
3+3	Right	13mm	94.12-S	94.26-S
2+3	Left	12mm	94.13-S	94.27-S
2+3	Right	12mm	94.14-S	94.28-S

2.0mm L-Plate with Bridge, 100°



Holes	Direction	Space	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	6mm	95.01-S	95.13-S
2+2	Right	6mm	95.02-S	95.14-S
2+2	Left	8mm	95.03-S	95.15-S
2+2	Right	8mm	95.04-S	95.16-S
2+2	Left	12mm	95.05-S	95.17-S
2+2	Right	12mm	95.06-S	95.18-S
3+3	Left	10mm	95.07-S	95.19-S
3+3	Right	10mm	95.08-S	95.20-S
3+4	Left	12mm	95.09-S	95.21-S
3+4	Right	12mm	95.10-S	95.22-S
4+6	Left	12mm	95.11-S	95.23-S
4+6	Right	12mm	95.12-S	95.24-S

2.0mm L-Plate, 130°



Holes	Direction	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	96.01-S	96.05-S
2+2	Right	96.02-S	96.06-S
2+4	Left	96.03-S	96.07-S
2+4	Right	96.04-S	96.08-S

2.0mm L-Plate with Bridge, Type-IV, 90°



Holes	Direction	Space	Codes (0.8mm)	Codes (1.0mm)
2+2	Left	6mm	97.01-S	97.07-S
2+2	Right	6mm	97.02-S	97.08-S
2+2	Left	8mm	97.03-S	97.09-S
2+2	Right	8mm	97.04-S	97.10-S
2+2	Left	10mm	97.05-S	97.11-S
2+2	Right	10mm	97.06-S	97.12-S


2.0mm Double X-Plate, Small



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
4+2	6mm	98.01	98.05
4+2	8mm	98.02	98.06
4+4	10mm	98.03	98.07
4+4	12mm	98.04	98.08

2.0mm H-Plate



Space	Codes (0.8mm)	Codes (1.0mm)
6+1	112.01	112.04
8+1	112.02	112.05
10+1	112.03	112.06

2.0mm Long T-Plate, Oblique Angle



Direction	Codes (0.8mm)	Codes (1.0mm)
Left	115.01	115.03
Right	115.02	115.04

2.0mm I-Plate



Space	Codes (0.8mm)	Codes (1.0mm)
6mm	113.01	113.04
9mm	113.02	113.05

2.0mm T-Plate with Bridge



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
2+2	6mm	114.01	114.08
2+2	8mm	114.02	114.09
2+2	12mm	114.03	114.10
2+4	6mm	114.04	114.11
2+3	8mm	114.05	114.12
2+3	15mm	114.06	114.13
2+6	6mm	114.07	114.14

2.0mm T-Plate



Direction	Codes (0.8mm)	Codes (1.0mm)
3+3	116.01	116.04
3+4	116.02	116.05
3+6	116.03	116.06

2.0mm T-Plate with Bridge



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
3+2	8mm	117.01	117.03
3+2	10mm	117.02	117.04

2.0mm Z-Platel



Space	Space	Codes (0.8mm)	Codes (1.0mm)
8mm	Left	120.01	120.05
8mm	Right	120.02	120.06
12mm	Left	120.03	120.07
12mm	Right	120.04	120.08

STERILE
2.0mm Double X-Plate, Small



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
4+2	6mm	98.01-S	98.05-S
4+2	8mm	98.02-S	98.06-S
4+4	10mm	98.03-S	98.07-S
4+4	12mm	98.04-S	98.08-S

2.0mm H-Plate



Space	Codes (0.8mm)	Codes (1.0mm)
6+1	112.01-S	112.04-S
8+1	112.02-S	112.05-S
10+1	112.03-S	112.06-S

2.0mm Long T-Plate, Oblique Angle



Direction	Codes (0.8mm)	Codes (1.0mm)
Left	115.01-S	115.03-S
Right	115.02-S	115.04-S

2.0mm I-Plate



Space	Codes (0.8mm)	Codes (1.0mm)
6mm	113.01-S	113.04-S
9mm	113.02-S	113.05-S

2.0mm T-Plate with Bridge



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
2+2	6mm	114.01-S	114.08-S
2+2	8mm	114.02-S	114.09-S
2+2	12mm	114.03-S	114.10-S
2+4	6mm	114.04-S	114.11-S
2+3	8mm	114.05-S	114.12-S
2+3	15mm	114.06-S	114.13-S
2+6	6mm	114.07-S	114.14-S

2.0mm T-Plate


Direction	Codes (0.8mm)	Codes (1.0mm)
3+3	116.01-S	116.04-S
3+4	116.02-S	116.05-S
3+6	116.03-S	116.06-S

2.0mm T-Plate with Bridge



Holes	Space	Codes (0.8mm)	Codes (1.0mm)
3+2	8mm	117.01-S	117.03-S
3+2	10mm	117.02-S	117.04-S

2.0mm Z-Platel



Space	Space	Codes (0.8mm)	Codes (1.0mm)
8mm	Left	120.01-S	120.05-S
8mm	Right	120.02-S	120.06-S
12mm	Left	120.03-S	120.07-S
12mm	Right	120.04-S	120.08-S


2.0mm Pseudo-Possessed Plate, Type-III



Holes	Codes (0.8mm)	Codes (1.0mm)
7	129.01	129.02

2.0mm Rectangular Grill Plate, Type-II



Holes	Codes (0.8mm)	Codes (1.0mm)
4x2	135.01	135.04
4x4	135.02	135.05
4x6	135.03	135.06

2.0mm Rectangular Grill Plate, Type-I


Holes	Codes (0.8mm)	Codes (1.0mm)
2x2	134.01	134.03
2x3	134.02	134.04

2.0mm Rectangular Grill Plate, Type -III



Holes	Codes (0.8mm)	Codes (1.0mm)
2x4	136.01	136.02

2.0mm Rectangular Grill Plate, Type -IV



Holes	Direction	Codes (0.8mm)	Codes (1.0mm)
2x4	Left	137.01	137.03
2x4	Right	137.02	137.04

STERILE
2.0mm Pseudo-Possessed Plate, Type-III



Holes	Codes (0.8mm)	Codes (1.0mm)
7	129.01-S	129.02-S

2.0mm Rectangular Grill Plate, Type-II


Holes	Codes (0.8mm)	Codes (1.0mm)
4x2	135.01-S	135.04-S
4x4	135.02-S	135.05-S
4x6	135.03-S	135.06-S

2.0mm Rectangular Grill Plate, Type-I


Holes	Codes (0.8mm)	Codes (1.0mm)
2x2	134.01-S	134.03-S
2x3	134.02-S	134.04-S

2.0mm Rectangular Grill Plate, Type-III


Holes	Codes (0.8mm)	Codes (1.0mm)
2x4	136.01-S	136.02-S

2.0mm Rectangular Grill Plate, Type-IV


Holes	Direction	Codes (0.8mm)	Codes (1.0mm)
2x4	Left	137.01-S	137.03-S
2x4	Right	137.02-S	137.04-S



LOCKING SYSTEM

2.0mm Locking Straight Plate with Bridge, Type-I



Holes	Codes (1.0mm)
4	143.01
6	143.02
8	143.03
12	143.04
14	143.05

2.0mm Locking T-Gap Plate, Straight




Holes	Space	Codes (1.0mm)
4	10mm	144.01
4	12mm	144.02
6	10mm	144.03
6	12mm	144.04

2.0mm Locking U-Plate




Holes	Codes (1.0mm)
3+3	145.01
4+4	145.02



2.0mm Locking U-Plate


Holes	Codes (1.0mm)
9	146.01

2.0mm Locking Square Frame Plate, Type-I




Holes	Codes (1.0mm)
2x4	147.01
2x6	147.02
2x8	147.03

2.0mm Locking Rectangular Grill Plate, Type-I



Holes	Codes (1.0mm)
4x2	148.01
4x6	148.02

2.0mm Locking Rectangular Grill Plate, Type-II



Holes	Codes (1.0mm)
2x3	149.01
2x4	149.02
2x5	149.03
2x6	149.04

LOCKING SYSTEM

STERILE

2.0mm Locking Straight Plate with Bridge, Type-I



Holes	Codes (1.0mm)
4	143.01-S
6	143.02-S
8	143.03-S
12	143.04-S
14	143.05-S

2.0mm Locking T-Gap Plate, Straight




Holes	Space	Codes (1.0mm)
4	10mm	144.01-S
4	12mm	144.02-S
6	10mm	144.03-S
6	12mm	144.04-S

2.0mm Locking U-Plate




Holes	Codes (1.0mm)
3+3	145.01-S
4+4	145.02-S



2.0mm Locking U-Plate


Holes	Codes (1.0mm)
9	146.01-S

2.0mm Locking Square Frame Plate, Type-I




Holes	Codes (1.0mm)
2x4	147.01-S
2x6	147.02-S
2x8	147.03-S

2.0mm Locking Rectangular Grill Plate, Type-I



Holes	Codes (1.0mm)
4x2	148.01-S
4x6	148.02-S


2.0mm Locking Rectangular Grill Plate, Type-II




Holes	Codes (1.0mm)
2x3	149.01-S
2x4	149.02-S
2x5	149.03-S
2x6	149.04-S


2.0mm Locking Straight Plate with Bridge, Type-II


Holes	Codes (1.2mm)
4	150.01

2.0mm Locking Rectangular Grill Plate, Type-III



Holes	Codes (1.2mm)
2x3	151.01
2x4	151.02
2x5	151.03
2x6	151.04

2.0mm Locking Straight Plate with Bridge, Type-III


Holes	Space	Codes (1.2mm)
4	10mm	152.01
4	12mm	152.02
6	10mm	152.03
6	12mm	152.04

2.0mm Locking Straight Plate, Type-II



Holes	Codes (1.2mm)
4	154.01
6	154.02
8	154.03
10	154.04
12	154.05
14	154.06
16	154.07

2.0mm Locking Straight Plate with Bridge, Type-IV



Holes	Space	Codes (1.2mm)
4	8mm	153.01
4	10mm	153.02
4	12mm	153.03
6	8mm	153.04
6	10mm	153.05
6	12mm	153.06

STERILE
2.0mm Locking Straight Plate with Bridge, Type-II


Holes	Codes (1.2mm)
4	150.01-S

2.0mm Locking Rectangular Grill Plate, Type-III



Holes	Codes (1.2mm)
2x3	151.01-S
2x4	151.02-S
2x5	151.03-S
2x6	151.04-S

2.0mm Locking Straight Plate with Bridge, Type-III



Holes	Space	Codes (1.2mm)
4	10mm	152.01-S
4	12mm	152.02-S
6	10mm	152.03-S
6	12mm	152.04-S

2.0mm Locking Straight Plate, Type-II



Holes	Codes (1.2mm)
4	154.01-S
6	154.02-S
8	154.03-S
10	154.04-S
12	154.05-S
14	154.06-S
16	154.07-S

2.0mm Locking Straight Plate with Bridge, Type-IV



Holes	Space	Codes (1.2mm)
4	8mm	153.01-S
4	10mm	153.02-S
4	12mm	153.03-S
6	8mm	153.04-S
6	10mm	153.05-S
6	12mm	153.06-S


2.0mm Locking Straight Plate with Bridge, Type-V



Holes	Space	Codes (1.5mm)
4	8mm	155.01
4	10mm	155.02
4	12mm	155.03
6	8mm	155.04

2.0mm Locking Straight Plate, Type-III


Holes	Space	Codes (1.5mm)
4	10mm	157.01
4	15mm	157.02

2.0mm One-Way Pressure Plate with Bridge



Holes	Space	Codes (1.5mm)
4	8mm	159.01
4	10mm	159.02
4	12mm	159.03
6	8mm	159.04
6	10mm	159.05

2.0mm One -Way Pressure Curved Plate


Holes	Direction	Space	Codes (1.5mm)
4	Left	6mm	162.01
4	Right	6mm	162.02
4	Left	10.5mm	162.03
4	Right	10.5mm	162.04
6	Left	10.5mm	162.05
6	Right	10.5mm	162.06
8	Left	10.5mm	162.07
8	Right	10.5mm	162.08

2.0mm Locking Straight Plate, Type-III



Holes	Codes (1.5mm)
6	156.01
8	156.02
26	156.03

2.0mm Locking Curved Plate, Type-I



Holes	Direction	Space	Codes (1.5mm)
4	Left	6mm	158.01
4	Right	6mm	158.02
4	Left	10.5mm	158.03
4	Right	10.5mm	158.04
6	Left	10.5mm	158.05
6	Right	10.5mm	158.06
8	Left	10.5mm	158.07
8	Right	10.5mm	158.08

2.0mm One -Way Pressure Plate



Holes	Codes (1.5mm)
6	160.01
8	160.02
26	160.03

2.0mm One -Way Pressure Plate



Holes	Space	Codes (1.5mm)
4	10mm	161.01
4	15mm	161.02

STERILE
2.0mm Locking Straight Plate with Bridge, Type-V



Holes	Space	Codes (1.5mm)
4	8mm	155.01-S
4	10mm	155.02-S
4	12mm	155.03-S
6	8mm	155.04-S

2.0mm Locking Straight Plate, Type-III


Holes	Space	Codes (1.5mm)
4	10mm	157.01-S
4	15mm	157.02-S

2.0mm One-Way Pressure Plate with Bridge



Holes	Space	Codes (1.5mm)
4	8mm	159.01-S
4	10mm	159.02-S
4	12mm	159.03-S
6	8mm	159.04-S
6	10mm	159.05-S

2.0mm One -Way Pressure Curved Plate


Holes	Direction	Space	Codes (1.5mm)
4	Left	6mm	162.01-S
4	Right	6mm	162.02-S
4	Left	10.5mm	162.03-S
4	Right	10.5mm	162.04-S
6	Left	10.5mm	162.05-S
6	Right	10.5mm	162.06-S
8	Left	10.5mm	162.07-S
8	Right	10.5mm	162.08-S

2.0mm Locking Straight Plate, Type-III



Holes	Codes (1.5mm)
6	156.01-S
8	156.02-S
26	156.03-S

2.0mm Locking Curved Plate, Type-I


Holes	Direction	Space	Codes (1.5mm)
4	Left	6mm	158.01-S
4	Right	6mm	158.02-S
4	Left	10.5mm	158.03-S
4	Right	10.5mm	158.04-S
6	Left	10.5mm	158.05-S
6	Right	10.5mm	158.06-S
8	Left	10.5mm	158.07-S
8	Right	10.5mm	158.08-S

2.0mm One -Way Pressure Plate


Holes	Codes (1.5mm)
6	160.01-S
8	160.02-S
26	160.03-S

2.0mm One -Way Pressure Plate


Holes	Space	Codes (1.5mm)
4	10mm	161.01-S
4	15mm	161.02-S

2.0mm Self Tapping Screws


Length	Codes
4mm	163.04
5mm	163.05
7mm	163.07
9mm	163.09
11mm	163.11
13mm	163.13
15mm	163.15
17mm	163.17
19mm	163.19

STERILE

Length	Codes
4mm	163.04-S
5mm	163.05-S
7mm	163.07-S
9mm	163.09-S
11mm	163.11-S
13mm	163.13-S
15mm	163.15-S
17mm	163.17-S
19mm	163.19-S


2.0mm Self Tapping & Self Drilling Screws


Length	Codes
4mm	164.04
5mm	164.05
7mm	164.07
9mm	164.09
11mm	164.11
13mm	164.13
15mm	164.15
17mm	164.17
19mm	164.19

STERILE

Length	Codes
4mm	164.04-S
5mm	164.05-S
7mm	164.07-S
9mm	164.09-S
11mm	164.11-S
13mm	164.13-S
15mm	164.15-S
17mm	164.17-S
19mm	164.19-S


2.0mm Locking Self Tapping Screws


Length	Codes
5mm	165.05
7mm	165.07
9mm	165.09
11mm	165.11
13mm	165.13

STERILE

Length	Codes
5mm	165.05-S
7mm	165.07-S
9mm	165.09-S
11mm	165.11-S
13mm	165.13-S



2.0mm Locking Self Tapping & Drilling Screws



Length	Codes
5mm	166.05
7mm	166.07
9mm	166.09
11mm	166.11
13mm	166.13

STERILE

Length	Codes
5mm	166.05-S
7mm	166.07-S
9mm	166.09-S
11mm	166.11-S
13mm	166.13-S



2.3mm Emergency Screws, Self Tapping




Length	Codes
5mm	167.05
7mm	167.07
9mm	167.09
11mm	167.11
13mm	167.13
15mm	167.15
17mm	167.17
19mm	167.19

STERILE

Length	Codes
5mm	167.05-S
7mm	167.07-S
9mm	167.09-S
11mm	167.11-S
13mm	167.13-S
15mm	167.15-S
17mm	167.17-S
19mm	167.19-S




7-134-32 Implant and Instrument Set for 2.0mm CMF System

0.8mm	Item	Qty.
84.01	2.0mm Straight Plate, 4 Holes	1
84.02	2.0mm Straight Plate, 5 Holes	1
84.03	2.0mm Straight Plate, 6 Holes	1
84.04	2.0mm Straight Plate, 8 Holes	1
84.05	2.0mm Straight Plate, 1 Holes	1
84.06	2.0mm Straight Plate, 12 Holes	1
84.07	2.0mm Straight Plate, 14 Holes	1
84.08	2.0mm Straight Plate, 16 Holes	1
90.01	2.0mm Orbital Plate with Bridge, Type-I, 4 Holes, Space 8.0mm	1
90.02	2.0mm Orbital Plate with Bridge, Type-I, 4 Holes, Space 8.0mm	1
90.03	2.0mm Orbital Plate with Bridge, Type-I, 6 Holes, Space 8.0mm	1
90.04	2.0mm Orbital Plate with Bridge, Type-I, 6 Holes, Space 8.0mm	1
90.05	2.0mm Orbital Plate with Bridge, Type-I, 6 Holes, Space 8.0mm	1
92.01	2.0mm Orbital Plate, Space 6.0mm	1
92.02	2.0mm Orbital Plate, Space 8.0mm	1
92.03	2.0mm Orbital Plate, Space 12.0mm	1
93.01	2.0mm L-Plate, Type-II, 90°, 2+2 Holes - Left	1
93.03	2.0mm L-Plate, Type-II, 90°, 2+4 Holes - Left	1
93.02	2.0mm L-Plate, Type-II, 90°, 2+2 Holes - Right	1

0.8mm	Item	Qty.
93.04	2.0mm L-Plate, Type-II, 90°, 2+4 Holes - Right	1
94.01	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 8.0mm - Left	1
94.03	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 10.0mm - Left	1
94.05	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 12.0mm - Left	1
94.07	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 15.0mm - Left	1
94.09	2.0mm L-Plate with Bridge, Type-III, 90°, 3+3 Holes, Space 9.0mm - Left	1
94.11	2.0mm L-Plate with Bridge, Type-III, 90°, 3+3 Holes, Space 13.0mm - Left	1
94.13	2.0mm L-Plate with Bridge, Type-III, 90°, 2+3 Holes, Space 12.0mm - Left	1
94.02	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 8.0mm - Right	1
94.04	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 10.0mm - Right	1
94.06	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 12.0mm - Right	1
94.08	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 15.0mm - Right	1
94.10	2.0mm L-Plate with Bridge, Type-III, 90°, 3+3 Holes, Space 9.0mm - Right	1
94.12	2.0mm L-Plate with Bridge, Type-III, 90°, 3+3 Holes, Space 13.0mm - Right	1
94.14	2.0mm L-Plate with Bridge, Type-III, 90°, 2+3 Holes, Space 12.0mm - Right	1
98.01	2.0mm Double X-Plate, Small, 4+2 Holes, Space 6.0mm	1
98.02	2.0mm Double X-Plate, Small, 4+2 Holes, Space 8.0mm	1
98.03	2.0mm Double X-Plate, Small, 4+4 Holes, Space 10.0mm	1
98.04	2.0mm Double X-Plate, Small, 4+4 Holes, Space 12.0mm	1
112.01	2.0mm H-Plate, 6+1 Holes	1
112.02	2.0mm H-Plate, 8+1 Holes	1
112.03	2.0mm H-Plate, 10+1 Holes	1
113.01	2.0mm I-Plate, Space 6.0mm	1
113.02	2.0mm I-Plate, Space 9.0mm	1
114.01	2.0mm T-Plate with Bridge, 2+2 Holes, Space 6.0mm	1
114.02	2.0mm T-Plate with Bridge, 2+2 Holes, Space 8.0mm	1
114.03	2.0mm T-Plate with Bridge, 2+2 Holes, Space 12.0mm	1
114.04	2.0mm T-Plate with Bridge, 2+4 Holes, Space 6.0mm	1
114.05	2.0mm T-Plate with Bridge, 2+3 Holes, Space 8.0mm	1
114.06	2.0mm T-Plate with Bridge, 2+3 Holes, Space 15.0mm	1
114.07	2.0mm T-Plate with Bridge, 2+6 Holes, Space 6.0mm	1
116.01	2.0mm T-Plate, 3+3 Holes	1
116.02	2.0mm T-Plate, 3+4 Holes	1
116.03	2.0mm T-Plate, 3+6 Holes	1
134.01	2.0mm Rectangular Grill Plate, Type-I, 2 x 2 Holes	1
134.02	2.0mm Rectangular Grill Plate, Type-I, 2 x 3 Hole	1
135.01	2.0mm Rectangular Grill Plate, Type-II, 4 x 2 Holes	1
135.02	2.0mm Rectangular Grill Plate, Type-II, 4 x 4 Holes	1
135.03	2.0mm Rectangular Grill Plate, Type-II, 4 x 6 Holes	1

0.8mm	Item	Qty.
84.09	2.0mm Straight Plate, 4 Holes	1
84.10	2.0mm Straight Plate, 5 Holes	1
84.11	2.0mm Straight Plate, 6 Holes	1
84.12	2.0mm Straight Plate, 8 Holes	1
84.13	2.0mm Straight Plate, 1 Holes	1
84.14	2.0mm Straight Plate, 12 Holes	1
84.15	2.0mm Straight Plate, 14 Holes	1
84.16	2.0mm Straight Plate, 16 Holes	1
90.06	2.0mm Orbital Plate with Bridge, Type-I, 4 Holes, Space 8.0mm	1
90.07	2.0mm Orbital Plate with Bridge, Type-I, 4 Holes, Space 8.0mm	1
90.08	2.0mm Orbital Plate with Bridge, Type-I, 6 Holes, Space 8.0mm	1
90.09	2.0mm Orbital Plate with Bridge, Type-I, 6 Holes, Space 8.0mm	1
90.10	2.0mm Orbital Plate with Bridge, Type-I, 6 Holes, Space 8.0mm	1
92.04	2.0mm Orbital Plate, Space 6.0mm	1
92.05	2.0mm Orbital Plate, Space 8.0mm	1
92.06	2.0mm Orbital Plate, Space 12.0mm	1
93.05	2.0mm L-Plate, Type-II, 90°, 2+2 Holes - Left	1
93.07	2.0mm L-Plate, Type-II, 90°, 2+4 Holes - Left	1
93.06	2.0mm L-Plate, Type-II, 90°, 2+2 Holes - Right	1
93.08	2.0mm L-Plate, Type-II, 90°, 2+4 Holes - Right	1
94.15	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 8.0mm - Left	1
94.17	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 10.0mm - Left	1
94.19	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 12.0mm - Left	1
94.21	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 15.0mm - Left	1
94.23	2.0mm L-Plate with Bridge, Type-III, 90°, 3+3 Holes, Space 9.0mm - Left	1
94.25	2.0mm L-Plate with Bridge, Type-III, 90°, 3+3 Holes, Space 13.0mm - Left	1
94.27	2.0mm L-Plate with Bridge, Type-III, 90°, 2+3 Holes, Space 12.0mm - Left	1
94.16	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 8.0mm - Right	1
94.18	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 10.0mm - Right	1
94.20	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 12.0mm - Right	1
94.22	2.0mm L-Plate with Bridge, Type-III, 90°, 2+2 Holes, Space 15.0mm - Right	1
94.24	2.0mm L-Plate with Bridge, Type-III, 90°, 3+3 Holes, Space 9.0mm - Right	1
94.26	2.0mm L-Plate with Bridge, Type-III, 90°, 3+3 Holes, Space 13.0mm - Right	1
94.28	2.0mm L-Plate with Bridge, Type-III, 90°, 2+3 Holes, Space 12.0mm - Right	1
98.05	2.0mm Double X-Plate, Small, 4+2 Holes, Space 6.0mm	1
98.06	2.0mm Double X-Plate, Small, 4+2 Holes, Space 8.0mm	1
98.07	2.0mm Double X-Plate, Small, 4+4 Holes, Space 10.0mm	1
98.08	2.0mm Double X-Plate, Small, 4+4 Holes, Space 12.0mm	1
112.04	2.0mm H-Plate, 6+1 Holes	1
112.05	2.0mm H-Plate, 8+1 Holes	1
112.06	2.0mm H-Plate, 10+1 Holes	1
113.03	2.0mm I-Plate, Space 6.0mm	1
113.04	2.0mm I-Plate, Space 9.0mm	1
114.08	2.0mm T-Plate with Bridge, 2+2 Holes, Space 6.0mm	1

0.8mm	Item	Qty.
114.09	2.0mm T-Plate with Bridge, 2+2 Holes, Space 8.0mm	1
114.10	2.0mm T-Plate with Bridge, 2+2 Holes, Space 12.0mm	1
114.11	2.0mm T-Plate with Bridge, 2+4 Holes, Space 6.0mm	1
114.12	2.0mm T-Plate with Bridge, 2+3 Holes, Space 8.0mm	1
114.13	2.0mm T-Plate with Bridge, 2+3 Holes, Space 15.0mm	1
114.14	2.0mm T-Plate with Bridge, 2+6 Holes, Space 6.0mm	1
116.04	2.0mm T-Plate, 3+3 Holes	1
116.05	2.0mm T-Plate, 3+4 Holes	1
116.06	2.0mm T-Plate, 3+6 Holes	1
134.03	2.0mm Rectangular Grill Plate, Type-I, 2 x 2 Holes	1
134.04	2.0mm Rectangular Grill Plate, Type-I, 2 x 3 Holes	1
135.04	2.0mm Rectangular Grill Plate, Type-II, 4 x 2 Holes	1
135.05	2.0mm Rectangular Grill Plate, Type-II, 4 x 4 Holes	1
135.06	2.0mm Rectangular Grill Plate, Type-II, 4 x 6 Holes	1
143.01	2.0mm Locking Straight Plate with Bridge, Type-I, 4 Holes	1
143.02	2.0mm Locking Straight Plate with Bridge, Type-I, 6 Holes	1
143.03	2.0mm Locking Straight Plate with Bridge, Type-I, 8 Holes	1
143.04	2.0mm Locking Straight Plate with Bridge, Type-I, 12 Holes	1
143.05	2.0mm Locking Straight Plate with Bridge, Type-I, 14 Holes	1
144.01	2.0mm Locking T-Gap Plate, Straight, 4 Holes, Space 10.0mm	1
144.02	2.0mm Locking T-Gap Plate, Straight, 4 Holes, Space 12.0mm	1
144.03	2.0mm Locking T-Gap Plate, Straight, 6 Holes, Space 10.0mm	1
144.04	2.0mm Locking T-Gap Plate, Straight, 6 Holes, Space 12.0mm	1
145.01	2.0mm Locking U-Plate, 3+3 Holes	1
145.02	2.0mm Locking U-Plate, 4+4 Holes	1
146.01	2.0mm Locking U-Plate, 9 Holes	1
148.01	2.0mm Locking Rectangular Grill Plate, Type-I, 4 x 2 Holes	1
148.02	2.0mm Locking Rectangular Grill Plate, Type-I, 4 x 6 Holes	1

1.2mm	Item	Qty.
150.01	2.0mm Locking Straight Plate with Bridge, Type-II, 4 Holes	1
151.01	2.0mm Locking Rectangular Grill Plate, Type-III, 2 x 3 Holes	1
151.02	2.0mm Locking Rectangular Grill Plate, Type-III, 2 x 4 Holes	1
151.03	2.0mm Locking Rectangular Grill Plate, Type-III, 2 x 5 Holes	1
151.04	2.0mm Locking Rectangular Grill Plate, Type-III, 2 x 6 Holes	1
152.01	2.0mm Locking Straight Plate with Bridge, Type-III, 4 Holes, Space 10.0mm	1
152.02	2.0mm Locking Straight Plate with Bridge, Type-III, 4 Holes, Space 12.0mm	1
152.03	2.0mm Locking Straight Plate with Bridge, Type-III, 6 Holes, Space 10.0mm	1
152.04	2.0mm Locking Straight Plate with Bridge, Type-III, 6 Holes, Space 12.0mm	1
153.01	2.0mm Locking Straight Plate with Bridge, Type-IV, 4 Holes, Space 8.0mm	1
153.02	2.0mm Locking Straight Plate with Bridge, Type-IV, 4 Holes, Space 10.0mm	1
153.03	2.0mm Locking Straight Plate with Bridge, Type-IV, 4 Holes, Space 12.0mm	1
153.04	2.0mm Locking Straight Plate with Bridge, Type-IV, 6 Holes, Space 8.0mm	1
153.05	2.0mm Locking Straight Plate with Bridge, Type-IV, 6 Holes, Space 10.0mm	1
153.06	2.0mm Locking Straight Plate with Bridge, Type-IV, 6 Holes, Space 12.0mm	1
154.01	2.0mm Locking Straight Plate, Type-II, 4 Holes	1
154.02	2.0mm Locking Straight Plate, Type-II, 6 Holes	1
154.03	2.0mm Locking Straight Plate, Type-II, 8 Holes	1
154.04	2.0mm Locking Straight Plate, Type-II, 10 Holes	1
154.05	2.0mm Locking Straight Plate, Type-II, 12 Holes	1
154.06	2.0mm Locking Straight Plate, Type-II, 14 Holes	1
154.07	2.0mm Locking Straight Plate, Type-II, 16 Holes	1

1.5mm	Item	Qty.
155.01	2.0mm Locking Straight Plate with Bridge, Type-V, 4 Holes, Space 8.0mm	1
155.02	2.0mm Locking Straight Plate with Bridge, Type-V, 4 Holes, Space 10.0mm	1
155.03	2.0mm Locking Straight Plate with Bridge, Type-V, 4 Holes, Space 12.0mm	1
155.04	2.0mm Locking Straight Plate with Bridge, Type-V, 6 Holes, Space 8.0mm	1
156.01	2.0mm Locking Straight Plate, Type-III, 6 Holes	1
156.02	2.0mm Locking Straight Plate, Type-III, 8 Holes	1
156.03	2.0mm Locking Straight Plate, Type-III, 26 Holes	1
157.01	2.0mm Locking Orbital Plate, 4 Holes, Space 10.0mm	1
157.02	2.0mm Locking Orbital Plate, 4 Holes, Space 15.0mm	1
159.01	2.0mm One -Way Pressure Plate with Bridge, 4 Holes, Space 8.0mm	1
159.02	2.0mm One -Way Pressure Plate with Bridge, 4 Holes, Space 10.0mm	1
159.03	2.0mm One -Way Pressure Plate with Bridge, 4 Holes, Space 12.0mm	1
159.04	2.0mm One -Way Pressure Plate with Bridge, 6 Holes, Space 8.0mm	1
159.05	2.0mm One -Way Pressure Plate with Bridge, 6 Holes, Space 10.0mm	1
160.01	2.0mm One -Way Pressure Plate, 6 Holes	1
160.02	2.0mm One -Way Pressure Plate, 8 Holes	1
160.03	2.0mm One -Way Pressure Plate, 26 Holes	1
161.01	2.0mm One -Way Pressure Orbital Plate 4 Holes, Space 10.0mm	1
161.02	2.0mm One -Way Pressure Orbital Plate 4 Holes, Space 15.0mm	1
162.01	2.0mm One -Way Pressure Curved Plate, 4 Holes, Space 6.0mm - Left	1
162.03	2.0mm One -Way Pressure Curved Plate, 4 Holes, Space 10.5mm - Left	1
162.05	2.0mm One -Way Pressure Curved Plate, 6 Holes, Space 10.5mm - Left	1
162.07	2.0mm One -Way Pressure Curved Plate, 8 Holes, Space 10.5mm - Left	1
162.02	2.0mm One -Way Pressure Curved Plate, 4 Holes, Space 6.0mm - Right	1
162.04	2.0mm One -Way Pressure Curved Plate, 4 Holes, Space 10.5mm - Right	1
162.06	2.0mm One -Way Pressure Curved Plate, 6 Holes, Space 10.5mm - Right	1
162.08	2.0mm One -Way Pressure Curved Plate, 8 Holes, Space 10.5mm - Right	1

Code	2.0mm Screw System	Qty.
163.04	2.0mm Self Tapping Screws, Length 4mm	13
163.05	2.0mm Self Tapping Screws, Length 5mm	13
163.07	2.0mm Self Tapping Screws, Length 7mm	13
163.09	2.0mm Self Tapping Screws, Length 9mm	13
163.11	2.0mm Self Tapping Screws, Length 11mm	13
163.13	2.0mm Self Tapping Screws, Length 13mm	13
163.15	2.0mm Self Tapping Screws, Length 15mm	13
163.17	2.0mm Self Tapping Screws, Length 17mm	13
163.19	2.0mm Self Tapping Screws, Length 19mm	13
164.09	2.0mm Self Tapping & Self Drilling Screws, Length 9mm	26
167.05	2.3mm Emergency Screws, Self Tapping, Length 5mm	13
167.07	2.3mm Emergency Screws, Self Tapping, Length 7mm	13
167.09	2.3mm Emergency Screws, Self Tapping, Length 9mm	13
167.11	2.3mm Emergency Screws, Self Tapping, Length 11mm	13
167.13	2.3mm Emergency Screws, Self Tapping, Length 13mm	13
167.15	2.3mm Emergency Screws, Self Tapping, Length 15mm	13
167.17	2.3mm Emergency Screws, Self Tapping, Length 17mm	13
167.19	2.3mm Emergency Screws, Self Tapping, Length 19mm	13

Bone Plate & Screw




1.0mm

1.5mm

2.0mm

2.4mm



Neuro
Facial
Mandible
Orthognathic


2.4mm Mandible Reconstruction Plate, Straight

Holes	Codes (2.0mm)
12	168.01
15	168.02
17	168.03
20	168.04

2.4mm Locking Two-way Compression Mandible Reconstruction Plate, Straight

Holes	Codes (2.0mm)
12	170.01
15	170.02
17	170.03
20	170.04

2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.0mm)
5+10	Left	171.01L
5+10	Right	171.01R
5+15	Left	171.02L
5+15	Right	171.02R
5+20	Left	171.03L
5+20	Right	171.03R

2.4mm Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.0mm)
5+10	Left	169.01L
5+10	Right	169.01R
5+15	Left	169.02L
5+15	Right	169.02R
5+20	Left	169.03L
5+20	Right	169.03R

2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight

Holes	Codes (2.0mm)
12	172.01
15	172.02
17	172.03
20	172.04

2.4mm Locking Two-Way Compression Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.0mm)
5+10	Left	173.01L
5+10	Right	173.01R
5+15	Left	173.02L
5+15	Right	173.02R
5+20	Left	173.03L
5+20	Right	173.03R

STERILE
2.4mm Mandible Reconstruction Plate, Straight

Holes	Codes (2.0mm)
12	168.01-S
15	168.02-S
17	168.03-S
20	168.04-S

2.4mm Locking Two-way Compression Mandible Reconstruction Plate, Straight

Holes	Codes (2.0mm)
12	170.01-S
15	170.02-S
17	170.03-S
20	170.04-S

2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.0mm)
5+10	Left	171.01L-S
5+10	Right	171.01R-S
5+15	Left	171.02L-S
5+15	Right	171.02R-S
5+20	Left	171.03L-S
5+20	Right	171.03R-S

2.4mm Mandible Reconstruction Plate, Angled


Holes	Direction	Codes (2.0mm)
5+10	Left	169.01L-S
5+10	Right	169.01R-S
5+15	Left	169.02L-S
5+15	Right	169.02R-S
5+20	Left	169.03L-S
5+20	Right	169.03R-S

2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight

Holes	Codes (2.0mm)
12	172.01-S
15	172.02-S
17	172.03-S
20	172.04-S


2.4mm Locking Two-Way Compression Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.0mm)
5+10	Left	173.01L-S
5+10	Right	173.01R-S
5+15	Left	173.02L-S
5+15	Right	173.02R-S
5+20	Left	173.03L-S
5+20	Right	173.03R-S




2.4mm Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.0mm)
5+25+5	174.01




2.4mm Locking Round Holes Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.0mm)
5+25+5	175.01



2.4mm Locking Two-Way Compression Mandible Reconstruction Plate, Double Angled


Holes	Codes (2.0mm)
5+25+5	176.01



STERILE

2.4mm Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.0mm)
5+25+5	174.01-S



STERILE

2.4mm Locking Round Holes Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.0mm)
5+25+5	175.01-S



STERILE

2.4mm Locking Two-Way Compression Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.0mm)
5+25+5	176.01-S



2.4mm Mandible Reconstruction Plate, Straight

Holes	Codes (2.5mm)
12	177.01
15	177.02
17	177.03
20	177.04

2.4mm Locking Two-Way Mandible Reconstruction Plate, Straight

Holes	Codes (2.5mm)
12	178.01
15	178.02
17	178.03
20	178.04

2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.5mm)
5+10	Left	179.01L
5+10	Right	179.01R
5+15	Left	179.02L
5+15	Right	179.02R
5+20	Left	179.03L
5+20	Right	179.03R

2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight

Holes	Codes (2.5mm)
12	180.01
15	180.02
17	180.03
20	180.04

2.4mm Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.5mm)
5+10	Left	181.01L
5+10	Right	181.01R
5+15	Left	181.02L
5+15	Right	181.02R
5+20	Left	181.03L
5+20	Right	181.03R

2.4mm Locking Two-Way Compression Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.5mm)
5+10	Left	182.01L
5+10	Right	182.01R
5+15	Left	182.02L
5+15	Right	182.02R
5+20	Left	182.03L
5+20	Right	182.03R

STERILE
2.4mm Mandible Reconstruction Plate, Straight

Holes	Codes (2.5mm)
12	177.01-S
15	177.02-S
17	177.03-S
20	177.04-S

2.4mm Locking Two-Way Mandible Reconstruction Plate, Straight

Holes	Codes (2.5mm)
12	178.01-S
15	178.02-S
17	178.03-S
20	178.04-S

2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.5mm)
5+10	Left	179.01L-S
5+10	Right	179.01R-S
5+15	Left	179.02L-S
5+15	Right	179.02R-S
5+20	Left	179.03L-S
5+20	Right	179.03R-S

2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight


Holes	Codes (2.5mm)
12	180.01-S
15	180.02-S
17	180.03-S
20	180.04-S

2.4mm Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.5mm)
5+10	Left	181.01L-S
5+10	Right	181.01R-S
5+15	Left	181.02L-S
5+15	Right	181.02R-S
5+20	Left	181.03L-S
5+20	Right	181.03R-S


2.4mm Locking Two-Way Compression Mandible Reconstruction Plate, Angled

Holes	Direction	Codes (2.5mm)
5+10	Left	182.01L-S
5+10	Right	182.01R-S
5+15	Left	182.02L-S
5+15	Right	182.02R-S
5+20	Left	182.03L-S
5+20	Right	182.03R-S



2.4mm Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.5mm)
5+25+5	183.01




2.4mm Locking Round Holes Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.5mm)
5+25+5	184.01



2.4mm Locking Two-Way Compression Mandible Reconstruction Plate, Double Angled


Holes	Codes (2.5mm)
5+25+5	185.01



STERILE

2.4mm Mandible Reconstruction Plate, Double Angled


Holes	Codes (2.5mm)
5+25+5	183.01-S



STERILE

2.4mm Locking Round Holes Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.5mm)
5+25+5	184.01-S



STERILE

2.4mm Locking Two-Way Compression Mandible Reconstruction Plate, Double Angled

Holes	Codes (2.5mm)
5+25+5	185.01-S

2.4mm Self Tapping Screws


Length	Codes
7mm	186.07
9mm	186.09
11mm	186.11
13mm	186.13
15mm	186.15
17mm	186.17
19mm	186.19

STERILE

Length	Codes
7mm	186.07-S
9mm	186.09-S
11mm	186.11-S
13mm	186.13-S
15mm	186.15-S
17mm	186.17-S
19mm	186.19-S


2.4mm Locking Self Tapping Screws


Length	Codes
7mm	187.07
9mm	187.09
10mm	187.10
12mm	187.12
13mm	187.13
14mm	187.14
16mm	187.16
18mm	187.18
20mm	187.20

STERILE

Length	Codes
7mm	187.07-S
9mm	187.09-S
10mm	187.10-S
12mm	187.12-S
13mm	187.13-S
14mm	187.14-S
16mm	187.16-S
18mm	187.18-S
20mm	187.20-S



2.7mm Emergency Screws, Self Tapping


Length	Codes
7mm	188.07
9mm	188.09
10mm	188.10
12mm	186.12
13mm	188.13
14mm	188.14
16mm	188.16
18mm	188.18
20mm	188.20

STERILE

Length	Codes
7mm	188.07-S
9mm	188.09-S
10mm	188.10-S
12mm	186.12-S
13mm	188.13-S
14mm	188.14-S
16mm	188.16-S
18mm	188.18-S
20mm	188.20-S




7-134-46 Implant and Instrument Set for 2.4mm CMF System

Code	Item	Qty.
168.01	2.4mm Mandible Reconstruction Plate, Straight, 12 Holes	1
168.02	2.4mm Mandible Reconstruction Plate, Straight, 15 Holes	1
168.03	2.4mm Mandible Reconstruction Plate, Straight, 17 Holes	1
168.04	2.4mm Mandible Reconstruction Plate, Straight, 20 Holes	1
169.01L	2.4mm Mandible Reconstruction Plate, Angled, 5+10 Holes - Left	1
169.02L	2.4mm Mandible Reconstruction Plate, Angled, 5+15 Holes - Left	1
169.03L	2.4mm Mandible Reconstruction Plate, Angled, 5+20 Holes - Left	1
169.01R	2.4mm Mandible Reconstruction Plate, Angled, 5+10 Holes -Right	1
169.02R	2.4mm Mandible Reconstruction Plate, Angled, 5+15 Holes -Right	1
169.03R	2.4mm Mandible Reconstruction Plate, Angled, 5+20 Holes -Right	1
171.01L	2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled, 5+10 Holes - Left	1
171.02L	2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled, 5+15 Holes - Left	1
171.03L	2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled, 5+20 Holes - Left	1
171.01R	2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled, 5+10 Holes - Right	1
171.02R	2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled, 5+15 Holes - Right	1
171.03R	2.4mm Locking Round Holes Mandible Reconstruction Plate, Angled, 5+20 Holes - Right	1
172.01	2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight, 12 Holes	1
172.02	2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight, 15 Holes	1
172.03	2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight, 17 Holes	1

172.04	2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight, 20 Holes	1
174.01	2.4mm Mandible Reconstruction Plate, Double Angled, 5+25+5 Holes	1
175.01	2.4mm Locking Round Holes Mandible Reconstruction Plate, Double Angled, 5+25+5 Holes	1
177.01	2.4mm Mandible Reconstruction Plate, Straight, 12 Holes	1
177.02	2.4mm Mandible Reconstruction Plate, Straight, 15 Holes	1
177.03	2.4mm Mandible Reconstruction Plate, Straight, 17 Holes	1
177.04	2.4mm Mandible Reconstruction Plate, Straight, 20 Holes	1
180.01	2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight, 12 Holes	1
180.02	2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight, 15 Holes	1
180.03	2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight, 17 Holes	1
180.04	2.4mm Locking Round Holes Mandible Reconstruction Plate, Straight, 20 Holes	1

Code	2.4mm Screw System	Qty.
186.07	2.4mm Self Tapping Screws, Length 7mm	18
186.09	2.4mm Self Tapping Screws, Length 9mm	18
186.11	2.4mm Self Tapping Screws, Length 11mm	18
186.13	2.4mm Self Tapping Screws, Length 13mm	18
186.15	2.4mm Self Tapping Screws, Length 15mm	18
186.17	2.4mm Self Tapping Screws, Length 17mm	18
186.19	2.4mm Self Tapping Screws, Length 19mm	18
188.07	2.7mm Emergency Screws, Self Tapping, Length 7mm	9
188.09	2.7mm Emergency Screws, Self Tapping, Length 9mm	9
188.10	2.7mm Emergency Screws, Self Tapping, Length 10mm	9
188.12	2.7mm Emergency Screws, Self Tapping, Length 12mm	9
188.13	2.7mm Emergency Screws, Self Tapping, Length 13mm	9
188.14	2.7mm Emergency Screws, Self Tapping, Length 14mm	9
188.16	2.7mm Emergency Screws, Self Tapping, Length 16mm	9
188.18	2.7mm Emergency Screws, Self Tapping, Length 18mm	9
188.20	2.7mm Emergency Screws, Self Tapping, Length 20mm	9

Bone Screw



Lag screw

Kirschner Screw

IMF Screw

Anchor Screw



Neuro
Facial
Mandible
Orthognathic

Screw System	Matched Drill Bit	Pin Emergency Screw	Name	Indication
1.5mm 2.0mm 2.4mm	1.2mm 1.7mm 2.0mm	None	Kirschner Wire, Lag Screw	Applicable for fracture of ankle bump, fracture fixation in hard fixation cases by bone plates.
1.5mm 2.0mm 2.4mm	Self Tapping & Drilling Screw	None	Intermaxillary Fixation Screw, Anchorage Screw, Traction Screw	Keep normal occlusion relationship, applicable for fracture of maxillary, mandible and orthognathic surgery. Features: Screw self-tapping & self-drilling directly, ligating by titanium wire.
1.3mm 1.9mm 2.4mm	1.3mm 1.9mm 2.4mm	None	Anchor Screw	Applicable for connection between bone and tendon in head of joint because of injury, diseases, sport injury.

Anchor Screws



Length	Codes
Ø1.3mmx3.5mm	200.01
Ø1.3mmx4.5mm	200.02
Ø1.9mmx4.5mm	200.03
Ø1.9mmx5.5mm	200.04
Ø2.4mmx8.5mm	200.05
Ø2.4mmx9.5mm	200.06

STERILE

Length	Codes
Ø1.3mmx3.5mm	200.01-S
Ø1.3mmx4.5mm	200.02-S
Ø1.9mmx4.5mm	200.03-S
Ø1.9mmx5.5mm	200.04-S
Ø2.4mmx8.5mm	200.05-S
Ø2.4mmx9.5mm	200.06-S

1.5mm Kirschner Screws


Length	Codes
18mm	189.18
20mm	189.20
22mm	189.22
24mm	189.24

STERILE

Length	Codes
18mm	189.18-S
20mm	189.20-S
22mm	189.22-S
24mm	189.24-S


1.9mm Emergency Screws


Length	Codes
18mm	190.18
20mm	190.20
22mm	190.22
24mm	190.24

STERILE

Length	Codes
18mm	190.18-S
20mm	190.20-S
22mm	190.22-S
24mm	190.24-S


2.0mm Kirschner Screws


Length	Codes
18mm	191.18
20mm	191.20
22mm	191.22
24mm	191.24

STERILE

Length	Codes
18mm	191.18-S
20mm	191.20-S
22mm	191.22-S
24mm	191.24-S


2.3mm Emergency Screws


Length	Codes
18mm	192.18
20mm	192.20
22mm	192.22
24mm	192.24

STERILE

Length	Codes
18mm	192.18-S
20mm	192.20-S
22mm	192.22-S
24mm	192.24-S


2.3mm Emergency Screws


Length	Codes
18mm	193.18
20mm	193.20
22mm	193.22
24mm	193.24

STERILE

Length	Codes
18mm	193.18-S
20mm	193.20-S
22mm	193.22-S
24mm	193.24-S



1.5mm Self Tapping & Self Drilling IMF Screws, Type-I

STERILE

Length	Codes
5mm	194.05
7mm	194.07
9mm	194.09
11mm	194.11
13mm	194.13

Length	Codes
5mm	194.05-S
7mm	194.07-S
9mm	194.09-S
11mm	194.11-S
13mm	194.13-S


1.5mm Self Tapping & Self Drilling IMF Screws, Type-II

STERILE

Length	Codes
5mm	195.05
7mm	195.07
9mm	195.09
11mm	195.11
13mm	195.13

Length	Codes
5mm	195.05-S
7mm	195.07-S
9mm	195.09-S
11mm	195.11-S
13mm	195.13-S


1.5mm Self Tapping & Self Drilling IMF Screws, Type-III

STERILE

Length	Codes
5mm	196.05
7mm	196.07
9mm	196.09
11mm	196.11
13mm	196.13

Length	Codes
5mm	196.05-S
7mm	196.07-S
9mm	196.09-S
11mm	196.11-S
13mm	196.13-S



2.0mm Self Tapping & Self Drilling IMF Screws, Type-I

STERILE

Length	Codes
5mm	197.05
7mm	197.07
9mm	197.09
11mm	197.11
13mm	197.13

Length	Codes
5mm	197.05-S
7mm	197.07-S
9mm	197.09-S
11mm	197.11-S
13mm	197.13-S


2.0mm Self Tapping & Self Drilling IMF Screws, Type-II

STERILE

Length	Codes
5mm	198.05
7mm	198.07
9mm	198.09
11mm	198.11
13mm	198.13

Length	Codes
5mm	198.05-S
7mm	198.07-S
9mm	198.09-S
11mm	198.11-S
13mm	198.13-S


2.0mm Self Tapping & Self Drilling IMF Screws, Type-III

STERILE

Length	Codes
5mm	199.05
7mm	199.07
9mm	199.09
11mm	199.11
13mm	199.13

Length	Codes
5mm	199.05-S
7mm	199.07-S
9mm	199.09-S
11mm	199.11-S
13mm	199.13-S





Instrumentation System

1.0mm
1.5mm
2.0mm
2.4mm

Neuro
Facial
Mandible
Orthognathic

Drill Bit with Stop Series



Stryker Interface	General Interface	Description
CS50.0803	CG50.0803	0.8mm x 50 x 3mm
CS50.0805	CG50.0805	0.8mm x 50 x 5mm
CS50.0807	CG50.0807	0.8mm x 50 x 7mm
CS50.0815	CG50.0815	0.8mm x 50 x 15mm
CS50.1205	CG50.1205	1.2mm x 50 x 5mm
CS50.1207	CG50.1207	1.2mm x 50 x 7mm
CS50.1209	CG50.1209	1.2mm x 50 x 9mm
CS50.1211	CG50.1211	1.2mm x 50 x 11mm
CS50.1213	CG50.1213	1.2mm x 50 x 13mm
CS50.1215	CG50.1215	1.2mm x 50 x 15mm
CS50.1705	CG50.1705	1.7mm x 50 x 5mm
CS50.1707	CG50.1707	1.7mm x 50 x 7mm
CS50.1709	CG50.1709	1.7mm x 50 x 9mm
CS50.1711	CG50.1711	1.7mm x 50 x 11mm
CS50.1713	CG50.1713	1.7mm x 50 x 13mm
CS50.1715	CG50.1715	1.7mm x 50 x 15mm
CS50.1717	CG50.1717	1.7mm x 50 x 17mm
CS50.1719	CG50.1719	1.7mm x 50 x 19mm
CS50.1721	CG50.1721	1.7mm x 50 x 21mm
CS50.1723	CG50.1723	1.7mm x 50 x 23mm
CS50.2005	CG50.2005	2.0mm x 50 x 5mm
CS50.2007	CG50.2007	2.0mm x 50 x 7mm
CS50.2009	CG50.2009	2.0mm x 50 x 9mm
CS50.2011	CG50.2011	2.0mm x 50 x 11mm
CS50.2013	CG50.2013	2.0mm x 50 x 13mm
CS50.2015	CG50.2015	2.0mm x 50 x 15mm

Stryker Interface	General Interface	Description
CS105.0803	CG105.0803	0.8mm x 105 x 3mm
CS105.0805	CG105.0805	0.8mm x 105 x 5mm
CS105.0807	CG105.0807	0.8mm x 105 x 7mm
CS105.0815	CG105.0815	0.8mm x 105 x 15mm
CS105.1205	CG105.1205	1.2mm x 105 x 5mm
CS105.1207	CG105.1207	1.2mm x 105 x 7mm
CS105.1209	CG105.1209	1.2mm x 105 x 9mm
CS105.1211	CG105.1211	1.2mm x 105 x 11mm
CS105.1213	CG105.1213	1.2mm x 105 x 13mm
CS105.1215	CG105.1215	1.2mm x 105 x 15mm
CS105.1705	CG105.1705	1.7mm x 105 x 5mm
CS105.1707	CG105.1707	1.7mm x 105 x 7mm
CS105.1709	CG105.1709	1.7mm x 105 x 9mm
CS105.1711	CG105.1711	1.7mm x 105 x 11mm
CS105.1713	CG105.1713	1.7mm x 105 x 13mm
CS105.1715	CG105.1715	1.7mm x 105 x 15mm
CS105.1717	—	1.7mm x 105 x 17mm
CS105.1719	—	1.7mm x 105 x 19mm
CS105.1721	—	1.7mm x 105 x 21mm
CS105.1723	—	1.7mm x 105 x 23mm
CS105.2005	CG105.2005	2.0mm x 105 x 5mm
CS105.2007	CG105.2007	2.0mm x 105 x 7mm
CS105.2009	CG105.2009	2.0mm x 105 x 9mm
CS105.2011	CG105.2011	2.0mm x 105 x 11mm
CS105.2013	CG105.2013	2.0mm x 105 x 13mm
CS105.2015	CG105.2015	2.0mm x 105 x 15mm

7-134-12 Drill Guide Ø0.8mm for CMF System



7-134-13 Drill Bit, Ø0.8mm x Length 50mm for CMF System



7-134-14 Drill Bit, Ø0.8mm x Length 4mm for CMF System



7-134-15 Drill Bit, Ø0.8mm x Length 5mm for CMF System



7-134-16 Drill Bit, Ø0.8mm x Length 6mm for CMF System



**Implant and Instrument Set
1.0mm CMF System**

7-134-17 Drill Bit, Ø0.8mm x Length 7mm for CMF System




7-134-18 Self Retaining Crosshead Screwdriver Shaft for 1.0mm Screw



7-134-19 1.0mm Screw Caddy for CMF System



7-134-20 Container with Tray & Lid for 1.0mm CMF Set




7-134-11

Implant and Instrument Set for 1.0mm CMF System

Code	Item	Qty.
7-134-11	Implant and Instrument Set for 1.0mm CMF System	1
7-134-12	Drill Guide Ø0.8mm for CMF System	1
7-134-13	Drill Bit, Ø0.8mm x Length 50mm for CMF System	2
7-134-14	Drill Bit, Ø0.8mm x Length 4mm for CMF System	1
7-134-15	Drill Bit, Ø0.8mm x Length 5mm for CMF System	1
7-134-16	Drill Bit, Ø0.8mm x Length 6mm for CMF System	1
7-134-17	Drill Bit, Ø0.8mm x Length 7mm for CMF System	1
7-134-18	Self Retaining Crosshead Screwdriver Shaft for 1.0mm Screw	2
7-134-19	1.0mm Screw Caddy for CMF System	1
7-134-20	Container with Tray & Lid for 1.0mm CMF Set	1

7-134-22 Drill Guide Ø1.1mm for CMF System



7-134-23 Drill Bit, Ø1.1mm x Length 50mm for CMF System



7-134-24 Drill Bit, Ø1.1mm x Length 4mm for CMF System



7-134-25 Drill Bit, Ø1.1mm x Length 5mm for CMF System



7-134-26 Drill Bit, Ø1.1mm x Length 6mm for CMF System



Implant and Instrument Set 1.5mm CMF System

7-134-27 Drill Bit, Ø1.1mm x Length 7mm for CMF System



7-134-28 Self Retaining Crosshead Screwdriver Shaft for 1.5mm Screw




7-134-29 Self Retaining Crosshead Screwdriver Shaft for 1.3mm & 1.9mm Emergency Screw



7-134-30 1.5mm Screw Caddy for CMF System



7-134-31 Container with Tray & Lid for 1.5mm CMF Set




7-134-21 Implant and Instrument Set for 1.5mm CMF System

Code	Item	Qty.
7-134-21	Implant and Instrument Set for 1.5mm CMF System	1
7-134-22	Drill Guide Ø1.1mm for CMF System	1
7-134-23	Drill Bit, Ø1.1mm x Length 50mm for CMF System	2
7-134-24	Drill Bit, Ø1.1mm x Length 4mm for CMF System	1
7-134-25	Drill Bit, Ø1.1mm x Length 5mm for CMF System	1
7-134-26	Drill Bit, Ø1.1mm x Length 6mm for CMF System	1
7-134-27	Drill Bit, Ø1.1mm x Length 7mm for CMF System	1
7-134-28	Self Retaining Crosshead Screwdriver Shaft for 1.5mm Screw	1
7-134-29	Self Retaining Crosshead Screwdriver Shaft for 1.3mm & 1.9mm Emergency Screw	1
7-134-30	1.5mm Screw Caddy for CMF System	1
7-134-31	Container with Tray & Lid for 1.5mm CMF Set	1

7-134-33 Double Drill Guide Ø1.4/2.0mm for CMF System



7-134-34 Drill Bit, Ø1.4mm x Length 60mm for CMF System



7-134-35 Drill Sleeve, Ø1.4mm for CMF System



7-134-36 Drill Bit, Ø1.4mm x Length 5mm for CMF System



7-134-37 Drill Bit, Ø1.4mm x Length 7mm for CMF System



7-134-38 Drill Bit, Ø1.4mm x Length 9mm for CMF System



7-134-39 Drill Bit, Ø1.4mm x Length 11mm for CMF System



7-134-40 Drill Bit, Ø1.4mm x Length 13mm for CMF System



7-134-41 Tap, Ø2.0mm for CMF System



7-134-42 Self Retaining Crosshead Screwdriver Shaft for 2.0mm Screw



Implant and Instrument Set 2.0mm CMF System

7-134-43


Self Retaining Crosshead Screwdriver Shaft for 2.3mm Emergency Screw



7-134-44

2.0mm Screw Caddy for CMF System


7-134-45

Container with Tray & Lid for 2.0mm CMF Set




7-134-32 Implant and Instrument Set for 2.0mm CMF System

Code	Item	Qty.
7-134-32	Implant and Instrument Set for 2.0mm CMF System	1
7-134-33	Double Drill Guide Ø1.4/2.0mm for CMF System	1
7-134-34	Drill Bit, Ø1.4mm x Length 60mm for CMF System	2
7-134-35	Drill Sleeve, Ø1.4mm for CMF System	2
7-134-36	Drill Bit, Ø1.4mm x Length 5mm for CMF System	1
7-134-37	Drill Bit, Ø1.4mm x Length 7mm for CMF System	1
7-134-38	Drill Bit, Ø1.4mm x Length 9mm for CMF System	1
7-134-39	Drill Bit, Ø1.4mm x Length 11mm for CMF System	1
7-134-40	Drill Bit, Ø1.4mm x Length 13mm for CMF System	1
7-134-41	Tap, Ø2.0mm for CMF System	1
7-134-42	Self Retaining Crosshead Screwdriver Shaft for 2.0mm Screw	2
7-134-43	Self Retaining Crosshead Screwdriver Shaft for 2.3mm Emergency Screw	2
7-134-44	2.0mm Screw Caddy for CMF System	1
7-134-45	Container with Tray & Lid for 2.0mm CMF Set	1

7-134-47 Double Drill Guide Ø1.6/2.4mm for CMF System



7-134-48 Drill Bit, Ø1.6mm x Length 60mm for CMF System



7-134-49 Drill Bit, Ø1.6mm x Length 70mm for CMF System



7-134-50 Drill Sleeve, Ø1.6mm for CMF System



7-134-51 Tap, Ø2.4mm for CMF System



7-134-52 Self Retaining Crosshead Screwdriver Shaft for 2.4mm Screw



7-134-53 Straight Template for CMF System




7-134-54 L-Type Template for CMF System




7-134-55 2.4mm Screw Caddy for CMF System



7-134-56 Container with Tray & Lid for 2.4mm CMF Set




7-134-46 Implant and Instrument Set for 2.4mm CMF System

Code	Item	Qty.
7-134-46	Implant and Instrument Set for 2.4mm CMF System	1
7-134-47	Double Drill Guide Ø1.6/2.4mm for CMF System	1
7-134-48	Drill Bit, Ø1.6mm x Length 60mm for CMF System	2
7-134-49	Drill Bit, Ø1.6mm x Length 70mm for CMF System	2
7-134-50	Drill Sleeve, Ø1.6mm for CMF System	2
7-134-51	Tap, Ø2.4mm for CMF System	1
7-134-52	Self Retaining Crosshead Screwdriver Shaft for 2.4mm Screw	2
7-134-53	Straight Template for CMF System	1
7-134-54	L-Type Template for CMF System	1
7-134-55	2.4mm Screw Caddy for CMF System	1
7-134-56	Container with Tray & Lid for 2.4mm CMF Set	1

7-134-01 Quick Coupling Handle for Screwdrivers (for CMF System)



7-134-02 Quick Coupling Handle for Drill bits (for CMF System)



7-134-03 Plate Holding Forceps for CMF System



7-134-04 Position Ball tip Awl for CMF System



7-134-05 Large Depth Gauge measuring upto 65mm for CMF System



7-134-06 Bender for Mandible Plates



7-134-07 Bending Pliers with nose for CMF System



7-134-08 Bender for Midface Plates



7-134-09 Cutting Pliers for CMF System




7-134-10 Small Depth Gauge measuring upto 65mm for CMF System




7-134-57

Implant Tray for Cranio Maxillofacial Instrument Set




7-134-58

Instrument Tray for Cranio Maxillofacial Instrument Set



7-134-59

Container for Cranio Maxillofacial Instrument Set



7-134 Cranio Maxillofacial Instrument Set



7-134 Cranio Maxillofacial Instrument Set

Code	Item	Qty.
7-134-01	Quick Coupling Handle for Screwdrivers (for CMF System)	1
7-134-02	Quick Coupling Handle for Drill bits (for CMF System)	1
7-134-03	Plate Holding Forceps for CMF System	1
7-134-04	Position Ball tip Awl for CMF System	1
7-134-05	Large Depth Gauge measuring upto 65mm for CMF System	1
7-134-06	Bender for Mandible Plates	2
7-134-07	Bending Pliers with nose for CMF System	1
7-134-08	Bender for Midface Plates	2
7-134-09	Cutting Pliers for CMF System	1
7-134-10	Small Depth Gauge measuring upto 65mm for CMF System	1
7-134-57	Implant Tray for Cranio Maxillofacial Instrument Set	1
7-134-58	Instrument Tray for Cranio Maxillofacial Instrument Set	1
7-134-59	Container for Cranio Maxillofacial Instrument Set	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com