

HERBERT SCREW

HERBERT SCREW

2.5mm | 3.0mm



about us

Auxin Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements



Orthopaedic Implants

INDICATIONS


For fixation of bone fractures or for bone reconstruction.

Examples Include:

- Fixation of Osteotomies for Hallux Valgus treatment (such as Scarf, Chevron, etc.).
- Mono or Bi-cortical Osteotomies in the Foot or Hand.
- Distal or Proximal Metatarsal or Metacarpal Osteotomies.
- Arthrodesis in Hand or Foot surgery.
- Fixation of small bone fragments, in long bones or small bones fractures.

1.0mm Guide wire

Screw	Compression
Dia 2.5mm	Bicortical compression
Dia 3.0mm	Bicortical compression



Self drill &
Self tapping flutes



Star Head



- Totally Intra-Osseous Implants.
- Dual threads.
- Material: Titanium Alloy

2.5mm Cannulated Herbert Screw, Short Thread, (Star Head)



STERILE

Code	Length
HS25-110	10mm
HS25-112	12mm
HS25-114	14mm
HS25-116	16mm
HS25-118	18mm
HS25-120	20mm
HS25-122	22mm
HS25-124	24mm
HS25-126	26mm
HS25-128	28mm
HS25-130	30mm
HS25-132	32mm
HS25-134	34mm
HS25-136	36mm
HS25-138	38mm
HS25-140	40mm

Code	Length
HS25-110-S	10mm
HS25-112-S	12mm
HS25-114-S	14mm
HS25-116-S	16mm
HS25-118-S	18mm
HS25-120-S	20mm
HS25-122-S	22mm
HS25-124-S	24mm
HS25-126-S	26mm
HS25-128-S	28mm
HS25-130-S	30mm
HS25-132-S	32mm
HS25-134-S	34mm
HS25-136-S	36mm
HS25-138-S	38mm
HS25-140-S	40mm



2.5mm Cannulated Herbert Screw, Long Thread, (Star Head)



STERILE

Code	Length
HS25-216	16mm
HS25-218	18mm
HS25-220	20mm
HS25-222	22mm
HS25-224	24mm
HS25-226	26mm
HS25-228	28mm
HS25-230	30mm
HS25-232	32mm
HS25-234	34mm
HS25-236	36mm
HS25-238	38mm
HS25-240	40mm

Code	Length
HS25-216-S	16mm
HS25-218-S	18mm
HS25-220-S	20mm
HS25-222-S	22mm
HS25-224-S	24mm
HS25-226-S	26mm
HS25-228-S	28mm
HS25-230-S	30mm
HS25-232-S	32mm
HS25-234-S	34mm
HS25-236-S	36mm
HS25-238-S	38mm
HS25-240-S	40mm



3.0mm Cannulated Herbert Screw, Short Thread, (Star Head)



STERILE

Code	Length
HS30-110	10mm
HS30-112	12mm
HS30-114	14mm
HS30-116	16mm
HS30-118	18mm
HS30-120	20mm
HS30-122	22mm
HS30-124	24mm
HS30-126	26mm
HS30-128	28mm
HS30-130	30mm
HS30-132	32mm
HS30-134	34mm
HS30-136	36mm
HS30-138	38mm
HS30-140	40mm

Code	Length
HS30-110-S	10mm
HS30-112-S	12mm
HS30-114-S	14mm
HS30-116-S	16mm
HS30-118-S	18mm
HS30-120-S	20mm
HS30-122-S	22mm
HS30-124-S	24mm
HS30-126-S	26mm
HS30-128-S	28mm
HS30-130-S	30mm
HS30-132-S	32mm
HS30-134-S	34mm
HS30-136-S	36mm
HS30-138-S	38mm
HS30-140-S	40mm



3.0mm Cannulated Herbert Screw, Long Thread, (Star Head)



STERILE

Code	Length
HS30-216	16mm
HS30-218	18mm
HS30-220	20mm
HS30-222	22mm
HS30-224	24mm
HS30-226	26mm
HS30-228	28mm
HS30-230	30mm
HS30-232	32mm
HS30-234	34mm
HS30-236	36mm
HS30-238	38mm
HS30-240	40mm

Code	Length
HS30-216-S	16mm
HS30-218-S	18mm
HS30-220-S	20mm
HS30-222-S	22mm
HS30-224-S	24mm
HS30-226-S	26mm
HS30-228-S	28mm
HS30-230-S	30mm
HS30-232-S	32mm
HS30-234-S	34mm
HS30-236-S	36mm
HS30-238-S	38mm
HS30-240-S	40mm



HS101-19 Screw Caddy for 2.5/3.0mm Cannulated Herbert Screw System



2.5mm Cannulated Herbert Screw, Short Thread, (Star Head)

Code	Length	Unit
HS25-110	10mm	2
HS25-112	12mm	2
HS25-114	14mm	2
HS25-116	16mm	2
HS25-118	18mm	2
HS25-120	20mm	2
HS25-122	22mm	2
HS25-124	24mm	2
HS25-126	26mm	2
HS25-128	28mm	2
HS25-130	30mm	2
HS25-132	32mm	2
HS25-134	34mm	2
HS25-136	36mm	2
HS25-138	38mm	2
HS25-140	40mm	2

2.5mm Cannulated Herbert Screw, Long Thread, (Star Head)

Code	Length	Unit
HS25-216	16mm	2
HS25-218	18mm	2
HS25-220	20mm	2
HS25-222	22mm	2
HS25-224	24mm	2
HS25-226	26mm	2
HS25-228	28mm	2
HS25-230	30mm	2
HS25-232	32mm	2
HS25-234	34mm	2
HS25-236	36mm	2
HS25-238	38mm	2
HS25-240	40mm	2

3.0mm Cannulated Herbert Screw, Short Thread, (Star Head)

Code	Length	Unit
HS30-110	10mm	2
HS30-112	12mm	2
HS30-114	14mm	2
HS30-116	16mm	2
HS30-118	18mm	2
HS30-120	20mm	2
HS30-122	22mm	2
HS30-124	24mm	2
HS30-126	26mm	2
HS30-128	28mm	2
HS30-130	30mm	2
HS30-132	32mm	2
HS30-134	34mm	2
HS30-136	36mm	2
HS30-138	38mm	2
HS30-140	40mm	2

3.0mm Cannulated Herbert Screw, Long Thread, (Star Head)

Code	Length	Unit
HS30-216	16mm	2
HS30-218	18mm	2
HS30-220	20mm	2
HS30-222	22mm	2
HS30-224	24mm	2
HS30-226	26mm	2
HS30-228	28mm	2
HS30-230	30mm	2
HS30-232	32mm	2
HS30-234	34mm	2
HS30-236	36mm	2
HS30-238	38mm	2
HS30-240	40mm	2

Herbert Screw Instruments

HS101-01 Screw Holding Forcep for Herbert Screw



HS101-02 Holding Sleeve for 3.0mm Herbert Screw



HS101-03 Holding Sleeve for 2.5mm Herbert Screw



HS101-04 Holding Sleeve Handle for Herbert Screw



HS101-05 Guide Sleeve, Ø4.9mm, for Herbert Screw



HS101-06 Double Drill Guide, Ø2.0/1.1mm, for Herbert Screw



HS101-07 Guide Wire, Ø0.8mm x Length 150mm



HS101-08 Guide Wire, Ø1.0mm x Length 150mm



HS101-09 Sharp Hook for Herbert Screw



HS101-10 Depth Gauge for Herbert Screw



HS101-11 Double Periosteal Elevator for Herbert Screw



HS101-12 Cannulated Screwdriver Shaft, T8, for Herbert Screw



HS101-13 Screwdriver Shaft, T8, for Herbert Screw



HS101-14 Cannulated Drill Bit with Quick Coupling End, Ø2.0mm x Length 150mm, for Herbert Screw



HS101-15 Cannulated Reamer Bit with Quick Coupling End, Ø4.8mm x Length 150mm, for Herbert Screw



HS101-16 Handle With Quick Coupling for Herbert Screw




HS101-18 Cleaning Stylet Ø1.0mm for Herbert Screw




HS101-19 Screw Caddy for 2.5/3.0mm Cannulated Herbert Screw System




HS101-21 Instrument Tray for Herbert Screw Instrument Set



HS101-20 Container for Herbert Screw Instrument Set



HS101 Herbert Screw Instrument Set



HS101 Herbert Screw Instrument Set

Code	Set Consisting of	Units
HS101-01	Screw Holding Forcep for Herbert Screw	1
HS101-02	Holding Sleeve for 3.0mm Herbert Screw	1
HS101-03	Holding Sleeve for 2.5mm Herbert Screw	1
HS101-04	Holding Sleeve Handle for Herbert Screw	1
HS101-05	Guide Sleeve, Ø4.9mm, for Herbert Screw	1
HS101-06	Double Drill Guide, Ø2.0/1.1mm, for Herbert Screw	1
HS101-07	Guide Wire, Ø0.8mm x Length 150mm	4
HS101-08	Guide Wire, Ø1.0mm x Length 150mm	4
HS101-09	Sharp Hook for Herbert Screw	1
HS101-10	Depth Gauge for Herbert Screw	1
HS101-11	Double Periosteal Elevator for Herbert Screw	1
HS101-12	Cannulated Screwdriver Shaft, T8, for Herbert Screw	1
HS101-13	Screwdriver Shaft, T8, for Herbert Screw	1
HS101-14	Cannulated Drill Bit with Quick Coupling End, Ø2.0mm x Length 150mm, for Herbert Screw	1
HS101-15	Cannulated Reamer Bit with Quick Coupling End, Ø4.8mm x Length 150mm, for Herbert Screw	1
HS101-16	Handle With Quick Coupling for Herbert Screw	1
HS101-18	Cleaning Stylet Ø1.0mm for Herbert Screw	1
HS101-19	Screw Caddy for 2.5/3.0mm Cannulated Herbert Screw System	1
HS101-21	Instrument Tray for Herbert Screw Instrument Set	1
HS101-20	Container for Herbert Screw Instrument Set	1



USA Auxein Inc. 1500 Nw 89th Court, Suite 107-108 Doral, Florida 33172 Tel: +1 305 395 6062 E Fax: +1 305 395 6262 Email: USoffice@auxein.com	MEXICO Auxein México, S.A. de C.V. Tepic 139 int 801, Colonia Roma Sur, Alcaldía Cuauhtémoc, CDMX, México, C.P. 06760 Tel: +521 55 7261 0318 Email: info@auxein.mx	INDIA Auxein Medical Pvt. Ltd. Plot No. 168-169-170, Phase-4, Kundli Industrial Area, HSIIDC, Sector-57, Sonapat - 131028, Haryana Tel: +91 99106 43638 Fax: +91 86077 70197 Email: info@auxein.com
--	---	--