

HAND SYSTEM

# LOW Profile HAND SYSTEM

0.8mm | 1.3mm





# about us

Auxin Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

## Our Achievements



# **Orthopaedic Implants**

## Product Overview

Auxein is a global leader of innovative orthopaedic and medical solutions. We are dedicated to developing products, service methods, and approaches that improve patient care.



## The 0.8mm | 1.3mm Low Profile Hand System

The 0.8mm | 1.3mm Low Profile Hand System is designed to provide both standard and fracture-specific fixation for metacarpal and phalangeal fractures, as well as fixation for fusions and osteotomies. This comprehensive system contains plates for fractures of the metacarpal neck, fractures of the base of the first metacarpal, avulsion fractures, and rotational malunions. Additionally, the system contains standard-shaped, cut-to-length, and bend-to-fit plates and Low profile Screw & Wise-Lock Low profile Screw, Star Head for less complicated fractures.

Low-profile plates and screws and a rounded-edge plate cutter are designed to minimize soft tissue irritation. Versatile screws, customizable plates, and dedicated instrumentation offer a comprehensive system to streamline the surgical experience.

### Indications for Use:

The 0.8mm | 1.3mm Low Profile Hand System is designed for the management of fractures, fusions, and osteotomies of the distal, middle, and proximal phalanges and metacarpals and other bones of appropriate size for the devices.

**0.8mm Wise-Lock Low Profile Avulsion Fracture Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-001-02TI-S</b>	<b>12-001-02TI</b>	2	Titanium



**0.8mm Wise-Lock Low Profile Straight Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-002-03TI-S</b>	<b>12-002-03TI</b>	3	Titanium
<b>12-002-04TI-S</b>	<b>12-002-04TI</b>	4	Titanium
<b>12-002-05TI-S</b>	<b>12-002-05TI</b>	5	Titanium
<b>12-002-06TI-S</b>	<b>12-002-06TI</b>	6	Titanium
<b>12-002-07TI-S</b>	<b>12-002-07TI</b>	7	Titanium
<b>12-002-10TI-S</b>	<b>12-002-10TI</b>	10	Titanium



**0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-003-09TI-S</b>	<b>12-003-09TI</b>	9	Titanium



**0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate, Short**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-004-04LTI-S</b>	<b>12-004-04LTI</b>	4	Titanium
<b>12-004-04RTI-S</b>	<b>12-004-04RTI</b>	4	Titanium
<b>12-004-05LTI-S</b>	<b>12-004-05LTI</b>	5	Titanium
<b>12-004-05RTI-S</b>	<b>12-004-05RTI</b>	5	Titanium
<b>12-004-06LTI-S</b>	<b>12-004-06LTI</b>	6	Titanium
<b>12-004-06RTI-S</b>	<b>12-004-06RTI</b>	6	Titanium



**0.8mm Wise-Lock Low Profile Compression Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-005-04TI-S</b>	<b>12-005-04TI</b>	4	Titanium
<b>12-005-05TI-S</b>	<b>12-005-05TI</b>	5	Titanium
<b>12-005-06TI-S</b>	<b>12-005-06TI</b>	6	Titanium



**0.8mm Wise-Lock Low Profile T-Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-006-05TI-S</b>	<b>12-006-05TI</b>	5	Titanium
<b>12-006-06TI-S</b>	<b>12-006-06TI</b>	6	Titanium
<b>12-006-08TI-S</b>	<b>12-006-08TI</b>	8	Titanium
<b>12-006-012TI-S</b>	<b>12-006-012TI</b>	12	Titanium



**0.8mm Wise-Lock Low Profile L-Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-007-04LTI-S</b>	<b>12-007-04LTI</b>	4	Titanium
<b>12-007-04RTI-S</b>	<b>12-007-04RTI</b>	4	Titanium
<b>12-007-06LTI-S</b>	<b>12-007-06LTI</b>	6	Titanium
<b>12-007-06RTI-S</b>	<b>12-007-06RTI</b>	6	Titanium
<b>12-007-08LTI-S</b>	<b>12-007-08LTI</b>	8	Titanium
<b>12-007-08RTI-S</b>	<b>12-007-08RTI</b>	8	Titanium



**0.8mm Wise-Lock Low Profile Offset Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-008-04TI-S</b>	<b>12-008-04TI</b>	4	Titanium
<b>12-008-06TI-S</b>	<b>12-008-06TI</b>	6	Titanium
<b>12-008-10TI-S</b>	<b>12-008-10TI</b>	10	Titanium




1.5mm Wise-Lock Low Profile Screw, Star Head



**STERILE**

Code	Length (mm)	
12-015-05TI	5	Titanium
12-015-06TI	6	Titanium
12-015-07TI	7	Titanium
12-015-08TI	8	Titanium
12-015-09TI	9	Titanium
12-015-10TI	10	Titanium
12-015-11TI	11	Titanium
12-015-12TI	12	Titanium
12-015-13TI	13	Titanium
12-015-14TI	14	Titanium
12-015-16TI	16	Titanium
12-015-18TI	18	Titanium
12-015-20TI	20	Titanium

Code	Length (mm)	
12-015-05TI-S	5	Titanium
12-015-06TI-S	6	Titanium
12-015-07TI-S	7	Titanium
12-015-08TI-S	8	Titanium
12-015-09TI-S	9	Titanium
12-015-10TI-S	10	Titanium
12-015-11TI-S	11	Titanium
12-015-12TI-S	12	Titanium
12-015-13TI-S	13	Titanium
12-015-14TI-S	14	Titanium
12-015-16TI-S	16	Titanium
12-015-18TI-S	18	Titanium
12-015-20TI-S	20	Titanium



1.5mm Low Profile Screw, Star Head



**STERILE**

Code	Length (mm)	
12-016-05TI	5	Titanium
12-016-06TI	6	Titanium
12-016-07TI	7	Titanium
12-016-08TI	8	Titanium
12-016-09TI	9	Titanium
12-016-10TI	10	Titanium
12-016-11TI	11	Titanium
12-016-12TI	12	Titanium
12-016-13TI	13	Titanium
12-016-14TI	14	Titanium
12-016-16TI	16	Titanium
12-016-18TI	18	Titanium
12-016-20TI	20	Titanium

Code	Length (mm)	
12-016-05TI-S	5	Titanium
12-016-06TI-S	6	Titanium
12-016-07TI-S	7	Titanium
12-016-08TI-S	8	Titanium
12-016-09TI-S	9	Titanium
12-016-10TI-S	10	Titanium
12-016-11TI-S	11	Titanium
12-016-12TI-S	12	Titanium
12-016-13TI-S	13	Titanium
12-016-14TI-S	14	Titanium
12-016-16TI-S	16	Titanium
12-016-18TI-S	18	Titanium
12-016-20TI-S	20	Titanium






**2.3mm Wise-Lock Low Profile Screw, Star Head**



**STERILE**

Code	Length (mm)	
12-017-05TI	5	Titanium
12-017-06TI	6	Titanium
12-017-07TI	7	Titanium
12-017-08TI	8	Titanium
12-017-09TI	9	Titanium
12-017-10TI	10	Titanium
12-017-11TI	11	Titanium
12-017-12TI	12	Titanium
12-017-13TI	13	Titanium
12-017-14TI	14	Titanium
12-017-16TI	16	Titanium
12-017-18TI	18	Titanium
12-017-20TI	20	Titanium

Code	Length (mm)	
12-017-05TI-S	5	Titanium
12-017-06TI-S	6	Titanium
12-017-07TI-S	7	Titanium
12-017-08TI-S	8	Titanium
12-017-09TI-S	9	Titanium
12-017-10TI-S	10	Titanium
12-017-11TI-S	11	Titanium
12-017-12TI-S	12	Titanium
12-017-13TI-S	13	Titanium
12-017-14TI-S	14	Titanium
12-017-16TI-S	16	Titanium
12-017-18TI-S	18	Titanium
12-017-20TI-S	20	Titanium



**2.3mm Low Profile Screw, Star Head**



**STERILE**

Code	Length (mm)	
12-018-05TI	5	Titanium
12-018-06TI	6	Titanium
12-018-07TI	7	Titanium
12-018-08TI	8	Titanium
12-018-09TI	9	Titanium
12-018-10TI	10	Titanium
12-018-11TI	11	Titanium
12-018-12TI	12	Titanium
12-018-13TI	13	Titanium
12-018-14TI	14	Titanium
12-018-16TI	16	Titanium
12-018-18TI	18	Titanium
12-018-20TI	20	Titanium

Code	Length (mm)	
12-018-05TI-S	5	Titanium
12-018-06TI-S	6	Titanium
12-018-07TI-S	7	Titanium
12-018-08TI-S	8	Titanium
12-018-09TI-S	9	Titanium
12-018-10TI-S	10	Titanium
12-018-11TI-S	11	Titanium
12-018-12TI-S	12	Titanium
12-018-13TI-S	13	Titanium
12-018-14TI-S	14	Titanium
12-018-16TI-S	16	Titanium
12-018-18TI-S	18	Titanium
12-018-20TI-S	20	Titanium



7-113-01 Drill Bit with Quick Coupling Ø1.1mm x Length 90mm



7-113-02 Drill Bit with AO Coupling Ø1.1mm x Length 90mm



7-113-03 Threaded Drill Guide Ø1.1mm



7-095-22 Stardrive Shaft, T6




7-095-23 Plate Positioner 0.8/1.3mm



7-113-04 Depth Gauge for Ø1.5mm Screw



**7-113 0.8 mm Phalange Implant and Instrument Case**



**7-113 0.8 mm Phalange Implant and Instrument Case**
**Plates**

Code	Set Consisting of	Holes	Units
<b>12-001-02TI</b>	0.8mm Wise-Lock Low Profile Avulsion Fracture Plate	2	1
<b>12-002-03TI</b>	0.8mm Wise-Lock Low Profile Straight Plate	3	1
<b>12-002-04TI</b>	0.8mm Wise-Lock Low Profile Straight Plate	4	1
<b>12-002-05TI</b>	0.8mm Wise-Lock Low Profile Straight Plate	5	1
<b>12-002-06TI</b>	0.8mm Wise-Lock Low Profile Straight Plate	6	1
<b>12-002-07TI</b>	0.8mm Wise-Lock Low Profile Straight Plate	7	1
<b>12-002-10TI</b>	0.8mm Wise-Lock Low Profile Straight Plate	10	1
<b>12-003-09TI</b>	0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate	9	1
<b>12-004-04LTI</b>	0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate, Short, left	4	1
<b>12-004-04RTI</b>	0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate, Short, right	4	1
<b>12-004-05LTI</b>	0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate, Short, left	5	1
<b>12-004-05RTI</b>	0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate, Short, right	5	1
<b>12-004-06LTI</b>	0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate, Short, left	6	1
<b>12-004-06RTI</b>	0.8mm Wise-Lock Low Profile Curved Medial/Lateral Plate, Short, right	6	1
<b>12-005-04TI</b>	0.8mm Wise-Lock Low Profile Compression Plate	4	1
<b>12-005-05TI</b>	0.8mm Wise-Lock Low Profile Compression Plate	5	1
<b>12-005-06TI</b>	0.8mm Wise-Lock Low Profile Compression Plate	6	1
<b>12-006-05TI</b>	0.8mm Wise-Lock Low Profile T-Plate	5	1
<b>12-006-06TI</b>	0.8mm Wise-Lock Low Profile T-Plate	6	1
<b>12-006-08TI</b>	0.8mm Wise-Lock Low Profile T-Plate	8	1
<b>12-006-012TI</b>	0.8mm Wise-Lock Low Profile T-Plate	12	1
<b>12-007-04LTI</b>	0.8mm Wise-Lock Low Profile L-Plate, left	4	1
<b>12-007-04RTI</b>	0.8mm Wise-Lock Low Profile L-Plate, right	4	1
<b>12-007-06LTI</b>	0.8mm Wise-Lock Low Profile L-Plate, left	6	1
<b>12-007-06RTI</b>	0.8mm Wise-Lock Low Profile L-Plate, right	6	1
<b>12-007-08LTI</b>	0.8mm Wise-Lock Low Profile L-Plate, left	8	1
<b>12-007-08RTI</b>	0.8mm Wise-Lock Low Profile L-Plate, right	8	1
<b>12-008-04TI</b>	0.8mm Wise-Lock Low Profile Offset Plate	4	1
<b>12-008-06TI</b>	0.8mm Wise-Lock Low Profile Offset Plate	6	1
<b>12-008-10TI</b>	0.8mm Wise-Lock Low Profile Offset Plate	10	1

**Screws**

Set Consisting of	Length	Units
1.5mm Wise-Lock Low Profile Screw, (Star Head)	5mm to 20mm	6
1.5mm Low Profile Screw, (Star Head)	5mm to 20mm	4
2.3mm Wise-Lock Low Profile Screw, (Star Head)	5mm to 20mm	6
2.3mm Low Profile Screw, (Star Head)	5mm to 20mm	4

**Instruments**

Code	Set Consisting of	Units
<b>7-113-01</b>	Drill Bit with Quick Coupling Ø1.1mm x Length 90mm	1
<b>7-113-02</b>	Drill Bit with AO Coupling Ø1.1mm x Length 90mm	2
<b>7-113-03</b>	Threaded Drill Guide Ø1.1mm	1
<b>7-095-22</b>	Stardrive Shaft, T6	2
<b>7-095-23</b>	Plate Positioner 0.8/1.3mm	1
<b>7-113-04</b>	Depth Gauge for Ø1.5mm Screw	1

**1.3mm Wise-Lock Low Profile Metacarpal Neck Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-009-06LTI-S</b>	<b>12-009-06LTI</b>	6	Titanium
<b>12-009-06RTI-S</b>	<b>12-009-06RTI</b>	6	Titanium



**1.3mm Wise-Lock Low Profile Compression Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-010-04TI-S</b>	<b>12-010-04TI</b>	4	Titanium
<b>12-010-05TI-S</b>	<b>12-010-05TI</b>	5	Titanium
<b>12-010-06TI-S</b>	<b>12-010-06TI</b>	6	Titanium



**1.3mm Wise-Lock Low Profile Straight Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-011-03TI-S</b>	<b>12-011-03TI</b>	3	Titanium
<b>12-011-04TI-S</b>	<b>12-011-04TI</b>	4	Titanium
<b>12-011-05TI-S</b>	<b>12-011-05TI</b>	5	Titanium
<b>12-011-06TI-S</b>	<b>12-011-06TI</b>	6	Titanium
<b>12-011-07TI-S</b>	<b>12-011-07TI</b>	7	Titanium
<b>12-011-10TI-S</b>	<b>12-011-10TI</b>	10	Titanium



**1.3mm Wise-Lock Low Profile T-Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-012-04TI-S</b>	<b>12-012-04TI</b>	4	Titanium
<b>12-012-05TI-S</b>	<b>12-012-05TI</b>	5	Titanium
<b>12-012-06TI-S</b>	<b>12-012-06TI</b>	6	Titanium
<b>12-012-07TI-S</b>	<b>12-012-07TI</b>	7	Titanium
<b>12-012-08TI-S</b>	<b>12-012-08TI</b>	8	Titanium
<b>12-012-11TI-S</b>	<b>12-012-11TI</b>	11	Titanium



**1.3mm Wise-Lock Low Profile Rolando Fracture Hook Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-013-07TI-S</b>	<b>12-013-07TI</b>	7	Titanium



**1.3mm Wise-Lock Low Profile Rotational Correction Plate**

**STERILE**



Sterile	Non Sterile	Holes	
<b>12-014-06TI-S</b>	<b>12-014-06TI</b>	6	Titanium



**1.5mm Wise-Lock Low Profile Screw, Star Head**


Code	Length (mm)	
12-015-05TI	5	Titanium
12-015-06TI	6	Titanium
12-015-07TI	7	Titanium
12-015-08TI	8	Titanium
12-015-09TI	9	Titanium
12-015-10TI	10	Titanium
12-015-11TI	11	Titanium
12-015-12TI	12	Titanium
12-015-13TI	13	Titanium
12-015-14TI	14	Titanium
12-015-16TI	16	Titanium
12-015-18TI	18	Titanium
12-015-20TI	20	Titanium

**STERILE**

Code	Length (mm)	
12-015-05TI-S	5	Titanium
12-015-06TI-S	6	Titanium
12-015-07TI-S	7	Titanium
12-015-08TI-S	8	Titanium
12-015-09TI-S	9	Titanium
12-015-10TI-S	10	Titanium
12-015-11TI-S	11	Titanium
12-015-12TI-S	12	Titanium
12-015-13TI-S	13	Titanium
12-015-14TI-S	14	Titanium
12-015-16TI-S	16	Titanium
12-015-18TI-S	18	Titanium
12-015-20TI-S	20	Titanium


**1.5mm Low Profile Screw, Star Head**


Code	Length (mm)	
12-016-05TI	5	Titanium
12-016-06TI	6	Titanium
12-016-07TI	7	Titanium
12-016-08TI	8	Titanium
12-016-09TI	9	Titanium
12-016-10TI	10	Titanium
12-016-11TI	11	Titanium
12-016-12TI	12	Titanium
12-016-13TI	13	Titanium
12-016-14TI	14	Titanium
12-016-16TI	16	Titanium
12-016-18TI	18	Titanium
12-016-20TI	20	Titanium

**STERILE**

Code	Length (mm)	
12-016-05TI-S	5	Titanium
12-016-06TI-S	6	Titanium
12-016-07TI-S	7	Titanium
12-016-08TI-S	8	Titanium
12-016-09TI-S	9	Titanium
12-016-10TI-S	10	Titanium
12-016-11TI-S	11	Titanium
12-016-12TI-S	12	Titanium
12-016-13TI-S	13	Titanium
12-016-14TI-S	14	Titanium
12-016-16TI-S	16	Titanium
12-016-18TI-S	18	Titanium
12-016-20TI-S	20	Titanium





**2.3mm Wise-Lock Low Profile Screw, Star Head**


Code	Length (mm)	
<b>12-017-05TI</b>	5	Titanium
<b>12-017-06TI</b>	6	Titanium
<b>12-017-07TI</b>	7	Titanium
<b>12-017-08TI</b>	8	Titanium
<b>12-017-09TI</b>	9	Titanium
<b>12-017-10TI</b>	10	Titanium
<b>12-017-11TI</b>	11	Titanium
<b>12-017-12TI</b>	12	Titanium
<b>12-017-13TI</b>	13	Titanium
<b>12-017-14TI</b>	14	Titanium
<b>12-017-16TI</b>	16	Titanium
<b>12-017-18TI</b>	18	Titanium
<b>12-017-20TI</b>	20	Titanium

**STERILE**


Code	Length (mm)	
<b>12-017-05TI-S</b>	5	Titanium
<b>12-017-06TI-S</b>	6	Titanium
<b>12-017-07TI-S</b>	7	Titanium
<b>12-017-08TI-S</b>	8	Titanium
<b>12-017-09TI-S</b>	9	Titanium
<b>12-017-10TI-S</b>	10	Titanium
<b>12-017-11TI-S</b>	11	Titanium
<b>12-017-12TI-S</b>	12	Titanium
<b>12-017-13TI-S</b>	13	Titanium
<b>12-017-14TI-S</b>	14	Titanium
<b>12-017-16TI-S</b>	16	Titanium
<b>12-017-18TI-S</b>	18	Titanium
<b>12-017-20TI-S</b>	20	Titanium


**2.3mm Low Profile Screw, Star Head**


Code	Length (mm)	
<b>12-018-05TI</b>	5	Titanium
<b>12-018-06TI</b>	6	Titanium
<b>12-018-07TI</b>	7	Titanium
<b>12-018-08TI</b>	8	Titanium
<b>12-018-09TI</b>	9	Titanium
<b>12-018-10TI</b>	10	Titanium
<b>12-018-11TI</b>	11	Titanium
<b>12-018-12TI</b>	12	Titanium
<b>12-018-13TI</b>	13	Titanium
<b>12-018-14TI</b>	14	Titanium
<b>12-018-16TI</b>	16	Titanium
<b>12-018-18TI</b>	18	Titanium
<b>12-018-20TI</b>	20	Titanium

**STERILE**

Code	Length (mm)	
<b>12-018-05TI-S</b>	5	Titanium
<b>12-018-06TI-S</b>	6	Titanium
<b>12-018-07TI-S</b>	7	Titanium
<b>12-018-08TI-S</b>	8	Titanium
<b>12-018-09TI-S</b>	9	Titanium
<b>12-018-10TI-S</b>	10	Titanium
<b>12-018-11TI-S</b>	11	Titanium
<b>12-018-12TI-S</b>	12	Titanium
<b>12-018-13TI-S</b>	13	Titanium
<b>12-018-14TI-S</b>	14	Titanium
<b>12-018-16TI-S</b>	16	Titanium
<b>12-018-18TI-S</b>	18	Titanium
<b>12-018-20TI-S</b>	20	Titanium




7-114-01 Drill Bit with Quick Coupling Ø2.0mm x Length 90mm



7-114-02 Drill Bit with AO Coupling Ø2.0mm x Length 90mm



7-114-03 Threaded Drill Guide Ø2.0mm



7-095-22 Stardrive Shaft, T6



7-095-23 Plate Positioner 0.8/1.3mm



7-114-04 Depth Gauge for Ø2.3mm Screw



7-114 1.3 mm Metacarpal Implant and Instrument Case



**7-114 1.3 mm Metacarpal Implant and Instrument Case**
**Plates**

Code	Set Consisting of	Holes	Units
<b>12-009-06LTI</b>	1.3mm Wise-Lock Low Profile Metacarpal Neck Plate, Left	6	1
<b>12-009-06RTI</b>	1.3mm Wise-Lock Low Profile Metacarpal Neck Plate, Right	6	1
<b>12-010-04TI</b>	1.3mm Wise-Lock Low Profile Compression Plate	4	1
<b>12-010-05TI</b>	1.3mm Wise-Lock Low Profile Compression Plate	5	1
<b>12-010-06TI</b>	1.3mm Wise-Lock Low Profile Compression Plate	6	1
<b>12-011-03TI</b>	1.3mm Wise-Lock Low Profile Straight Plate	3	1
<b>12-011-04TI</b>	1.3mm Wise-Lock Low Profile Straight Plate	4	1
<b>12-011-05TI</b>	1.3mm Wise-Lock Low Profile Straight Plate	5	1
<b>12-011-06TI</b>	1.3mm Wise-Lock Low Profile Straight Plate	6	1
<b>12-011-07TI</b>	1.3mm Wise-Lock Low Profile Straight Plate	7	1
<b>12-011-10TI</b>	1.3mm Wise-Lock Low Profile Straight Plate	10	1
<b>12-012-04TI</b>	1.3mm Wise-Lock Low Profile T-Plate	4	1
<b>12-012-05TI</b>	1.3mm Wise-Lock Low Profile T-Plate	5	1
<b>12-012-06TI</b>	1.3mm Wise-Lock Low Profile T-Plate	6	1
<b>12-012-07TI</b>	1.3mm Wise-Lock Low Profile T-Plate	7	1
<b>12-012-08TI</b>	1.3mm Wise-Lock Low Profile T-Plate	8	1
<b>12-012-11TI</b>	1.3mm Wise-Lock Low Profile T-Plate	11	1
<b>12-013-07TI</b>	1.3mm Wise-Lock Low Profile Rolando Fracture Hook Plate	7	1
<b>12-014-06TI</b>	1.3mm Wise-Lock Low Profile Rotational Correction Plate	6	1

**Screws**

Set Consisting of	Length	Units
1.5mm Wise-Lock Low Profile screw	5mm to 20mm	6
1.5mm Low Profile screw	5mm to 20mm	4
2.3mm Wise-Lock Low Profile screw	5mm to 20mm	6
2.3mm Low Profile screw	5mm to 20mm	4

**Instruments**

Code	Set Consisting of	Units
<b>7-114-01</b>	Ø2.0mm Drill Bit x Length 90mm with Quick Coupling	1
<b>7-114-02</b>	Ø2.0mm Drill Bit x Length 90mm with AO Coupling	2
<b>7-114-03</b>	Ø2.0mm Drill Guide	1
<b>7-095-22</b>	T6 Stardrive Shaft	2
<b>7-095-23</b>	0.8/1.3mm Plate Positioner	1
<b>7-114-04</b>	Depth Gauge for Ø2.3mm Screw	1

7-095-01

Bending Pliers for Low Profile Hand System




7-095-02

Double Periosteal Elevator, Flat, for Low Profile Hand System




7-095-03

Countersink for Low Profile Hand System



7-095-04

Position Pin for Low Profile Hand System



7-095-05

Retractor for Low Profile Hand System



7-095-06

Cleaning Hook for Low Profile Hand System



**7-095-07** Bone Forcep 1.1/2.0 for Low Profile Hand System



**7-095-08** Reduction Forcep with Points, for Low Profile Hand System



**7-095-09** Reduction Forcep for Low Profile Hand System



**7-095-10** Osteotomy Guide for Low Profile Hand System



**7-095-11** Distractor for Low Profile Hand System



**7-095-12** Kirschner Wire with Threaded Trocar Tip, Ø1.0mm x Thread Length 19mm x Length 150mm



**7-095-13** Kirschner Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 19mm x Length 150mm



**7-095-14** Kirschner Wire with Threaded Trocar Tip, Ø1.6mm x Thread Length 19mm x Length 150mm



**7-095-15** Kirschner Wire, Ø1.0mm x Length 150mm



**7-095-16** Kirschner Wire, Ø1.2mm x Length 150mm



**7-095-17** Kirschner Wire, Ø1.6mm x Length 150mm




**7-095-18** Handle with Mini Quick Coupling for Low Profile Hand System




**7-095-19** Mini Quick Coupling (AO) for Low Profile Hand System




**7-095-20** Instrument Trays for 0.8-1.3mm Low Profile Hand Instrument Set



**7-095-21** Container for 0.8-1.3mm Low Profile Hand Instrument Set





**7-095 0.8-1.3mm Low Profile Hand System Instrument Set**


Item Code	Item Name	Unit
7-095-01	Bending Pliers for Low Profile Hand System	2
7-095-02	Double Periosteal Elevator, Flat, for Low Profile Hand System	1
7-095-03	Countersink for Low Profile Hand System	2
7-095-04	Position Pin for Low Profile Hand System	4
7-095-05	Retractor for Low Profile Hand System	2
7-095-06	Cleaning Hook for Low Profile Hand System	1
7-095-07	Bone Forcep 1.1/2.0 for Low Profile Hand System	1
7-095-08	Reduction Forcep with Points, for Low Profile Hand System	1
7-095-09	Reduction Forcep for Low Profile Hand System	1
7-095-10	Osteotomy Guide for Low Profile Hand System	1
7-095-11	Distractor for Low Profile Hand System	1
7-095-12	Kirschner Wire with Threaded Trocar Tip, Ø1.0mm x Thread Length 19mm x Length 150mm	2
7-095-13	Kirschner Wire with Threaded Trocar Tip, Ø1.2mm x Thread Length 19mm x Length 150mm	2
7-095-14	Kirschner Wire with Threaded Trocar Tip, Ø1.6mm x Thread Length 19mm x Length 150mm	2
7-095-15	Kirschner Wire, Ø1.0mm x Length 150mm	2
7-095-16	Kirschner Wire, Ø1.2mm x Length 150mm	2
7-095-17	Kirschner Wire, Ø1.6mm x Length 150mm	2
7-095-18	Handle with Mini Quick Coupling for Low Profile Hand System	1
7-095-19	Mini Quick Coupling (AO) for Low Profile Hand System	1
7-095-20	Instrument Trays for 0.8-1.3mm Low Profile Hand Instrument Set	1
7-095-21	Container for 0.8-1.3mm Low Profile Hand Instrument Set	1







**USA**

Auxein Inc.  
1500 Nw 89th Court, Suite 107-108  
Doral, Florida 33172  
Tel: +1 305 395 6062  
E Fax: +1 305 395 6262  
Email: USoffice@auxein.com

**MEXICO**

Auxein México, S.A. de C.V.  
Tepic 139 int 801, Colonia Roma Sur,  
Alcaldía Cuauhtémoc, CDMX,  
México, C.P. 06760  
Tel: +521 55 7261 0318  
Email: info@auxein.mx

**INDIA**

Auxein Medical Pvt. Ltd.  
Plot No. 168-169-170, Phase-4,  
Kundli Industrial Area,  
HSIIDC, Sector-57, Sonapat - 131028, Haryana  
Tel: +91 99106 43638 | Fax: +91 86077 70197  
Email: info@auxein.com