


ANKLE ARTHRODESIS NAILING SYSTEM

ANKLE

ARTHRODESIS NAILING SYSTEM



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements




Orthopaedic Implants

Indications

1. Distal tibial fracture and non-unions.
2. Severe foot and ankle deformity.
3. Osteoarthritis.
4. Arthritis.
5. Instability and skeletal defects, including but not limited to after tumor resection.
6. Avascular necrosis of the Talus.
7. Failed joint replacement or failed ankle fusion.
8. Rheumatoid arthritis and pseudoarthrosis.






Design Features



1. Provide Static and Dynamic Proximal Locking option for clinical choice.
2. Single piece of Aiming Shaft ensure locking accuracy.
3. Distal sliding hole allows maximal 5mm compression distance.
4. Distal multi-planer screws enhance angular stability.
5. Available in different sizes.
6. Left and Right Configuration.

Design Features



Ankle Arthrodesis Nail
STERILE

Ø Dia	Length	Item Code
10mm	160mm	1-076-160LTI-S
10mm	200mm	1-076-200LTI-S
10mm	250mm	1-076-250LTI-S
11.5mm	160mm	1-083-160LTI-S
11.5mm	200mm	1-083-200LTI-S
11.5mm	250mm	1-083-250LTI-S



Ø Dia	Length	Item Code
10mm	160mm	1-076-160LTI
10mm	200mm	1-076-200LTI
10mm	250mm	1-076-250LTI
11.5mm	160mm	1-083-160LTI
11.5mm	200mm	1-083-200LTI
11.5mm	250mm	1-083-250LTI

STERILE

Ø Dia	Length	Item Code
10mm	160mm	1-076-160RTI-S
10mm	200mm	1-076-200RTI-S
10mm	250mm	1-076-250RTI-S
11.5mm	160mm	1-083-160RTI-S
11.5mm	200mm	1-083-200RTI-S
11.5mm	250mm	1-083-250RTI-S






Ø Dia	Length	Item Code
10mm	160mm	1-076-160RTI
10mm	200mm	1-076-200RTI
10mm	250mm	1-076-250RTI
11.5mm	160mm	1-083-160RTI
11.5mm	200mm	1-083-200RTI
11.5mm	250mm	1-083-250RTI





5.0mm Locking Bolt for Ankle Arthrodesis Nail (Hex Head)

		
Dia	Length	Item Code
5.0mm	26mm	1-086-26TI
5.0mm	28mm	1-086-28TI
5.0mm	30mm	1-086-30TI
5.0mm	32mm	1-086-32TI
5.0mm	34mm	1-086-34TI
5.0mm	36mm	1-086-36TI
5.0mm	38mm	1-086-38TI
5.0mm	40mm	1-086-40TI
5.0mm	42mm	1-086-42TI
5.0mm	44mm	1-086-44TI
5.0mm	46mm	1-086-46TI
5.0mm	48mm	1-086-48TI
5.0mm	50mm	1-086-50TI
5.0mm	52mm	1-086-52TI
5.0mm	54mm	1-086-54TI
5.0mm	56mm	1-086-56TI
5.0mm	58mm	1-086-58TI
5.0mm	60mm	1-086-60TI
5.0mm	62mm	1-086-62TI
5.0mm	64mm	1-086-64TI
5.0mm	66mm	1-086-66TI
5.0mm	68mm	1-086-68TI
5.0mm	70mm	1-086-70TI
5.0mm	72mm	1-086-72TI
5.0mm	74mm	1-086-74TI
5.0mm	76mm	1-086-76TI
5.0mm	78mm	1-086-78TI
5.0mm	80mm	1-086-80TI
5.0mm	85mm	1-086-85TI
5.0mm	90mm	1-086-90TI
5.0mm	95mm	1-086-95TI
5.0mm	100mm	1-086-100TI
5.0mm	105mm	1-086-105TI
5.0mm	110mm	1-086-110TI



5.0mm Locking Bolt for Ankle Arthrodesis Nail (Hex Head)
STERILE

Dia	Length	Item Code
5.0mm	26mm	1-086-26TI-S
5.0mm	28mm	1-086-28TI-S
5.0mm	30mm	1-086-30TI-S
5.0mm	32mm	1-086-32TI-S
5.0mm	34mm	1-086-34TI-S
5.0mm	36mm	1-086-36TI-S
5.0mm	38mm	1-086-38TI-S
5.0mm	40mm	1-086-40TI-S
5.0mm	42mm	1-086-42TI-S
5.0mm	44mm	1-086-44TI-S
5.0mm	46mm	1-086-46TI-S
5.0mm	48mm	1-086-48TI-S
5.0mm	50mm	1-086-50TI-S
5.0mm	52mm	1-086-52TI-S
5.0mm	54mm	1-086-54TI-S
5.0mm	56mm	1-086-56TI-S
5.0mm	58mm	1-086-58TI-S
5.0mm	60mm	1-086-60TI-S
5.0mm	62mm	1-086-62TI-S
5.0mm	64mm	1-086-64TI-S
5.0mm	66mm	1-086-66TI-S
5.0mm	68mm	1-086-68TI-S
5.0mm	70mm	1-086-70TI-S
5.0mm	72mm	1-086-72TI-S
5.0mm	74mm	1-086-74TI-S
5.0mm	76mm	1-086-76TI-S
5.0mm	78mm	1-086-78TI-S
5.0mm	80mm	1-086-80TI-S
5.0mm	85mm	1-086-85TI-S
5.0mm	90mm	1-086-90TI-S
5.0mm	95mm	1-086-95TI-S
5.0mm	100mm	1-086-100TI-S
5.0mm	105mm	1-086-105TI-S
5.0mm	110mm	1-086-110TI-S




End cap for Ankle Arthrodesis Nail

STERILE

Length	Item Code
0 mm	1-084-01TI-S
2.5mm	1-084-02TI-S
5.0mm	1-084-03TI-S
7.5mm	1-084-04TI-S
10.0mm	1-084-05TI-S



Length	Item Code
0 mm	1-084-01TI
2.5mm	1-084-02TI
5.0mm	1-084-03TI
7.5mm	1-084-04TI
10.0mm	1-084-05TI



Compression Screw for Ankle Arthrodesis Nail (Hex Head), Titanium

STERILE

Item Code

1-085-01TI-S


Item Code

1-085-01TI


7-174-01 Radiographic Ruler for Ankle Arthrodesis Nail



7-174-02 Guide Wire with Threaded Tip Ø3.2mm x Thread Length 10.0mm X Length 300mm



7-174-03 Guide Wire with Threaded Tip Ø3.2mm x Thread Length 10.0mm X Length 345mm



7-174-04 Entry Reamer Protection Sleeve Ø15/12.5mm for Ankle Arthrodesis Nail



7-174-05 Guide Wire Sleeve Ø12.5/3.2mm for Ankle Arthrodesis Nail



7-174-06 Cannulated Awl for Ankle Arthrodesis Nail



7-174-07 Cannulated Entry Reamer Ø12.5 / Ø3.2mm for Ankle Arthrodesis Nail



7-174-08 Adaptor with Jacob Coupling for Ankle Arthrodesis Nail



7-174-09 T- Handle with Jacob Coupling for Ankle Arthrodesis Nail



7-174-10 Reduction Device for Ankle Arthrodesis Nail



7-174-11 Guide Wire Interducer for Ankle Arthrodesis Nail




7-174-12 Guide Wire with Ball Tip Ø2.5/4.0mm x Length 680mm, Nitinol for Ankle Arthrodesis Nail



7-174-13 Flexible Reamer with Shaft Ø9mm for Ankle Arthrodesis Nail




7-174-14 Flexible Reamer with Shaft Ø10mm for Ankle Arthrodesis Nail



7-174-15 Flexible Reamer with Shaft Ø11mm for Ankle Arthrodesis Nail




7-174-16 Flexible Reamer with Shaft Ø12mm for Ankle Arthrodesis Nail



7-174-17 Flexible Reamer with Shaft Ø13mm for Ankle Arthrodesis Nail



7-174-18 Insertion Handle for Ankle Arthrodesis Nail



7-174-19

Nail Holding Bolt (M8x1) for Ankle Arthrodesis Nail


7-174-20

Wrench SW11 for Ankle Arthrodesis Nail


7-174-21

Cannulated Impactor-Extractor Rod (M8x1) for Ankle Arthrodesis Nail


7-174-22

Sliding Hammer for Ankle Arthrodesis Nail


7-174-23

Outer Protection Sleeve Ø10/8.0mm, Length 155mm for Ankle Arthrodesis Nail


7-174-24

Drill Protection Sleeve Ø8.0 / Ø4.0mm, Length 170mm for Ankle Arthrodesis Nail



7-174-25 Trocar Ø4.0 x 186mm Length for Ankle Arthrodesis Nail



7-174-26 Drill Bit with Jacob Coupling, Ø4.0mm x Length 300mm for Ankle Arthrodesis Nail



7-174-27 Stopper for Ø4.0mm Drill bit



7-174-28 Allen Key, Hex 3mm for Ankle Arthrodesis Nail



7-174-29 Depth Gauge upto 110mm for Ankle Arthrodesis Nail



7-174-30 T-Handle Screwdriver, SW4.5 for Ankle Arthrodesis Nail



7-174-31 Inner Threaded Bar, M3 for T-Handle Screwdriver




7-174-32 Hex Screwdriver Shaft, SW4.5 with Jacob Coupling for Ankle Arthrodesis Nail




7-174-33 Hex Screwdriver, SW3.5 for Ankle Arthrodesis Nail




7-174-34 Implant Tray for Ankle Arthrodesis Nail



7-174-35 Instrument Tray-I for Ankle Arthrodesis Nail




7-174-36 Instrument Tray-II for Ankle Arthrodesis Nail



7-174-37

Instrument Box with Cover for Ankle Arthrodesis Nail



7-174 Ankle Arthrodesis Nailing Instruments Set


7-174 Arthrodesis Nailing Instruments Set

Code	Set Consisting of	Units
7-174-01	Radiographic Ruler for Ankle Arthrodesis Nail	1
7-174-02	Guide Wire with Threaded Tip Ø3.2mm x Thread Length 10.0mm X Length 300mm	2
7-174-03	Guide Wire with Threaded Tip Ø3.2mm x Thread Length 10.0mm X Length 345mm	2
7-174-04	Entry Reamer Protection Sleeve Ø15/12.5mm for Ankle Arthrodesis Nail	1
7-174-05	Guide Wire Sleeve ØØ12.5/3.2mm for Ankle Arthrodesis Nail	1
7-174-06	Cannulated Awl for Ankle Arthrodesis Nail	1
7-174-07	Cannulated Entry Reamer Ø12.5 / Ø3.2mm for Ankle Arthrodesis Nail	1
7-174-08	Adaptor with Jacob Coupling for Ankle Arthrodesis Nail	1
7-174-09	T- Handle with Jacob Coupling for Ankle Arthrodesis Nail	1
7-174-10	Reduction Device for Ankle Arthrodesis Nail	1
7-174-11	Guide Wire Interducer for Ankle Arthrodesis Nail	1
7-174-12	Guide Wire with Ball Tip Ø2.5/4.0mm x Length 680mm, Nitinol for Ankle Arthrodesis Nail	2
7-174-13	Flexible Reamer with Shaft Ø9mm for Ankle Arthrodesis Nail	1
7-174-14	Flexible Reamer with Shaft Ø10mm for Ankle Arthrodesis Nail	1
7-174-15	Flexible Reamer with Shaft Ø11mm for Ankle Arthrodesis Nail	1
7-174-16	Flexible Reamer with Shaft Ø12mm for Ankle Arthrodesis Nail	1
7-174-17	Flexible Reamer with Shaft Ø13mm for Ankle Arthrodesis Nail	1
7-174-18	Insertion Handle for Ankle Arthrodesis Nail	1
7-174-19	Nail Holding Bolt (M8x1) for Ankle Arthrodesis Nail	2
7-174-20	Wrench SW11 for Ankle Arthrodesis Nail	1
7-174-21	Cannulated Impactor-Extractor Rod (M8x1) for Ankle Arthrodesis Nail	1
7-174-22	Sliding Hammer for Ankle Arthrodesis Nail	1
7-174-23	Outer Protection Sleeve Ø10/8.0mm, Length 155mm for Ankle Arthrodesis Nail	2
7-174-24	Drill Sleeve Ø8.0/4.0mm, Length 170mm for Ankle Arthrodesis Nail	2
7-174-25	Trocar Ø4.0 x 186mm Length for Ankle Arthrodesis Nail	2
7-174-26	Drill Bit with Jacob Coupling, Ø4.0mm x Length 300mm for Ankle Arthrodesis Nail	1
7-174-27	Stopper for Ø4.0mm Drill bit	1
7-174-28	Allen Key, Hex 3mm for Ankle Arthrodesis Nail	1
7-174-29	Depth Gauge upto 110mm for Ankle Arthrodesis Nail	1
7-174-30	T-Handle Screwdriver, SW4.5 for Ankle Arthrodesis Nail	1
7-174-31	Inner Threaded Bar, M3 for T-Handle Screwdriver	1
7-174-32	Hex Screwdriver Shaft, SW4.5 with Jacob Coupling for Ankle Arthrodesis Nail	1
7-174-33	Hex Screwdriver, SW3.5 for Ankle Arthrodesis Nail	1
7-174-34	Implant Tray for Ankle Arthrodesis Nail	1
7-174-35	Instrument Tray-I for Ankle Arthrodesis Nail	1
7-174-36	Intrument Tray-II for Ankle Arthrodesis Nail	1
7-174-37	Instrument Box with Cover for Ankle Arthrodesis Nail	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com