

# JIN TYPE - HUMERUS NAILING SYSTEM

## JIN TYPE - HUMERUS NAILING SYSTEM





# about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

## Our Achievements



# **Orthopaedic Implants**

**JIN Type - Humerus Nail**


Ø Dia	Length	Stainless Steel	Titanium
6mm	180mm	<b>1296-6-180</b>	<b>TI-1296-6-180</b>
6mm	200mm	<b>1296-6-200</b>	<b>TI-1296-6-200</b>
6mm	220mm	<b>1296-6-220</b>	<b>TI-1296-6-220</b>
6mm	240mm	<b>1296-6-240</b>	<b>TI-1296-6-240</b>
6mm	260mm	<b>1296-6-260</b>	<b>TI-1296-6-260</b>
6mm	280mm	<b>1296-6-280</b>	<b>TI-1296-6-280</b>
6mm	300mm	<b>1296-6-300</b>	<b>TI-1296-6-300</b>
6mm	320mm	<b>1296-6-320</b>	<b>TI-1296-6-320</b>
7mm	180mm	<b>1296-7-180</b>	<b>TI-1296-7-180</b>
7mm	200mm	<b>1296-7-200</b>	<b>TI-1296-7-200</b>
7mm	220mm	<b>1296-7-220</b>	<b>TI-1296-7-220</b>
7mm	240mm	<b>1296-7-240</b>	<b>TI-1296-7-240</b>
7mm	260mm	<b>1296-7-260</b>	<b>TI-1296-7-260</b>
7mm	280mm	<b>1296-7-280</b>	<b>TI-1296-7-280</b>
7mm	300mm	<b>1296-7-300</b>	<b>TI-1296-7-300</b>
7mm	320mm	<b>1296-7-320</b>	<b>TI-1296-7-320</b>
8mm	180mm	<b>1296-8-180</b>	<b>TI-1296-8-180</b>
8mm	200mm	<b>1296-8-200</b>	<b>TI-1296-8-200</b>
8mm	220mm	<b>1296-8-220</b>	<b>TI-1296-8-220</b>
8mm	240mm	<b>1296-8-240</b>	<b>TI-1296-8-240</b>
8mm	260mm	<b>1296-8-260</b>	<b>TI-1296-8-260</b>
8mm	280mm	<b>1296-8-280</b>	<b>TI-1296-8-280</b>
8mm	300mm	<b>1296-8-300</b>	<b>TI-1296-8-300</b>
8mm	320mm	<b>1296-8-320</b>	<b>TI-1296-8-320</b>
9mm	180mm	<b>1296-9-180</b>	<b>TI-1296-9-180</b>
9mm	200mm	<b>1296-9-200</b>	<b>TI-1296-9-200</b>
9mm	220mm	<b>1296-9-220</b>	<b>TI-1296-9-220</b>
9mm	240mm	<b>1296-9-240</b>	<b>TI-1296-9-240</b>
9mm	260mm	<b>1296-9-260</b>	<b>TI-1296-9-260</b>
9mm	280mm	<b>1296-9-280</b>	<b>TI-1296-9-280</b>
9mm	300mm	<b>1296-9-300</b>	<b>TI-1296-9-300</b>
9mm	320mm	<b>1296-9-320</b>	<b>TI-1296-9-320</b>
10mm	180mm	<b>1296-10-180</b>	<b>TI-1296-10-180</b>
10mm	200mm	<b>1296-10-200</b>	<b>TI-1296-10-200</b>
10mm	220mm	<b>1296-10-220</b>	<b>TI-1296-10-220</b>
10mm	240mm	<b>1296-10-240</b>	<b>TI-1296-10-240</b>
10mm	260mm	<b>1296-10-260</b>	<b>TI-1296-10-260</b>
10mm	280mm	<b>1296-10-280</b>	<b>TI-1296-10-280</b>
10mm	300mm	<b>1296-10-300</b>	<b>TI-1296-10-300</b>
10mm	320mm	<b>1296-10-320</b>	<b>TI-1296-10-320</b>



Stainless Steel

Titanium

**JIN Type - Humerus Nail**
**STERILE**

Ø Dia	Length	Stainless Steel	Titanium
6mm	180mm	<b>1296-6-180-S</b>	<b>TI-1296-6-180-S</b>
6mm	200mm	<b>1296-6-200-S</b>	<b>TI-1296-6-200-S</b>
6mm	220mm	<b>1296-6-220-S</b>	<b>TI-1296-6-220-S</b>
6mm	240mm	<b>1296-6-240-S</b>	<b>TI-1296-6-240-S</b>
6mm	260mm	<b>1296-6-260-S</b>	<b>TI-1296-6-260-S</b>
6mm	280mm	<b>1296-6-280-S</b>	<b>TI-1296-6-280-S</b>
6mm	300mm	<b>1296-6-300-S</b>	<b>TI-1296-6-300-S</b>
6mm	320mm	<b>1296-6-320-S</b>	<b>TI-1296-6-320-S</b>
7mm	180mm	<b>1296-7-180-S</b>	<b>TI-1296-7-180-S</b>
7mm	200mm	<b>1296-7-200-S</b>	<b>TI-1296-7-200-S</b>
7mm	220mm	<b>1296-7-220-S</b>	<b>TI-1296-7-220-S</b>
7mm	240mm	<b>1296-7-240-S</b>	<b>TI-1296-7-240-S</b>
7mm	260mm	<b>1296-7-260-S</b>	<b>TI-1296-7-260-S</b>
7mm	280mm	<b>1296-7-280-S</b>	<b>TI-1296-7-280-S</b>
7mm	300mm	<b>1296-7-300-S</b>	<b>TI-1296-7-300-S</b>
7mm	320mm	<b>1296-7-320-S</b>	<b>TI-1296-7-320-S</b>
8mm	180mm	<b>1296-8-180-S</b>	<b>TI-1296-8-180-S</b>
8mm	200mm	<b>1296-8-200-S</b>	<b>TI-1296-8-200-S</b>
8mm	220mm	<b>1296-8-220-S</b>	<b>TI-1296-8-220-S</b>
8mm	240mm	<b>1296-8-240-S</b>	<b>TI-1296-8-240-S</b>
8mm	260mm	<b>1296-8-260-S</b>	<b>TI-1296-8-260-S</b>
8mm	280mm	<b>1296-8-280-S</b>	<b>TI-1296-8-280-S</b>
8mm	300mm	<b>1296-8-300-S</b>	<b>TI-1296-8-300-S</b>
8mm	320mm	<b>1296-8-320-S</b>	<b>TI-1296-8-320-S</b>
9mm	180mm	<b>1296-9-180-S</b>	<b>TI-1296-9-180-S</b>
9mm	200mm	<b>1296-9-200-S</b>	<b>TI-1296-9-200-S</b>
9mm	220mm	<b>1296-9-220-S</b>	<b>TI-1296-9-220-S</b>
9mm	240mm	<b>1296-9-240-S</b>	<b>TI-1296-9-240-S</b>
9mm	260mm	<b>1296-9-260-S</b>	<b>TI-1296-9-260-S</b>
9mm	280mm	<b>1296-9-280-S</b>	<b>TI-1296-9-280-S</b>
9mm	300mm	<b>1296-9-300-S</b>	<b>TI-1296-9-300-S</b>
9mm	320mm	<b>1296-9-320-S</b>	<b>TI-1296-9-320-S</b>
10mm	180mm	<b>1296-10-180-S</b>	<b>TI-1296-10-180-S</b>
10mm	200mm	<b>1296-10-200-S</b>	<b>TI-1296-10-200-S</b>
10mm	220mm	<b>1296-10-220-S</b>	<b>TI-1296-10-220-S</b>
10mm	240mm	<b>1296-10-240-S</b>	<b>TI-1296-10-240-S</b>
10mm	260mm	<b>1296-10-260-S</b>	<b>TI-1296-10-260-S</b>
10mm	280mm	<b>1296-10-280-S</b>	<b>TI-1296-10-280-S</b>
10mm	300mm	<b>1296-10-300-S</b>	<b>TI-1296-10-300-S</b>
10mm	320mm	<b>1296-10-320-S</b>	<b>TI-1296-10-320-S</b>



Stainless Steel

Titanium

### End Cap For JIN Type - Humerus Nail



Stainless Steel	Titanium
1296-001	TI-1296-001



**STERILE**

Stainless Steel	Titanium
1296-001-S	TI-1296-001-S

### Ø3.5mm Locking Bolt, Self-Tapping, For JIN Type - Humerus Nail




Length	Stainless Steel	Titanium
20mm	1296-3.5-20	TI-1296-3.5-20
25mm	1296-3.5-25	TI-1296-3.5-25
30mm	1296-3.5-30	TI-1296-3.5-30
35mm	1296-3.5-35	TI-1296-3.5-35
40mm	1296-3.5-40	TI-1296-3.5-40
45mm	1296-3.5-45	TI-1296-3.5-45
50mm	1296-3.5-50	TI-1296-3.5-50
55mm	1296-3.5-55	TI-1296-3.5-55
60mm	1296-3.5-60	TI-1296-3.5-60
65mm	1296-3.5-65	TI-1296-3.5-65
70mm	1296-3.5-70	TI-1296-3.5-70
75mm	1296-3.5-75	TI-1296-3.5-75
80mm	1296-3.5-80	TI-1296-3.5-80


**STERILE**

Length	Stainless Steel	Titanium
20mm	1296-3.5-20-S	TI-1296-3.5-20-S
25mm	1296-3.5-25-S	TI-1296-3.5-25-S
30mm	1296-3.5-30-S	TI-1296-3.5-30-S
35mm	1296-3.5-35-S	TI-1296-3.5-35-S
40mm	1296-3.5-40-S	TI-1296-3.5-40-S
45mm	1296-3.5-45-S	TI-1296-3.5-45-S
50mm	1296-3.5-50-S	TI-1296-3.5-50-S
55mm	1296-3.5-55-S	TI-1296-3.5-55-S
60mm	1296-3.5-60-S	TI-1296-3.5-60-S
65mm	1296-3.5-65-S	TI-1296-3.5-65-S
70mm	1296-3.5-70-S	TI-1296-3.5-70-S
75mm	1296-3.5-75-S	TI-1296-3.5-75-S
80mm	1296-3.5-80-S	TI-1296-3.5-80-S

**1-017 Implant Box for JIN Type - Humerus Nail**






**Tray 1**



Dia	Length	Qty
Ø6mm	180mm to 320mm	1pcs

*Stainless Steel or Titanium Implants can be placed as per requirement*

**Tray 2**



Dia	Length	Qty
Ø7mm	180mm to 320mm	1pcs

*Stainless Steel or Titanium Implants can be placed as per requirement*

**Tray 3**



Dia	Length	Qty
Ø8mm	180mm to 320mm	1pcs

*Stainless Steel or Titanium Implants can be placed as per requirement*

**Tray 4**



Dia	Length	Qty
Ø9mm	180mm to 320mm	1pcs

*Stainless Steel or Titanium Implants can be placed as per requirement*

**Tray 5**


Dia	Length	Qty
Ø10mm	180mm to 320mm	1pcs

*Stainless Steel or Titanium Implants can be placed as per requirement*



### Ø3.5mm Locking Bolt, Self-Tapping, For JIN Type - Humerus Nail

Dia	Length	Qty
Ø3.5mm	20mm to 80mm	6pcs


*Stainless Steel or Titanium Implants  
can be placed as per requirement*

### End Cap For JIN Type - Humerus Nail

Qty
2pcs

*Stainless Steel or Titanium Implants  
can be placed as per requirement*

**INS-1296-01** Proximal Targeting Arm for JIN Type - Humerus Nail



**INS-1296-02** Distal Targeting Arm for JIN Type - Humerus Nail



**INS-1296-03** Bolt for Distal Targeting Device M6 - Hex 5mm, for JIN Type - Humerus Nail



**INS-1296-04** Positioning Rod Ø5mm for JIN Type - Humerus Nail



**INS-1296-05** Positioning Block for JIN Type - Humerus Nail




**INS-1296-06** Trocar Ø6mm for JIN Type - Humerus Nail



**INS-1296-08** Drill Sleeve 8/5.7mm for Drill Bit Ø5.5mm - JIN Type - Humerus Nail



**INS-1296-10** Drill Sleeve Ø8/6mm for JIN Type - Humerus Nail



**INS-1296-11** Depth Gauge Drill Sleeve Ø6mm for Drill Bit Ø2.9/3.5mm - JIN Type - Humerus Nail



**INS-1296-12** Insertion Handle for JIN Type - Humerus Nail



**INS-1296-13** Nail Holding Bolt for JIN Type - Humerus Nail



**INS-1296-15** Step Drill Bit Plain Shank/Jacob Chuck End Ø2.9/3.5mm x Length 250mm for JIN Type - Humerus Nail



**INS-1296-16** Distal Targeting Device for Length 180mm To 320mm - JIN Type - Humerus Nail



**INS-1296-07** Drill Bit Plain Shank/Jacob Chuck End Ø5.5mm x Length 300mm, for JIN Type - Humerus Nail



**INS-1296-17** Cannulated Awl for JIN Type - Humerus Nail



**INS-1296-18** Combination Wrench 12mm for JIN Type - Humerus Nail



**INS-1296-19** Handle For Reamer (AO Connection) for JIN Type - Humerus Nail



**INS-1296-20** Allen Key, Hex 5mm for JIN Type - Humerus Nail



**INS-1296-21** T-Handle Screwdriver, Hex 3.5mm, for JIN Type - Humerus Nail



**INS-1296-22** Ø5.5mm T-Handle Drill Bit with Flat End for JIN Type - Humerus Nail



**INS-1296-23** T-Handle Tap for Ø3.5mm Locking Bolt - JIN Type - Humerus Nail



**INS-1296-24** T-Handle Screwdriver, Hex 2.5mm, for JIN Type - Humerus Nail



**INS-1296-27** Fixed Reamer Ø8.0mm x Length 400mm with Quick Coupling End for JIN Type - Humerus Nail



**INS-1296-28**Fixed Reamer  $\varnothing 9.0\text{mm}$  x Length 400mm with Quick Coupling End for JIN Type - Humerus Nail**INS-1296-29**Fixed Reamer  $\varnothing 10.0\text{mm}$  x Length 400mm with Quick Coupling End for JIN Type - Humerus Nail**7-023-04**Flexible Reamer  $\varnothing 6.5\text{mm}$  for JIN Type - Humerus Nail**7-023-05**Flexible Reamer  $\varnothing 7.0\text{mm}$  for JIN Type - Humerus Nail**7-023-06**Flexible Reamer  $\varnothing 8.0\text{mm}$  for JIN Type - Humerus Nail**7-023-07**Flexible Reamer  $\varnothing 9.0\text{mm}$  for JIN Type - Humerus Nail

**7-023-13** Flexible Reamer Ø10.0mm for JIN Type - Humerus Nail



**7-023-01** Protection Sleeve for Entry Reamer - JIN Type - Humerus Nail



**7-023-03** Quick Coupling Handle for Jacob Connection for JIN Type - Humerus Nail



**7-023-08** Radiographic Ruler for JIN Type - Humerus Nail



**7-023-15** Guide Wire with Threaded Tip, Ø2.5mm x Thread Length 10mm x Length 200mm



**7-023-14** Extraction Handle for JIN Type- Humerus Nail






**INS-1296-30** Extraction Rod with Ram for JIN Type- Humerus Nail



**7-023-02** Entry Reamer for JIN Type- Humerus Nail



**7-023-09** Wire Pusher for JIN Type- Humerus Nail




**7-023-11** Reduction Device for JIN Type- Humerus Nail



**7-023-12** Plastic Exchange Tube for JIN Type- Humerus Nail




**INS-1296-32** Guide Wire with Ball Tip Ø2.5/3.2mm x Length 1000mm - Stainless Steel, for JIN Type- Humerus Nail




**7-023-16**

Guide Wire without Ball Tip Ø2.5mm x Length 1000mm - Stainless Steel, for JIN Type- Humerus Nail




**INS-1296T**

Instrument Trays For JIN Type - Humerus Nailing Set




**INS-1296B**

Container For JIN Type - Humerus Nailing Set



**INS-1296-000** JIN Type - Humerus Nailing Instrument Set



**INS-1296-000 JIN Type- Humerus Nailing Instrument Set**

Code	Set Consisting of	Units
<b>INS-1296-01</b>	Proximal Targeting Arm for JIN Type - Humerus Nail	1
<b>INS-1296-02</b>	Distal Targeting Arm for JIN Type - Humerus Nail	1
<b>INS-1296-03</b>	Bolt for Distal Targeting Device M6 - Hex 5mm, for JIN Type - Humerus Nail	2
<b>INS-1296-04</b>	Positioning Rod Ø5mm for JIN Type - Humerus Nail	1
<b>INS-1296-05</b>	Positioning Block for JIN Type - Humerus Nail	1
<b>INS-1296-06</b>	Trocar Ø6mm for JIN Type - Humerus Nail	1
<b>INS-1296-08</b>	Drill Sleeve 8/5.7mm for Drill Bit Ø5.5mm - JIN Type - Humerus Nail	1
<b>INS-1296-10</b>	Drill Sleeve Ø8/6mm for JIN Type - Humerus Nail	2
<b>INS-1296-11</b>	Depth Gauge Drill Sleeve Ø6mm for Drill Bit Ø2.9/3.5mm - JIN Type - Humerus Nail	1
<b>INS-1296-12</b>	Insertion Handle for JIN Type - Humerus Nail	1
<b>INS-1296-13</b>	Nail Holding Bolt for JIN Type - Humerus Nail	1
<b>INS-1296-15</b>	Step Drill Bit Plain Shank/Jacob Chuck End Ø2.9/3.5mm x Length 250mm for JIN Type - Humerus Nail	2
<b>INS-1296-16</b>	Distal Targeting Device for Length 180mm To 320mm - JIN Type - Humerus Nail	1
<b>INS-1296-07</b>	Drill Bit Plain Shank/Jacob Chuck End Ø5.5mm x Length 300mm, for JIN Type - Humerus Nail	1
<b>INS-1296-17</b>	Cannulated Awl for JIN Type - Humerus Nail	1
<b>INS-1296-18</b>	Combination Wrench 12mm for JIN Type - Humerus Nail	1
<b>INS-1296-19</b>	Handle For Reamer (AO Connection) for JIN Type - Humerus Nail	1
<b>INS-1296-20</b>	Allen Key, Hex 5mm for JIN Type - Humerus Nail	1
<b>INS-1296-21</b>	T-Handle Screwdriver, Hex 3.5mm, for JIN Type - Humerus Nail	1
<b>INS-1296-22</b>	Ø5.5mm T-Handle Drill Bit with Flat End for JIN Type - Humerus Nail	1
<b>INS-1296-23</b>	T-Handle Tap for Ø3.5mm Locking Bolt - JIN Type - Humerus Nail	1
<b>INS-1296-24</b>	T-Handle Screwdriver, Hex 2.5mm, for JIN Type - Humerus Nail	1
<b>INS-1296-27</b>	Fixed Reamer Ø8.0mm x Length 400mm with Quick Coupling End for JIN Type - Humerus Nail	1
<b>INS-1296-28</b>	Fixed Reamer Ø9.0mm x Length 400mm with Quick Coupling End for JIN Type - Humerus Nail	1
<b>INS-1296-29</b>	Fixed Reamer Ø10.0mm x Length 400mm with Quick Coupling End for JIN Type - Humerus Nail	1
<b>7-023-04</b>	Flexible Reamer Ø6.5mm for JIN Type - Humerus Nail	1
<b>7-023-05</b>	Flexible Reamer Ø7.0mm for JIN Type - Humerus Nail	1
<b>7-023-06</b>	Flexible Reamer Ø8.0mm for JIN Type - Humerus Nail	1
<b>7-023-07</b>	Flexible Reamer Ø9.0mm for JIN Type - Humerus Nail	1
<b>7-023-13</b>	Flexible Reamer Ø10.0mm for JIN Type - Humerus Nail	1
<b>7-023-01</b>	Protection Sleeve for Entry Reamer - JIN Type - Humerus Nail	1
<b>7-023-03</b>	Quick Coupling Handle for Jacob Connection for JIN Type - Humerus Nail	1
<b>7-023-08</b>	Radiographic Ruler for JIN Type - Humerus Nail	1
<b>7-023-15</b>	Guide Wire with Threaded Tip, Ø2.5mm x Thread Length 10mm x Length 200mm	3

Code	Set Consisting of	Units
<b>7-023-14</b>	Extraction Handle for JIN Type- Humerus Nail	1
<b>INS-1296-30</b>	Extraction Rod with Ram for JIN Type- Humerus Nail	1
<b>7-023-02</b>	Entry Reamer for JIN Type- Humerus Nail	1
<b>7-023-09</b>	Wire Pusher for JIN Type- Humerus Nail	1
<b>7-023-11</b>	Reduction Device for JIN Type- Humerus Nail	1
<b>7-023-12</b>	Plastic Exchange Tube for JIN Type- Humerus Nail	1
<b>INS-1296-32</b>	Guide Wire with Ball Tip Ø2.5/3.2mm x Length 1000mm - Stainless Steel, for JIN Type- Humerus Nail	1
<b>7-023-16</b>	Guide Wire without Ball Tip Ø2.5mm x Length 1000mm - Stainless Steel, for JIN Type- Humerus Nail	1
<b>INS-1296T</b>	Instrument Trays For JIN Type - Humerus Nailing Set	3
<b>INS-1296B</b>	Container For JIN Type - Humerus Nailing Set	1





**USA**

Auxein Inc.  
1500 Nw 89th Court, Suite 107-108  
Doral, Florida 33172  
Tel: +1 305 395 6062  
E Fax: +1 305 395 6262  
Email: USoffice@auxein.com

**MEXICO**

Auxein México, S.A. de C.V.  
Tepic 139 int 801, Colonia Roma Sur,  
Alcaldía Cuauhtémoc, CDMX,  
México, C.P. 06760  
Tel: +521 55 7261 0318  
Email: info@auxein.mx

**INDIA**

Auxein Medical Pvt. Ltd.  
Plot No. 168-169-170, Phase-4,  
Kundli Industrial Area,  
HSIIDC, Sector-57, Sonapat - 131028, Haryana  
Tel: +91 99106 43638 | Fax: +91 86077 70197  
Email: info@auxein.com