


METATARSAL

Shortening System



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75 countries worldwide.

Our Achievements





METATARSAL


Shortening System

INTRODUCTION

The Auxein's Metatarsal Shortening System was designed to provide fixation for an in-line, parallel shortening osteotomy of the lesser metatarsals. The system was designed to help eliminate potential complications that may arise from alternative osteotomies, such as post-operative floating toe, malunion of the metatarsal head, an altered distal articular surface, opening of the joint capsule, and/or the effects of an open capsulotomy.



Metatarsal Shortening Implant



Dia	Item Code
Ø2.7mm	12-019-2.7TI
Ø3.5mm	12-019-3.5TI
Ø4.2mm	12-019-4.2TI
Ø5.0mm	12-019-5.0TI

Dia	Item Code
Ø2.7mm	12-019-2.7TI-S
Ø3.5mm	12-019-3.5TI-S
Ø4.2mm	12-019-4.2TI-S
Ø5.0mm	12-019-5.0TI-S

2.0mm Cortical Screw, Self-Tapping, Star Head



Length (mm)	Item Code
6mm	TI-102.206
7mm	TI-102.207
8mm	TI-102.208
9mm	TI-102.209
10mm	TI-102.210
11mm	TI-102.211
12mm	TI-102.212
13mm	TI-102.213
14mm	TI-102.214
15mm	TI-102.215
16mm	TI-102.216
18mm	TI-102.218
20mm	TI-102.220
22mm	TI-102.222
24mm	TI-102.224
26mm	TI-102.226
28mm	TI-102.228
30mm	TI-102.230
32mm	TI-102.232
34mm	TI-102.234
36mm	TI-102.236
38mm	TI-102.238
40mm	TI-102.240

Length (mm)	Item Code
6mm	TI-102.206-S
7mm	TI-102.207-S
8mm	TI-102.208-S
9mm	TI-102.209-S
10mm	TI-102.210-S
11mm	TI-102.211-S
12mm	TI-102.212-S
13mm	TI-102.213-S
14mm	TI-102.214-S
15mm	TI-102.215-S
16mm	TI-102.216-S
18mm	TI-102.218-S
20mm	TI-102.220-S
22mm	TI-102.222-S
24mm	TI-102.224-S
26mm	TI-102.226-S
28mm	TI-102.228-S
30mm	TI-102.230-S
32mm	TI-102.232-S
34mm	TI-102.234-S
36mm	TI-102.236-S
38mm	TI-102.238-S
40mm	TI-102.240-S

7-161-01 Implant Inserter Ø2.7 mm for Metatarsal Shortening System



7-161-02 Implant Inserter Ø3.5 mm for Metatarsal Shortening System



7-161-03 Implant Inserter Ø4.2 mm for Metatarsal Shortening System




7-161-04 Implant Inserter Ø5.0 mm for Metatarsal Shortening System



7-161-05 Implant Sizers Ø2.7 mm for Metatarsal Shortening System




7-161-06 Implant Sizers Ø3.5 mm for Metatarsal Shortening System




7-161-07

Implant Sizers Ø4.2 mm for Metatarsal Shortening System


7-161-08

Implant Sizers Ø5.0 mm for Metatarsal Shortening System


7-161-09

Pilot Hole Wire Ø2.25 mm for Metatarsal Shortening System


7-161-10

Pilot Hole Wire Ø2.75 mm for Metatarsal Shortening System


7-161-11


K-wire Ø1.30 mm x 80mm for Metatarsal Shortening System


7-161-12


K-wire Ø1.10 mm x 80mm for Metatarsal Shortening System




7-161-13 Wire Guide for Metatarsal Shortening System




7-161-14 Depth Gauge for Metatarsal Shortening System



7-161-15 Cut Guide for 3 mm, 5 mm & 7 mm cuts




7-161-16 Cut Guide for 4 mm, 6 mm & 8 mm cuts



7-161-17 K-wire Ruler / Sizer for Metatarsal Shortening System



7-161-18 Shmogler Metatarsal Elevator for Metatarsal Shortening System



7-161-19 Slot Wrench for Metatarsal Shortening System



7-161-20 Mini Skin Hook for Metatarsal Shortening System




7-161-21 Ø2.7 mm Reamer for Metatarsal Shortening System



7-161-22 Screw Driver Shaft, T6 with AO Quick Coupling




7-161-23 Screw Holding Sleeve for 2.0mm Cortical Screw



7-161-24 Mini Pin Distractor for Metatarsal Shortening System




7-161-25 Mini Pin Compressor for Metatarsal Shortening System




7-161-26 Ratcheting Mini-AO Handle for Metatarsal Shortening System




7-161-27 Implant Caddy for Shortening Instrument System



7-161-28 Lid for Metatarsal Shortening Instrument System



7-161-29 Container for Metatarsal Shortening Instrument System



Optional

1138-39

Star Screwdriver for 2.0mm Screws - Hand & Foot System



7-164 : Intramedullary Lapidus Nail Instrument Set

Code	Description	Qty.
7-161-01	Implant Inserter Ø2.7 mm for Metatarsal Shortening System	1
7-161-02	Implant Inserter Ø3.5 mm for Metatarsal Shortening System	1
7-161-03	Implant Inserter Ø4.2 mm for Metatarsal Shortening System	1
7-161-04	Implant Inserter Ø5.0 mm for Metatarsal Shortening System	1
7-161-05	Implant Sizers Ø2.7 mm for Metatarsal Shortening System	1
7-161-06	Implant Sizers Ø3.5 mm for Metatarsal Shortening System	1
7-161-07	Implant Sizers Ø4.2 mm for Metatarsal Shortening System	1
7-161-08	Implant Sizers Ø5.0 mm for Metatarsal Shortening System	1
7-161-09	Pilot Hole Wire Ø2.25 mm for Metatarsal Shortening System	2
7-161-10	Pilot Hole Wire Ø2.75 mm for Metatarsal Shortening System	2
7-161-11	K-wire Ø1.30 mm x 80mm for Metatarsal Shortening System	4
7-161-12	K-wire Ø1.10 mm x 80mm for Metatarsal Shortening System	4
7-161-13	Wire Guide for Metatarsal Shortening System	1
7-161-14	Depth Gauge for Metatarsal Shortening System	1
7-161-15	Cut Guide for 3 mm, 5 mm & 7 mm cuts	1
7-161-16	Cut Guide for 4 mm, 6 mm & 8 mm cuts	1
7-161-17	K-wire Ruler / Sizer for Metatarsal Shortening System	1
7-161-18	Shmogler Metatarsal Elevator for Metatarsal Shortening System	1
7-161-19	Slot Wrench for Metatarsal Shortening System	1
7-161-20	Mini Skin Hook for Metatarsal Shortening System	1
7-161-21	Ø2.7 mm Reamer for Metatarsal Shortening System	1
7-161-22	Screw Driver Shaft, T6 with AO Quick Coupling	1
7-161-23	Screw Holding Sleeve for 2.0mm Cortical Screw	1
7-161-24	Mini Pin Distractor for Metatarsal Shortening System	1
7-161-25	Mini Pin Compressor for Metatarsal Shortening System	1
7-161-26	Ratcheting Mini-AO Handle for Metatarsal Shortening System	1
7-161-27	Implant Caddy for Shortening Instrument System	1
7-161-28	Lid for Metatarsal Shortening Instrument System	1
7-161-29	Container for Metatarsal Shortening Instrument System	
	Optional	
1138-39	Star Screwdriver for 2.0mm Screws - Hand & Foot System	1

**USA**

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com