


RETRO DISTAL TIBIA

RETRO DISTAL TIBIA

Nailing System



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 75+ countries worldwide.


Our Achievements



Orthopaedic Implants

INTRODUCTION

Auxein's Retro Distal Tibial Nail help to solve closed/open distal tibial fractures or one-fifth of distal tibial fractures and the metaphyseal fractures.



Type II + Blue
Left



Type II + Golden
Right



End Cap for
Retro Distal Tibia Nail



Ø3.5mm Cancellous Screw for
Retro Distal Tibia Nail



Ø4.0mm Cancellous Screw for
Retro Distal Tibia Nail



Ø3.5mm Cortical Screw for
Retro Distal Tibia Nail, Titanium



Ø4.0mm Cortical Screw for
Retro Distal Tibia Nail



Indication

- Closed and open distal tibial fractures or one-fifth of the distal tibial fractures and the metaphyseal fractures.
- Distal tibial fractures, simply involves the ankle joint

Features

- Three distal divergent locking screws, for enhanced support and fixation.
- Two proximal angled locking screws, to prevent fracture rotation and displacement.
- Special anatomical contour design for optimal nail position.
- Distal flat cutting design minimizes insertion challenges.
- Proximal Tapered end cut design, reduce irritation to soft tissues of medial malleolus


Retro Distal Tibia Nail




Code	Direction	Dia x Length
1-050-110LTI	Left	Ø7mm x 110mm
1-050-110RTI	Right	Ø7mm x 110mm
1-051-120LTI	Left	Ø8mm x 120mm
1-051-120RTI	Right	Ø8mm x 120mm
1-051-140LTI	Left	Ø8mm x 140mm
1-051-140RTI	Right	Ø8mm x 140mm

STERILE

Code	Direction	Dia x Length
1-050-110LTI-S	Left	Ø7mm x 110mm
1-050-110RTI-S	Right	Ø7mm x 110mm
1-051-120LTI-S	Left	Ø8mm x 120mm
1-051-120RTI-S	Right	Ø8mm x 120mm
1-051-140LTI-S	Left	Ø8mm x 140mm
1-051-140RTI-S	Right	Ø8mm x 140mm



Type II + Blue
Left



Type II + Golden
Right

End Cap for Retro Distal Tibia Nail



Code	Length
1-056-00TI	0mm
1-056-05TI	5mm

STERILE

Code	Length
1-056-00TI-S	0mm
1-056-05TI-S	5mm



Ø3.5mm Cancellous Screw for Retro Distal Tibia Nail



Code	Length
1-052-30TI	30mm
1-052-32TI	32mm
1-052-34TI	34mm
1-052-36TI	36mm
1-052-38TI	38mm
1-052-40TI	40mm
1-052-42TI	42mm
1-052-44TI	44mm
1-052-46TI	46mm
1-052-48TI	48mm
1-052-50TI	50mm
1-052-52TI	52mm

STERILE

Code	Length
1-052-30TI-S	30mm
1-052-32TI-S	32mm
1-052-34TI-S	34mm
1-052-36TI-S	36mm
1-052-38TI-S	38mm
1-052-40TI-S	40mm
1-052-42TI-S	42mm
1-052-44TI-S	44mm
1-052-46TI-S	46mm
1-052-48TI-S	48mm
1-052-50TI-S	50mm
1-052-52TI-S	52mm



Ø4.0mm Cancellous Screw for Retro Distal Tibia Nail



Code	Length
1-053-30TI	30mm
1-053-32TI	32mm
1-053-34TI	34mm
1-053-36TI	36mm
1-053-38TI	38mm
1-053-40TI	40mm
1-053-42TI	42mm
1-053-44TI	44mm
1-053-46TI	46mm
1-053-48TI	48mm
1-053-50TI	50mm
1-053-52TI	52mm

STERILE

Code	Length
1-053-30TI-S	30mm
1-053-32TI-S	32mm
1-053-34TI-S	34mm
1-053-36TI-S	36mm
1-053-38TI-S	38mm
1-053-40TI-S	40mm
1-053-42TI-S	42mm
1-053-44TI-S	44mm
1-053-46TI-S	46mm
1-053-48TI-S	48mm
1-053-50TI-S	50mm
1-053-52TI-S	52mm



Ø3.5mm Cortical Screw for Retro Distal Tibia Nail, Titanium



STERILE

Code	Length
1-054-28TI	28mm
1-054-30TI	30mm
1-054-32TI	32mm
1-054-34TI	34mm
1-054-36TI	36mm
1-054-38TI	38mm
1-054-40TI	40mm
1-054-42TI	42mm
1-054-44TI	44mm
1-054-46TI	46mm
1-054-48TI	48mm
1-054-50TI	50mm
1-054-52TI	52mm

Code	Length
1-054-28TI-S	28mm
1-054-30TI-S	30mm
1-054-32TI-S	32mm
1-054-34TI-S	34mm
1-054-36TI-S	36mm
1-054-38TI-S	38mm
1-054-40TI-S	40mm
1-054-42TI-S	42mm
1-054-44TI-S	44mm
1-054-46TI-S	46mm
1-054-48TI-S	48mm
1-054-50TI-S	50mm
1-054-52TI-S	52mm



Ø4.0mm Cortical Screw for Retro Distal Tibia Nail



STERILE

Code	Length
1-055-28TI	28mm
1-055-30TI	30mm
1-055-32TI	32mm
1-055-34TI	34mm
1-055-36TI	36mm
1-055-38TI	38mm
1-055-40TI	40mm
1-055-42TI	42mm
1-055-44TI	44mm
1-055-46TI	46mm
1-055-48TI	48mm
1-055-50TI	50mm
1-055-52TI	52mm

Code	Length
1-055-28TI-S	28mm
1-055-30TI-S	30mm
1-055-32TI-S	32mm
1-055-34TI-S	34mm
1-055-36TI-S	36mm
1-055-38TI-S	38mm
1-055-40TI-S	40mm
1-055-42TI-S	42mm
1-055-44TI-S	44mm
1-055-46TI-S	46mm
1-055-48TI-S	48mm
1-055-50TI-S	50mm
1-055-52TI-S	52mm



7-144-01

Aiming Jig Ø7.0mm X 110, for Retro Distal Tibia Nail



7-144-02

Aiming Jig Ø8.0mm X 120, for Retro Distal Tibia Nail




7-144-03

Aiming Jig Ø8.0mmX140, for Retro Distal Distal Tibia Nail



7-144-04

Positioning Arm Coupler for Retro Distal Tibia Nail



7-144-05

Positioning Arm for Retro Distal Tibia Nail



7-144-06

Outer Sleeve Ø8.0-5.6mm-I for Retro Distal Tibia Nail



7-144-07

Cannulated Sleeve for Drill Bit Ø2.8mm for Retro Distal Tibia Nail



7-144-08

Trocar Ø2.8mm x Length 140mm for Retro Distal Tibia Nail



7-144-09

Outer Sleeve Ø8.0-5.6mm-II for Retro Distal Tibia Nail



7-144-10

Cannulated Sleeve for Drill Bit Ø3.2mm for Retro Distal Tibia Nail



7-144-11

Trocar Ø3.2mm x Length 140mm for Retro Distal Tibia Nail



7-144-12

K-Wire Ø2.0mm x Length 180mm, for Retro Distal Tibia Nail



7-144-13 Guide Wire Holder for Retro Distal Tibia Nail



7-144-14 Entry Guide Wire Ø10.35-Ø2.2mm for Retro Distal Tibia Nail



7-144-15 Protection Sleeve Ø10.5mm for Retro Distal Tibia Nail



7-144-16 Drill Bit Ø2.8 x 220mm Long, Jacob Coupler for Retro Distal Tibia Nail



7-144-17 Drill Bit Ø3.2 x 220mm Long, Jacob Coupler for Retro Distal Tibia Nail



7-144-18 Reamer Ø7.5mm, Jacob Coupling for Retro Distal Tibia Nail



7-144-19

Reamer Ø8.5mm, Jacob Coupling for Retro Distal Tibia Nail



7-144-20

Cannulated Entry Reamer Ø7.5mm, Jacob Coupler for Retro Distal Tibia Nail



7-144-21

Cannulated Entry Reamer Ø8.5mm, Jacob Coupler for Retro Distal Tibia Nail



7-144-22

Trial for Ø7x110, Jacob Coupler for Retro Distal Tibia Nail



7-144-23

Trial for Ø8x120, Jacob Coupler for Retro Distal Tibia Nail



7-144-24

Trial for Ø8x140, Jacob Coupler for Retro Distal Tibia Nail



7-144-25


Depth Gauge, Curved for Retro Distal Tibia Nail



7-144-26 Depth Gauge, Straight for Retro Distal Tibia Nail




7-144-27 Nail Holding Bolt for Retro Distal Tibia Nail



7-144-28 Hammer Guide Shaft for Retro Distal Tibia Nail



7-144-29 Combined Hammer for Retro Distal Tibia Nail



7-144-30 L Handle Pin Wrench Ø4.5mm for Impactor - Retro Distal Tibia Nail



7-144-31 Screwdriver Shaft T20 Jacob for Retro Distal Tibia Nail



7-144-32 Screwdriver T20 for Retro Distal Tibia Nail




7-144-33

T-Handle Quick Coupling for Retro Distal Tibia Nail

**7-144-34**

Instrument Box for Retro Distal Tibia Nail

**1-069**

Retro Distal Tibia Nail Implants Box With Screw Caddy



7-144 Retro Distal Tibia Nailing Instrument Set

Codes	Set Consisting of:	Units
7-144-01	Aiming Jig Ø7.0mm X 110, for Retro Distal Tibia Nail	1
7-144-02	Aiming Jig Ø8.0mm X 120, for Retro Distal Tibia Nail	1
7-144-03	Aiming Jig Ø8.0mmX140, for Retro Distal Distal Tibia Nail	1
7-144-04	Positioning Arm Coupler for Retro Distal Tibia Nail	1
7-144-05	Positioning Arm for Retro Distal Tibia Nail	1
7-144-06	Outer Sleeve Ø8.0-5.6mm-I for Retro Distal Tibia Nail	1
7-144-07	Cannulated Sleeve for Drill Bit Ø2.8mm for Retro Distal Tibia Nail	1
7-144-08	Trocar Ø2.8mm x Length 140mm for Retro Distal Tibia Nail	1
7-144-09	Outer Sleeve Ø8.0-5.6mm-II for Retro Distal Tibia Nail	1
7-144-10	Cannulated Sleeve for Drill Bit Ø3.2mm for Retro Distal Tibia Nail	1
7-144-11	Trocar Ø3.2mm x Length 140mm for Retro Distal Tibia Nail	1
7-144-12	K-Wire Ø2.0mm x Length 180mm, for Retro Distal Tibia Nail	4
7-144-13	Guide Wire Holder for Retro Distal Tibia Nail	1
7-144-14	Entry Guide Wire Ø10.35-Ø2.2mm for Retro Distal Tibia Nail	1
7-144-15	Protection Sleeve Ø10.5 mm for Retro Distal Tibia Nail	1
7-144-16	Drill Bit Ø2.8 x 220mm Long, Jacob Coupler for Retro Distal Tibia Nail	1
7-144-17	Drill Bit Ø3.2 x 220mm Long, Jacob Coupler for Retro Distal Tibia Nail	1
7-144-18	Reamer Ø7.5mm, Jacob Coupling for Retro Distal Tibia Nail	1
7-144-19	Reamer Ø8.5mm, Jacob Coupling for Retro Distal Tibia Nail	1
7-144-20	Cannulated Entry Reamer Ø7.5mm, Jacob Coupler for Retro Distal Tibia Nail	1
7-144-21	Cannulated Entry Reamer Ø8.5mm, Jacob Coupler for Retro Distal Tibia Nail	1
7-144-22	Trial for Ø7x110, Jacob Coupler for Retro Distal Tibia Nail	1
7-144-23	Trial for Ø8x120, Jacob Coupler for Retro Distal Tibia Nail	1
7-144-24	Trial for Ø8x140, Jacob Coupler for Retro Distal Tibia Nail	1
7-144-25	Depth Gauge, Curved for Retro Distal Tibia Nail	1
7-144-26	Depth Gauge, Straight for Retro Distal Tibia Nail	1
7-144-27	Nail Holding Bolt for Retro Distal Tibia Nail	3
7-144-28	Hammer Guide Shaft for Retro Distal Tibia Nail	1
7-144-29	Combined Hammer for Retro Distal Tibia Nail	1
7-144-30	L Handle Pin Wrench Ø4.5mm for Impactor - Retro Distal Tibia Nail	1
7-144-31	Screwdriver Shaft T20 Jacob for Retro Distal Tibia Nail	1
7-144-32	Screwdriver T20 for Retro Distal Tibia Nail	1
7-144-33	T-Handle Quick Coupling for Retro Distal Tibia Nail	1
7-144-34	Instrument Box for Retro Distal Tibia Nail	1
1-069	Retro Distal Tibia Nail Implants Box With Screw Caddy	1



USA

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma Sur,
Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonipat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com