





AV-WISELOCK

AV-WISELOCK™

2.7/3.5mm Small Fragment



about us

Auxein Medical is an integrated, research based, orthopaedic Implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavor to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in india and serves customers in over 90+ countries worldwide.

Our Achievements



7-012-01 Conical Extraction Screw, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System



7-012-02 Drill Bit with Quick Coupling with Stopper, Ø2.8mm x Length 165mm, for 2.7/3.5mm AV-Wiselock System



7-012-03 Star Screwdriver Shaft, T8, for 2.7/3.5mm AV-Wiselock System



7-012-04 Depth Gauge measuring upto 60mm for 2.7/3.5mm AV-Wiselock System



7-012-05 Drill Bit with Quick Coupling, Ø2.0mm x Length 140mm, for 2.7/3.5mm AV-Wiselock System



7-012-06 Bone Tap Quick Coupling for Cortical Screws, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System



7-012-07 Screwdriver Shaft, T15, for 2.7/3.5mm AV-Wiselock System



7-012-08 Threaded Drill Sleeve for Ø2.8mm Drill Bit - 2.7/3.5mm AV-Wiselock System



7-012-09 Self-Centering Double Drill Guide, Ø2.5/3.5mm, for 2.7/3.5mm AV-Wiselock System



7-012-10 Torque Limiting Attachment, 1.2Nm, for 2.7/3.5mm AV-Wiselock System



7-012-11 Torque Limiting Attachment, 1.5Nm, for 2.7/3.5mm AV-Wiselock System



7-012-12 Drill Bit with Quick Coupling End, Ø2.5mm x Length 180mm - 2.7/3.5mm AV-Wiselock System



- 7-012-13** Drill Bit with Quick Coupling End, Ø3.5mm x Length 110mm - 2.7/3.5mm AV-Wiselock System



- 7-012-14** 2.0mm Variable Angle Guide for 2.7/3.5mm AV-Wiselock System



- 7-012-15** Self-Centering Double Drill Guide, Ø2.0/2.7mm, for 2.7/3.5mm AV-Wiselock System



- 7-012-16** 2.8mm Variable Angle Double Drill Guide



- 7-012-17** Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System



- 7-012-18** Handle for Torque Limiting Attachment, for 2.7/3.5mm AV-Wiselock System



7-012-19 2.5Nm Torque Limiting Handle for 2.7/3.5mm AV-Wiselock System



7-012-20 Depth Gauge measuring upto 100mm, for 2.7/3.5mm AV-Wiselock System



7-012-21 Percutaneous Depth Gauge Adaptor for 3.5mm Screws - 2.7/3.5mm AV-Wiselock System



7-012-22 Kirschner Wire with One Sided Trocar Tip Ø1.6mm x Length 150mm




7-012-23 Periarticular Reduction Forcep, Small, for 2.7/3.5mm AV-Wiselock System



7-012-24 Periarticular Reduction Forcep, Medium, for 2.7/3.5mm AV-Wiselock System



7-012-25 Large Compression Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System



7-012-26 Large Distraction Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System



7-012-27 3.5mm Post for AV-Wiselock Locking Hole



7-012-28 2.8mm Compression Wire, 10mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-29 2.8mm Compression Wire, 15mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-30 2.8mm Compression Wire, 20mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-31 2.8mm Compression Wire, 25mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-32 2.8mm Compression Wire, 30mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-33 2.8mm Compression Wire, 35mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-34 2.8mm Compression Wire, 40mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-35 2.8mm Compression Wire, 45mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-36 2.8mm Compression Wire, 50mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-37 2.8mm Compression Wire, 55mm Thread for 2.7/3.5mm AV-Wiselock System



7-012-38 2.8mm Compression Wire, 60mm Thread for 2.7/3.5mm AV-Wiselock System




7-012-39 Guide Wire, Ø2.0mm x Length 150mm



7-012-40 Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 10mm x Length 150mm



7-012-41 Threaded Drill Sleeve for Drill Bit 2.0mm - 2.7/3.5mm AV-Wiselock System



7-012-42 Tap, HA 2.7mm for 2.7/3.5mm AV-Wiselock System



7-012-43 Screwdriver Star Head, T8, for 2.7/3.5mm AV-Wiselock System



7-012-44 T-Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System



7-012-45 Self-Centering Bone Holding Forcep, Speed Lock, 190mm for 2.7/3.5mm AV-Wiselock System



7-012-46 Reduction Forcep, Serrated Jaws, Speed Lock, 160mm for 2.7/3.5mm AV-Wiselock System



7-012-47 Reduction Forcep, Pointed, Ratchet Lock, 180mm for 2.7/3.5mm AV-Wiselock System



7-012-48 Bending Iron 1 for 2.7/3.5mm AV-Wiselock System



7-012-49 Bending Iron 2 for 2.7/3.5mm AV-Wiselock System



7-012-50 Ø1.2mm Kirschner Wire for 2.7/3.5mm AV-Wiselock System



7-012-51 Hohmann Retractor, 6.5mm for 2.7/3.5mm AV-Wiselock System



7-012-52 Hohmann Retractor, 8.5mm for 2.7/3.5mm AV-Wiselock System



7-012-53 Hohmann Retractor, 15.5mm for 2.7/3.5mm AV-Wiselock System



7-012-54 Periosteal Elevator with Silicon Handle, Curved, 6mm for 2.7/3.5mm AV-Wiselock System



7-012-55 Screw Holding Forcep for 2.7/3.5mm AV-Wiselock System



7-012-56 Countersink for 2.7/3.5mm AV-Wiselock System



7-012-57 Screw Holding Sleeve for 2.7mm Cortical Screw



7-012-60 Screw Holding Sleeve for 2.7mm Locking Screw



7-012-61 Screw Holding Sleeve for 3.5mm Cortical Screw



7-012-62 Screw Holding Sleeve for 3.5mm Locking Screw




7-050-05


T- Handle for Threaded Drill Sleeve


7-012-58

Instrument Trays for 2.7/3.5mm AV-Wiselock System


7-012-59

Container for 2.7/3.5mm AV-Wiselock System Instrument Set



7-158-01 Pheelos Aiming Device with Nut for AV-Wiselock System

3300-02 Centering Sleeve, for Pheelos Aiming Device

7-158-02 Drill Sleeve, for Pheelos Aiming Device for AV-Wiselock System

3300-04 Centering Sleeve for Kirschner Wire Ø 1.6mm

3300-05 Pheelos Direct Measuring Device

451-1.6-150 Kirschner Wire with Trocar Tip Both End Ø1.6mm x Length 150mm


7-158 PHEELOS Instrument Set for AV-Wiselock System

Codes	Set consisting of:	Units
7-158-01	Pheelos Aiming Device with Nut for AV-Wiselock System	1
7-158-02	Drill Sleeve, for Pheelos Aiming Device for AV-Wiselock System	1
3300-02	Centering Sleeve, for Pheelos Aiming Device	1
3300-04	Centering Sleeve for Kirschner Wire Ø 1.6mm	1
3300-05	Pheelos Direct Measuring Device	1
451-1.6-150	Kirschner Wire with Trocar Tip Both End Ø1.6mm x Length 150mm	1

7-012: 2.7/3.5mm AV-Wiselock System Instrument Set



7-012: 2.7/3.5mm AV-Wiselock System Instrument Set

Code	Set Consisting of	Units
7-012-01	Conical Extraction Screw, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-02	Drill Bit with Quick Coupling with Stopper, Ø2.8mm x Length 165mm, for 2.7/3.5mm AV-Wiselock System	2
7-012-03	Star Screwdriver Shaft, T8, for 2.7/3.5mm AV-Wiselock System	2
7-012-04	Depth Gauge measuring upto 60mm for 2.7/3.5mm AV-Wiselock System	1
7-012-05	Drill Bit with Quick Coupling, Ø2.0mm x Length 140mm, for 2.7/3.5mm AV-Wiselock System	2
7-012-06	Bone Tap Quick Coupling for Cortical Screws, Ø3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-07	Screwdriver Shaft, T15, for 2.7/3.5mm AV-Wiselock System	2
7-012-08	Threaded Drill Sleeve for Ø2.8mm Drill Bit - 2.7/3.5mm AV-Wiselock System	3
7-012-09	Self-Centering Double Drill Guide, Ø2.5/3.5mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-10	Torque Limiting Attachment, 1.2Nm, for 2.7/3.5mm AV-Wiselock System	1
7-012-11	Torque Limiting Attachment, 1.5Nm, for 2.7/3.5mm AV-Wiselock System	1
7-012-12	Drill Bit with Quick Coupling End, Ø2.5mm x Length 180mm - 2.7/3.5mm AV-Wiselock System	2
7-012-13	Drill Bit with Quick Coupling End, Ø3.5mm x Length 110mm - 2.7/3.5mm AV-Wiselock System	1
7-012-14	2.0mm Variable Angle Guide for 2.7/3.5mm AV-Wiselock System	1
7-012-15	Self-Centering Double Drill Guide, Ø2.0/2.7mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-16	2.8mm Variable Angle Double Drill Guide	1
7-012-17	Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System	1
7-012-18	Handle for Torque Limiting Attachment, for 2.7/3.5mm AV-Wiselock System	1
7-012-19	2.5Nm Torque Limiting Handle for 2.7/3.5mm AV-Wiselock System	1
7-012-20	Depth Gauge measuring upto 100mm, for 2.7/3.5mm AV-Wiselock System	1
7-012-21	Percutaneous Depth Gauge Adaptor for 3.5mm Screws - 2.7/3.5mm AV-Wiselock System	1
7-012-22	Kirschner Wire with One Sided Trocar Tip Ø1.6mm x Length 150mm	5
7-012-23	Periarticular Reduction Forcep, Small, for 2.7/3.5mm AV-Wiselock System	1
7-012-24	Periarticular Reduction Forcep, Medium, for 2.7/3.5mm AV-Wiselock System	1
7-012-25	Large Compression Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System	1
7-012-26	Large Distraction Forcep with Speed Lock for 2.7/3.5mm AV-Wiselock System	1
7-012-27	3.5mm Post for AV-Wiselock Locking Hole	2
7-012-28	2.8mm Compression Wire, 10mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-29	2.8mm Compression Wire, 15mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-30	2.8mm Compression Wire, 20mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-31	2.8mm Compression Wire, 25mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-32	2.8mm Compression Wire, 30mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-33	2.8mm Compression Wire, 35mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-34	2.8mm Compression Wire, 40mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-35	2.8mm Compression Wire, 45mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-36	2.8mm Compression Wire, 50mm Thread for 2.7/3.5mm AV-Wiselock System	2

Code	Set Consisting of	Units
7-012-37	2.8mm Compression Wire, 55mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-38	2.8mm Compression Wire, 60mm Thread for 2.7/3.5mm AV-Wiselock System	2
7-012-39	Guide Wire, Ø2.0mm x Length 150mm	5
7-012-40	Guide Wire with Threaded Tip, Ø2.0mm x Thread Length 10mm x Length 150mm	2
7-012-41	Threaded Drill Sleeve for Drill Bit 2.0mm - 2.7/3.5mm AV-Wiselock System	3
7-012-42	Tap, HA 2.7mm for 2.7/3.5mm AV-Wiselock System	1
7-012-43	Screwdriver Star Head, T8, for 2.7/3.5mm AV-Wiselock System	1
7-012-44	T-Handle with Quick Coupling for 2.7/3.5mm AV-Wiselock System	1
7-012-45	Self-Centering Bone Holding Forcep, Speed Lock, 190mm for 2.7/3.5mm AV-Wiselock System	2
7-012-46	Reduction Forcep, Serrated Jaws, Speed Lock, 160mm for 2.7/3.5mm AV-Wiselock System	1
7-012-47	Reduction Forcep, Pointed, Ratchet Lock, 180mm for 2.7/3.5mm AV-Wiselock System	1
7-012-48	Bending Iron 1 for 2.7/3.5mm AV-Wiselock System	1
7-012-49	Bending Iron 2 for 2.7/3.5mm AV-Wiselock System	1
7-012-50	Ø1.2mm Kirschner Wire for 2.7/3.5mm AV-Wiselock System	3
7-012-51	Hohmann Retractor, 6.5mm for 2.7/3.5mm AV-Wiselock System	1
7-012-52	Hohmann Retractor, 8.5mm for 2.7/3.5mm AV-Wiselock System	1
7-012-53	Hohmann Retractor, 15.5mm for 2.7/3.5mm AV-Wiselock System	2
7-012-54	Periosteal Elevator with Silicon Handle, Curved, 6mm for 2.7/3.5mm AV-Wiselock System	1
7-012-55	Screw Holding Forcep for 2.7/3.5mm AV-Wiselock System	1
7-012-56	Countersink for 2.7/3.5mm AV-Wiselock System	1
7-012-57	Screw Holding Sleeve for 2.7mm Cortical Screw	1
7-012-60	Screw Holding Sleeve for 2.7mm Locking Screw	1
7-012-61	Screw Holding Sleeve for 3.5mm Cortical Screw	1
7-012-62	Screw Holding Sleeve for 3.5mm Locking Screw	1
7-050-05	T- Handle for Threaded Drill Sleeve	1
7-012-58	Instrument Trays for 2.7/3.5mm AV-Wiselock System	3
7-012-59	Container for 2.7/3.5mm AV-Wiselock System Instrument Set	1

**USA**

Auxein Inc.
1500 Nw 89th Court, Suite 107-108
Doral, Florida 33172
Tel: +1 305 395 6062
E Fax: +1 305 395 6262
Email: USoffice@auxein.com

MEXICO

Auxein México, S.A. de C.V.
Tepic 139 int 801, Colonia Roma
Sur, Alcaldía Cuauhtémoc, CDMX,
México, C.P. 06760
Tel: +521 55 7261 0318
Email: info@auxein.mx

INDIA

Auxein Medical Pvt. Ltd.
Plot No. 168-169-170, Phase-4,
Kundli Industrial Area,
HSIIDC, Sector-57, Sonapat - 131028, Haryana
Tel: +91 99106 43638 | Fax: +91 86077 70197
Email: info@auxein.com